

Product datasheet for **RG225772**

IRAK4 (NM_001114182) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IRAK4 (NM_001114182) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	IRAK4
Synonyms:	IMD67; IPD1; IRAK-4; NY-REN-64; REN64
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG225772 representing NM_001114182
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGAACAAACCCATAACACCATCAACATATGTGCGCTGCCTCAATGTTGGACTAATTAGGAAGCTGTCAG
 ATTTTATTGATCCTCAAGAAGGATGGAAGAAGTTAGCTGTAGCTATTAACCAACCATCTGGTGATGATAG
 ATACAATCAGTTTACATAAAGGAGATTTGAAGCATTACTTCAAACCTGGAAAAAGTCCCACTTCTGAATTA
 CTGTTTACTGGGGCACCACAAATTGCACAGTTGGTGATCTTGTGGATCTTTTATCCAAAATGAATTTT
 TTGCTCCTGCGAGTCTTTTGTCCAGATGCTGTTCCAAAACCTGCTAATACTACCTACCTTAAAGAAGC
 TATAACAGTTTCAACAAAAACAGATGCCTTTCTGTGACAAAGACAGGACATTGATGACACCTGTGCAGAA
 CTTGAACAAAGCTATATGCCACCTGACTCCTCAAGTCCAGAAAAAAAAGTTTGAAGTTAGTGATACAC
 GTTTTACAGTTTTTTCATTTTATGAATTGAAGAATGTCACAAATAACTTTGATGAACGCCCATTTCTGT
 TGGTGGTAATAAAATGGGAGAGGGAGGATTTGGAGTTGTATATAAAGGCTACGTAATAACACAACCTGTG
 GCAGTGAAGAAGCTTGCAGCAATGGTTGACACTACTGAAAGAACTGAAACAGCAGTTTATCAAGAAA
 TAAAAGTAATGGCAAAGTGTCAACATGAAAACCTAGTAGAACTACTTGGTTTCTCAAGTGTGGAGATGA
 CCTCTGCTTAGTATATGTTTACATGCCTAATGGTTCATTGCTAGACAGACTCTCTTGGTGGATGGTACT
 CCACCCTTTCTTGGCACATGAGATGCAAGATTGCTCAGGGTGCAGCTAATGGCATCAATTTTCTACATG
 AAAATCATCATATTCATAGAGATATTAAGAGTGCAAAATCTTACTGGATGAAGCTTTTACTGCTAAAAT
 ATCTGACTTTGGCCTTGCACGGGCTTCTGAGAAGTTGCCAGACAGTCTGACTAGCAGAATTGTGGGA
 ACAACAGCTTATATGGCACCAGAAGCTTTCGCTGGGAGAAAATAACCCCAAATCTGATATTTACAGCTTTG
 GTGTGGTTTTACTAGAAAATAAATACTGGACTTCCAGCTGTGGATGAACACCGTGAACCTCAGTTATTGCT
 AGATATTAAGAAGAAAATTGAAGATGAAGAAAAGACAATTGAAGATTATATTGATAAAAAGATGAATGAT
 GCTGATCCACTTCAGTTGAAGCTATGACTCTGTTGCTAGTCAATGCTGCTGCATGAAAAGAAAAATAAGA
 GACCAGACATTAAGAAGTTCAACAGCTGCTGCAAGAGATGACAGCTTCT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG225772 representing NM_001114182
 Red=Cloning site Green=Tags(s)

MNKPITPSTYVRCLNVGLIRKLSDFIDPQEGWKKLAVAIKKPSGDDRYNQFHRRFEALLQTGKSPTSEL
 LFDWGTNCTVGDLDVLLIQNEFFAPASLLLPDAVPKTANTLPSKEAITVQQKQMPFCDKDRTLMTVPVQN
 LEQSYMPPDSSSPENKSLEVSDFRHSFSFYELKNVTNNFDERPISVGGNKMGEVGGVGVVYKGYVNNNTV
 AVKKLAAMVDITTEELKQQFDQEIKVMAKQHENLVELLGFSSDGDCLVYVYMPNGSLLDRLSCLDGT
 PPLSWHMRCKIAQGAANGINFLHENHHIHRDIKSANILLDEAFTAKISDFGLARASEKFAQTVMTSRIVG
 TTAYMAPEALRGEITPKSDIYSFGVLLLEIITGLPAVDEHREPQLLLDIKEEIEDEEKTIEDYIDKKMND
 ADSTSVEAMYSVASQCLHEKKNRPDIKKVQQLQEMTAS

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001114182

ORF Size: 1380 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001114182.2](#), [NP_001107654.1](#)

RefSeq Size: 2915 bp
RefSeq ORF: 1383 bp
Locus ID: 51135
UniProt ID: [Q9NWZ3](#)
Cytogenetics: 12q12
Protein Families: Druggable Genome, Protein Kinase
Protein Pathways: Apoptosis, Neurotrophin signaling pathway, Toll-like receptor signaling pathway
Gene Summary: This gene encodes a kinase that activates NF-kappaB in both the Toll-like receptor (TLR) and T-cell receptor (TCR) signaling pathways. The protein is essential for most innate immune responses. Mutations in this gene result in IRAK4 deficiency and recurrent invasive pneumococcal disease. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]

Product images:



Circular map for RG225772