

Product datasheet for **RG225847**

SLC39A14 (NM_001135153) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC39A14 (NM_001135153) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SLC39A14
Synonyms:	cig19; HCIN; HMNDYT2; LZT-Hs4; NET34; ZIP14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG225847 representing NM_001135153
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAGCTGCTGCTGCTGCACCCGGCCTCCAGAGCTGCCTCCTGCTGACCTGCTTGCTTATGGAGAA
 CCACCCCTGAGGCTCACGCTTATCCCTGGGTGCCACAGCTATCAGCGCTGCCTCCTCCTCGAGGATCT
 AATACATCGGTATGGCGAGGGTGACAGCCTCACTCTGCAGCAGCTGAAGGCCCTACTCAACCCACTGGAT
 GTGGGAGTGGGCCGGGTAATGTCACCCAGCACGTGCAAGGACACAGAACCTCTCCACGTGCTTTAGTT
 CTGGAGACCTTCTACTGCCACAATTTAGCGAGCAGTCGCGGATTGGGAGCAGCGAGCTCCAGGAGTT
 CTGCCCCACCATCTCCAGCAGCTGGATTCCCGGGCCTGCACCTCGGAGAACCAGGAAAACGAGGAGAA
 GAGCAGACGGAGGGGGCGCCAAGCGCTGTTGAAGTGTGGGGATACGGTCTCCTCTGTGTACCGTCA
 TCTCCCTCTGCTCCCTCCTGGGGCCAGCGTGGTCCCTCATGAAGAAGACCTTTACAAGAGGCTGCT
 GCTCTACTTCATAGCTCTGGCGATTGGAACCTCTACTCCAACGCCCTTCCAGCTCATCCCGGAGGCA
 TTTGGTTTCAACCTCTGGAAGATTATTATGTCTCCAAGTCTGCAGTGGTGGTTGGGGGCTTTATCTTT
 TCTTTTTACAGAGAAGATCTTGAAGATTCTTCTTAAGCAGAAAAATGAGCATCATGACACAGCCA
 TTATGCCTCTGAGTCGTTCCCTCCAAGAAGACCAGGAGGGGGTGTGGAGAAGCTGCAGAACGGG
 GACCTGGACCACATGATTCCTCAGCACTGCAGCAGTGAAGTGGACGGCAAGGCGCCCATGGTGGACGAGA
 AGGTCAATGTGGGCTCGCTCTGTGTCAGGACCTGCAGGCTTCCAGAGTGGTGTACTGGCTGAAAGG
 TGTCCGCTACTCTGATATCGGCACTCTGGCCTGGATGATCACTCTGAGCGACGGCCTCCATAATTTATC
 GATGGCCTGGCCATCGGTGCTTCTTCTACTGTGTCAGTTTTCCAAGGCATCAGCACCTCGGTGGCCATCC
 TCTGTGAGGAGTTCCACATGAGCTAGGAGACTTGTGTCATCTGCTCAACGCTGGGATGAGCATCCAACA
 AGCTCTCTTCTTCAACTTCTTCTGCTGCTGCTACCTGGTCTGGCCTTTGGCATCCTGGCCGGC
 AGCCACTTCTCTGCCAAGTGGATTTTGCCTAGCTGGAGGAATGTTCTGTATATTTCTCTGGCTGATA
 TGTTCCTGAGATGAATGAGGTCTGTCAAGAGGATGAAAGGAAGGGCAGCATCTTGATTCCATTTATCAT
 CCAGAACCTGGGCTCCTGACTGGATTACCATCATGGTGGTCTCACCATGTATTCAGGACAGATCCAG
 ATTGGG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG225847 representing NM_001135153
 Red=Cloning site Green=Tags(s)

MKLLLLHPAFQSCLLLTLGLWRTTPEAHASSLGAPAI SAASFLQDLI HRYGEGDSLTLQQLKALLNHL
 VGVGRGNVTQHVGHRNLSTCFSSGDLFTAHNSEQSRIGSSELQEF CPTILQQLDSRACTSENQENEEN
 EQTEEGRPSAVEVWGYGLLCVTVISLCSLLGASVVPFMKKT FYKRLLLYFIALAIGTLYSNALFQLIPEA
 FGFNPLEDYVSKSAVVFGGFYLF FFFTEKILKILLKQKNEHHHGHSHYASESLPSKKDQEEGVMEKLQNG
 DLDMIPQHCSEL D GKAPMVDEKVI VGSLSVQDLQASQSACYWLKGVRYSDIGTLAWMITLSDGLHNF
 DGLAIGASFTVSVFQGI STSVAILCEEFPHELGD FVILLNAGMSIQQALFFNLSACCCYLGLAFGILAG
 SHFSANWIFALAGMFLYISLADMFPENVEVCQEDERKGSILIPFIIQNLGLLTGFTIMVVL TMSGQIQ
 IG

TRTRPLE – GFP Tag – V

Restriction Sites:


Sgfl-Mlul

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001135153

ORF Size: 1476 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135153.3](#)

RefSeq Size: 4740 bp

RefSeq ORF: 1479 bp

Locus ID: 23516

UniProt ID: [Q15043](#)

Cytogenetics: 8p21.3

Protein Families: Transmembrane

Gene Summary: This gene encodes a member of the the SLC39A family of divalent metal transporters that mediates the cellular uptake of manganese, zinc, iron, and cadmium. The encoded protein contains eight transmembrane domains, a histidine-rich motif, and a metalloprotease motif, and is expressed on the plasma membrane and the endocytic vesicle membrane. It is an important transporter of nontransferrin-bound iron and a critical regulator of manganese homeostasis. Naturally occurring mutations in this gene are associated with neurodegeneration with brain iron accumulation and early-onset parkinsonism-dystonia with hypermanganesemia. [provided by RefSeq, May 2017]