

Product datasheet for **RG226031**

SLC6A6 (NM_001134367) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC6A6 (NM_001134367) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SLC6A6
Synonyms:	TAUT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG226031 representing NM_001134367
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGCCACCAAGGAGAAGCTGCAGTGTCTGAAAGATTTCCACAAGGACATCCTGAAGCCCTCACCAGGGA
 AGAGCCCAGGCACGCGGCTGAGGACGAGGCTGAGGGAAAACCTCCGAGAGGGAGAAAGTGGTCTAGCAA
 GATCGACTTTGTCTCTGTGGCTGGCGGCTTCGTGGGCTTGGGCAACGTCTGGCGCTTCCCGTACCTC
 TGCTACAAGAATGGTGGAGGTGCGTTTCTCATACCGTATTTATTTTCTGTTTGGGAGCGGCTGCCTG
 TGTTTTTCTTGGAGATCATCATAGGCCAGTACACCTCTGAAGGGGGCATCACCTGCTGGGAAAAGATCTG
 CCCCTTGTCTCTGATCGGCTATGCCTCCGTTGTAATTGTGTCCTCCTGAATGTCTACTACATCGTC
 ATCCTGGCCTGGCCACATACTACCTGTCCAGTCTTCCAGAAGGAGCTGCCCTGGGCACACTGCAACC
 ACAGCTGGAACACACCTCACTGCATGGAGGACACCATGCGCAAGAACAAGAGTGTCTGGATCACCATCAG
 CTCCACCAACTTACCTCCCCTGTCATCGAGTTCTGGGAGCGCAACGTGTGAGCTTGTCCCTGGAATC
 GACCACCCAGGCTCTCTGAAATGGGACCTCGCTCTCGCTTCTTTAGTCTGGCTAGTGTGTTTCTTCT
 GCATCTGGAAGGCGTCAAGTCCACTGGGAAGGTGCTCTACTTCACAGCCACTTTTCCATTGCCATGCT
 CCTGGTGTCTGGTCCGAGGGCTGACGCTGCCGGGCGCGGGCGCAGGCATCAAGTTCTATCTGTATCCT
 GACATCACCCGCCTTGGAGACCCACAGGTGTGGATTGACGCTGGGACTCAGATATTCTTCTCTTATGCCA
 TCTGCCTGGGGCTATGACCTCGCTGGGAGCTACAACAAGTACAAGTATAACTCGTACAGGACTGTAT
 GCTGCTGGGATGCCTGAACAGTGGTACCAAGTTTGTGTCTGGCTTCGCAATTTTTTCCATCCTGGGCTTC
 ATGGCACAAGAGCAAGGGGTGGACATTGCTGATGTGGCTGAGTCAGGTCCTGGCCTGGCCTTCATTGCCT
 ACCCAAAAGCTGTGACAATGATGCCGCTGCCACATTTTGGTCCATTCTTTTTTTTATTATGCTTCTCTT
 GCTTGGACTGGATAGCCAGTTTGTGAAGTTGAAGGACAGATCACATCCTTGGTTGATCTTTACCCATCC
 TTCTTAAGGAAGGTTATCGTCGGGAAATCTTCATCGCCTTCGTGTGTAGCATCAGCTACCTGCTGGGGC
 TGACGATGGTGACGGAGGGTGGCATGTATGTGTTTCACTCTTTGACTACTATGCAGCTAGCGGTGTATG
 CCTTTTGTGGTTGCAATCTTTGAATGTTTGTATTGCTGGATATAGGAGGTGATAACCTTTATGAT
 GGTATTGAGGACATGATTGGCTATCGGCCCGGGCCCTGGATGAAGTACAGCTGGGCTGTGATCACTCCAG
 TTCTCTGTGTTGGATGTTTCATCTTCTCGCTCGTCAAGTACGTACCCCTGACCTACAACAAAACATACGT
 GTACCCCAACTGGGCCATTGGGCTGGGCTGGAGCCTGGCCCTTCTCCATGCTCTGCGTCCCTTGGTC
 ATCGTCATCCGCCTTGCAGACTGAGGGGCCGTTCTTGTGAGAGTCAAGTACCTGCTGACCCCAAGGG
 AACCCAAACCGCTGGGCTGTGGAGCGCGAGGGAGCCACACCTTACAACCTCTCGCACCGTCATGAACGGCGC
 TCTCGTGAACCGACCCACATCATTGTGGAGACCATGATG

ACGCGTACGCGGGCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG226031 representing NM_001134367
 Red=Cloning site Green=Tags(s)


MATKEKLQCLKDFHKDILKPSPGKSPGTRPEDEAEGKPPQREKWSSKIDFVLSVAGGFVGLGNVWRFPYL
 CYKNGGGAFLIPYFIFLFGSGLPVFFLEIIIGQYTSEGGITCWEKICPLFSGIGYASVVIIVSLNNVYIV
 ILAWATYYLFQSFQKELPWAHCNHSWNTPHCMEDTMRKNKSVWITISSTNFTSPVIEFWERNVLSLSPGI
 DHPGSLKWDALCLLLVWLVCFFCIWKGVIRSTGKVVYFTATFPFAMLLVLLVRGLTLPAGAGAIKFYLYP
 DITRLEDPOQWIDAGTQIFFSYAICLGAMTSLGSYNKYKNSYRDCMLLGCLNSGTSFVSGFAIFSILGF
 MAQEQQVDIADVAESGGLAFIAYPKAVTMMPLPTFWSILFFIMLLLLGLDSQFVEVEGQITSLVDLYPS
 FLRKGYRREIFIAFVCSISYLLGLTMVTEGMYVFQLFDYYAASGVCLLWVAFFECFVIAWIYGGDNLYD
 GIEDMIGYRPGPMKYSAVITPVLCVGCIFISLVKYVPLTYNKTYVYPNWAIGLWLSLSSMLCVPLV
 IVIRLCQTEGPFLLVRVKYLLTPREPNRWAVEREGATPYNSTVMNGALVKPHTIIIVETMM

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001134367

ORF Size: 1860 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001134367.2](#)

RefSeq Size: 6512 bp


RefSeq ORF: 2166 bp

Locus ID: 6533

UniProt ID: [P31641](#)

Cytogenetics:	3p25.1
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	This gene encodes a multi-pass membrane protein that is a member of a family of sodium and chloride-ion dependent transporters. The encoded protein transports taurine and beta-alanine. There is a pseudogene for this gene on chromosome 21. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013]

Product images:



Circular map for RG226031