

Product datasheet for **RG226044**

SHTN1 (NM_001127211) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SHTN1 (NM_001127211) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SHTN1
Synonyms:	KIAA1598; shootin-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG226044 representing NM_001127211
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAACAGCTCGGACGAAGAGAAGCAGCTGCAGCTCATTACCAGTCTGAAGGAGCAAGCAATAGGCGAAT
 ATGAAGACCTTAGAGCAGAGAACCAGAAAACAAAGGAGAAGTGTGACAAAATTAGGCAAGAACGAGATGA
 AGCCGTTAAAAAACTGGAAGAATTTAGAAAAATTTCTCACATGGTCATAGAGGAAGTTAATTTATGTCAG
 AACCATCTTGAAATAGAGAAGACTTGTGCGAGAAAGTGTGAAGCTTTGGCAACAAAGCTAAATAAAGAAA
 ATAAACGTTGAAAAGAATCAGCATGTTGTACATGGCCAAAGCTGGGACCAGATGTAATACTGAAGAGAT
 AAACATTGATGATGAAGATTCGACTACAGACACAGACGGTGCCGCCGAGACTTGTGTCTCAGTACAGTGT
 CAGAAGCAAAATTAAGAACTTCGAGATCAAATTGTATCTGTTGAGGAGGAAAAGAAGATTTTAGCCATTG
 AGCTGGAAAATCTCAAGAGCAAACTCGTAGAAGTAATTGAAGAAGTAAATAAAGTTAAACAAGAAAAGAC
 TGTTTTAAATTCAGAAGTTCTTGAACAGAGAAAAGTCTTAGAAAAATGCAATAGAGTGTCCATGTTAGCT
 GTAGAAGAGTATGAGGAGATGCAAGTAAACCTGGAGCTGGAGAAGGACCTTCGAAAGAAAGCAGAGTCAT
 TTGCACAAGAGATGTTTATTGAGCAAAAACAGCTAAAGAGACAAAGCCACCTTCTGCTGCAGAGCTCCAT
 CCCTGATCAGCAGCTTTTGAAAGCTTTAGACGAAAATGCAAACTCACCCAGCACTTGAAGAAGAGAGA
 ATTCAGCATCAACAAAAGGTCAAAGAATTAGAAGAGCACTAGAAAAAGAACTCCACAAAGAAATAC
 ACAACCTCAAACAGCACTGGAGCTTCTAGAGGAAGATAAAAAGGAATTGGAATTGAAATATCAGAATTC
 TGAAGAGAAAGCCAGAAATTTAAAGCACTCTGTTGATGAACTCCAGAAACGAGTGAACAGTCTGAGAAT
 TCAGTACCTCCACCACCTCCTCCTCCACCACCACTTCCCCTCCACCTCCCAATCCTATCCGATCCCTCA
 TGTCCATGATCCGGAACGATCCCAACCCAGTGGCAGTGGTGCTAAGAAAGAAAAGGCAACTCAACCAGA
 AACCACTGAAGAAGTACAGATCTAAAGAGGCAAGCAGTTGAAGAGATGATGGATAGAATTAAGGGA
 GTTCATCTTAGACCCGTTAATCAGACAGCCAGACCGAAGACAAAGCCAGAATCTTCGAAAGGCTGCGAAA
 GTGCAGTGGATGAACTAAAGGAATACTGGGACACTTAACAAATCCACTAGTTCAAGAAGCTTAAATC
 CCTTGACCCTGAAAACAGTGAACTGAGTTAGAAAGGATTTTGCCTGCGAGAAAGGTGACAGCAGAAGCA
 GATAGCAGTAGTCCAAGTGGATATTAGCCACCTCAGAGTCCAAATCCATGCCAGTGTGGGTTCTGTAT
 CCAGTGTAAACAAAACAGCCTTGAACAAGAAAACCTCTGGAGGCAGAAATTAACAGCCCGTCCCCCAAC
 ACCTGAGCCAGGTGAAGGGCCCCGTAATTTGGAAGGATGCACAAGTTCCAAGTTACGTTTCAGCCTCCC
 AGTAGCATTGGATGCAGGAAAAATACATTGACGGTGAAAAACAAGCCGAACCAAGTTGTAGTTTTAGATC
 CTGTTTCTACACATGAACCCCAACCAAGACAGGTTGCTGAAAAAGATCCAACCAACACAAGGAGGA
 TGAAGGCGAAATTCACCCAGAAAACAAAGAAGACAGCATTGAAAACGTGAGAGAGACAGACAGCTCCAAC
 TGC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:


>RG226044 representing NM_001127211
 Red=Cloning site Green=Tags(s)

MNSSDEEKQLQLITSLKEQAIGEYEDLRAENQKTEKCDKIRQERDEAVKKLEEFQKISHMVIEEVNFMQ
 NHLEIEKTCRESAEALATKLNKENKTLKRISMLYMAKLGPDVITEEINIDDEDSTTDGAAETCVSVQC
 QKQIKELRDQIVSVQEEKKILAELENLKSLEVEVIEVNKVKQEKTVLNSEVLEQRKVLEKCNVSMILA
 VEEYEEMQVNLELEKDLRKKAEFAQEMFIEQNKLRQSHLLQSSIPDQQLLKALDENAKLTQQLLEER
 IQHQQKVKELEELENELHKEIHNKQQLLELEDKKELELKYQNSEEKARNLKHVDELQKRVNQSEN
 SVPPPPPPPPPLPPPPPNPIRSLMSMIRKSHPSGSGAKKEKATQPETTEEVTDLKRQAVEEMMDRIKKG
 VHLRPVNQTARPKTKPESSKGCESAVDELKGILGTLNKTSSRSLKSLDPENSETELERILRRRKVTAEA
 DSSSPTGILATSESKSMPVLGSVSVTKTALNKKTLAEFNPSPPPTPEPGEGPRKLEGCTSSKVTFPQP
 SSIGCRKKYIDGKQAEPPVVLDPVSTHEPQTKDQVAEKDPTQHKEDGEIQPENKEDSIENVRETDSSN
 C

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001127211

ORF Size: 1893 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001127211.3](#)

RefSeq Size: 3950 bp

RefSeq ORF: 1896 bp


Locus ID: 57698

UniProt ID: [A0MZ66](#)

Cytogenetics: 10q25.3

Gene Summary: Involved in the generation of internal asymmetric signals required for neuronal polarization and neurite outgrowth. Mediates netrin-1-induced F-actin-substrate coupling or 'clutch engagement' within the axon growth cone through activation of CDC42, RAC1 and PAK1-dependent signaling pathway, thereby converting the F-actin retrograde flow into traction forces, concomitantly with filopodium extension and axon outgrowth. Plays a role in cytoskeletal organization by regulating the subcellular localization of phosphoinositide 3-kinase (PI3K) activity at the axonal growth cone. Plays also a role in regenerative neurite outgrowth. In the developing cortex, cooperates with KIF20B to promote both the transition from the multipolar to the bipolar stage and the radial migration of cortical neurons from the ventricular zone toward the superficial layer of the neocortex. Involved in the accumulation of phosphatidylinositol 3,4,5-trisphosphate (PIP3) in the growth cone of primary hippocampal neurons.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG226044