


Product datasheet for **RG226069**

CD105 (ENG) (NM_001114753) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CD105 (ENG) (NM_001114753) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CD105
Synonyms:	END; HHT1; ORW1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
Restriction Sites:	SgfI-RsrII
Cloning Scheme:	

Cloning sites used for ORF Shutting:



Cloning sites used for ORF Shutting:

EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI
 CTATAGGCGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGATCGCCGCCGCCAGATCT

HindIII NheI RsrII MluI NotI XhoI GFP Tag
 CAAGCTTAAGCTAGCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC ---

--- GAA GAA AGA GTT TAA ACGGCCGGCCGGAGCT
 - - E E R V Stop


ACCN:	NM_001114753
ORF Size:	1974 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001114753.3
RefSeq Size:	3072 bp
RefSeq ORF:	1977 bp
Locus ID:	2022
UniProt ID:	P17813
Cytogenetics:	9q34.11
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Gene Summary:	This gene encodes a homodimeric transmembrane protein which is a major glycoprotein of the vascular endothelium. This protein is a component of the transforming growth factor beta receptor complex and it binds to the beta1 and beta3 peptides with high affinity. Mutations in this gene cause hereditary hemorrhagic telangiectasia, also known as Osler-Rendu-Weber syndrome 1, an autosomal dominant multisystemic vascular dysplasia. This gene may also be involved in preeclampsia and several types of cancer. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2013]

Product images:



Circular map for RG226069