

Product datasheet for **RG226074**

PTPN12 (NM_001131008) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTPN12 (NM_001131008) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PTPN12
Synonyms:	PTP-PEST; PTPG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG226074 representing NM_001131008
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGATATGGGAGTATAATGTTGTGATCATTGTAATGGCCTGCCGAGAATTTGAGATGGGAAGGAAAAAT
GTGAGCGCTATTGGCCTTTGTATGGAGAAGACCCATAACGTTTGCACCATTTAAAATTTCTTGTGAGGA
TGAACAAGCAAGAACAGACTACTTCATCAGGACTCTTACTTGAATTTCAAATGAATCTCGTAGGCTG
TATCAGTTTCATTATGTGAAGTGGCCAGACCATGATGTTCTTCATCATTGATTCTATTCTGGACATGA
TAAGCTTAATGAGGAAATATCAAGAACATGAAGATGTTCTATTTGATTTCATTGCAGTGCAGGCTGTGG
AAGAACAGGTGCCATTTGTCCATAGATTATACGTGGAATTTACTAAAAGCTGGGAAAAATACCAGAGGAA
TTTAATGTATTTAATTAACAAGAAATGAGAACACAAAGGCATTCTGCAGTACAAACAAAGGAGCAAT
ATGAACCTGTTTCATAGAGCTATTGCCAACTGTTTAAAAACAGCTACAACATATGAAATTCATGGAGC
TCAGAAAATTGCTGATGGAGTGAATGAAATTAACACTGAAAACATGGTCAGCTCCATAGAGCCTGAAAA
CAAGATTCCTCCTCCTCCAAAACCAAGGACCCGCAGTTGCCTTGTGAAGGGGATGCTAAGAAGAAA
TACTGCAGCCACCGAACCTCATCCAGTCCACCCATCTTGACACCTTCTCCCCCTTCAGCTTTTCCAAC
AGTCACTACTGTGTGGCAGGACAATGATAGATACCATCCAAAGCCAGTGTTCATATGGTTTCATCAGAA
CAACATTCAGCAGACCTCAACAGAACTATAGTAAATCAACAGAACTCCAGGAAAAATGAATCAACAA
TTGAACAGATAGATAAAAAATTTGGAACGAAATTAAGTTTGGAGATTAAGAAGGTCCCTCTCAAGAGGG
ACCAAAAAGTTTTGATGGGAACACACTTTTGAATAGGGGACATGCAATTAATTAATCTGCTTCCACT
TGTATAGCTGATAAAATCTCTAAGCCACAGGAATTAAGTTCAGATCTAAATGTCGGTGATACTCCAGAA
ATTCTTGTGTGGACTGCAGTGAACACAATCAAACAAGTTTCAGTTACTCCACCAGAAGAATCCAGAA
TTCAGACACACCTCCAAGGCCAGACCGCTTGCCTTGTGATGAGAAAGGACATGTAACTGGTCAATTCAT
GGACCTGAAAATGCCATACCCATACCTGATTTATCTGAAGGCAATTCCTCAGATATCAACTATCAAACTA
GGAAAATGTGAGTTTAAACCAAGTCTACAACACAAGTTGAAACACCTGATCTTGTGGATCATGATAA
CACTTCACCCTTTCAGAACACCCCTCAGTTTACTAATCCACTTCACTCTGATGACTCAGACTCAGAT
GAAAGAACTCTGATGGTGTGTGACCCAGAATAAACTAATATTTCAACAGCAAGTCCACAGTTTCTG
CTGCCACTAGTACTGAAAGCATTCTACTAGGAAAGTATTGCCAATGTCCATTGCTAGACATAATATAGC
AGGAACAACACATTCAGGTGCTGAAAAAGATGTTGATGTTAGTGAAGATTCACCTCCTCCCTACCTGAA
AGAACTCCTGAATCGTTTGTGTTAGCAAGTGAACATAATACACCTGTAAGATCGGAATGGAGTGAACCTC
AAAGTCAGGAACGATCTGAACAAAAAAGTCTGAAGGCTTGATAACCTCTGAAAAATGAGAAATGTGATCA
TCCAGCGGGAGGTATTCATATGAAATGTGCATAGAATGTCCACCTACTTTCAGTGACAAGAGAGAACAA
ATATCAGAAAATCCAACAGAAGCCACAGATATTGGTTTTGGTAATCGATGTGAAAAACCAAGGACCAA
GAGATCCACCTTCAGAATGGACA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG226074 representing NM_001131008
 Red=Cloning site Green=Tags(s)


MIWEYNVVIIVMACREFEMGRKKCERYWPLYGEDPITFAPFKISCEDEQARTDYFIRTLLEFQNESRRL
 YQFHYVNWPDHDPSSFDSILDMISLMRKYQEHEDVPICIHCSAGCGRTGAICAIDYTWLLKAGKIPEE
 FNVFNLIQEMRTQRHSAVQTKEQYELVHRAIAQLFEKQLQLYEIHGAQKIADGVNEINTENMVSSIEPEK
 QDSPPPKPPRTRSCLEVGDAKEEILQPPEPHVPPILTPSPPSAFPTVTTVWQDNDRYHPKPVLMVSSSE
 QHSADLNRNYSKSTELPGKNESTIEQIDKKLERNLSFEIKKVPLQEGPKSFDGNTLLNRGHAIKIKSASP
 CIADKISKPQELSSDLNVGDTSQNSCVDCSVTQSNKVSVTPPEESQNSDTPPRPDRPLPDEKGHVTWSFH
 GPENAIPIPDLEGNSSDINYQTRKTVSLTPSPTTQVETPDLVDHNTSPLFRTPLSFTNPLHSDSDSD
 ERNSDGAVTQKNTNISTASATVSAATSTESISTRKVLPMIARHNIAGTTHSGAEKDVDVSEDSPPPLPE
 RTPESFVLASEHNTPVRESEWSELQSQRSEQKKSEGLITSENEKCDHPAGGIHYEMCIECPPTFSKREQ
 ISENPTEATDIGGNRCGKPKGPRDPPSEWT

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001131008

ORF Size: 1983 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001131008.2](#)

RefSeq Size: 3122 bp

RefSeq ORF: 1986 bp

Locus ID: 5782

UniProt ID: [Q05209](#)

Cytogenetics: 7q11.23

Protein Families: Druggable Genome, Phosphatase

Gene Summary: The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP contains a C-terminal PEST motif, which serves as a protein-protein interaction domain, and may regulate protein intracellular half-life. This PTP was found to bind and dephosphorylate the product of the oncogene c-ABL and thus may play a role in oncogenesis. This PTP was also shown to interact with, and dephosphorylate, various products related to cytoskeletal structure and cell adhesion, such as p130 (Cas), CAKbeta/PTK2B, PSTPIP1, and paxillin. This suggests it has a regulatory role in controlling cell shape and mobility. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Oct 2008]