

## Product datasheet for **RG226099**

### Lebercilin (LCA5) (NM\_001122769) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Lebercilin (LCA5) (NM_001122769) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LCA5
Synonyms:	C6orf152
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RG226099 representing NM\_001122769  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGGGAAAGAGCAGGAAGTCCAGGTACTGATCAAGAAAGAAAGGCAGGCAAACACCATTATTCTTACT  
TATCTGATTTTGAACGCCACAGTCTTCTGGCCGATCATCGTGGTCAGTTCTTACCTGCAAGTGTAG  
GAGAAAAATCCTAAAAGACAACTTCAGATGGCCAAGTACATCACCAAGCCCTCGAAAACCAAGCCCT  
AAGGGTCTACCAAACAGAAAGGGAGTCCGAGTGGGATTTTCGCTCCAGAGCCTCAATAGAGAGCCACTTC  
GGAAAGATACTGATCTTGTACAAAACGGATTCTGTCTGCAAGACTGCTAAAAATCAATGAGTTGCAGAA  
TGAAGTATCTGAACTCCAGGTCAAGTTAGCTGAGCTGCTAAAAGAAAATAAATCTTTGAAAAGGCTTCAG  
TACAGACAGGAGAAAGCCCTGAATAAGTTTGAAGATGCCGAAAATGAAATCTCACAACTTATATTTTCGTC  
ATAACAATGAGATTACAGCACTCAAAGAACGCTTAAGAAAATCTCAAGAGAAAGAACGGGCAACTGAGAA  
AAGGGTAAAAGATACAGAAAGTGAAGTATTTAGGACAAAATTTTCTTACAGAACTGAAAGAGATCTCT  
GAAGCTAGACACCTACCTGAACGAGATGATTTGGCAAAGAACTAGTTTACGAGAGTTAAAGTTAGATG  
ACACCGAGAGAAGAAATTAAGGAGCTATCGAAAACTTGAAGTACTAACAGTTTCCACGACAGTT  
GCTTGCTGAAAGGAAAAGGGCATATGAGGCTCATGATGAAAAATAAAGTTCTTCAAAGGAGGTACAGCGA  
CTATATCACAATAAAGGAAAAGGAGAGAACTGGATATAAAAAATATATATTCTAATCGTCTGCCAA  
AGTCTCTCCAAATAAAGAGAAAGAACTTGCAATTAAGAAAAATGCTGCATGCCAGAGTGAATTTGCAGA  
CCTGTGTACAAAAGGAGTACAAACCATGGAAGACTTCAAGCCAGAAGAAATATCCTTTAACTCCAGAAACA  
ATTATGTGTTACGAAAACAAATGGGAAGAACCAGGACATCTTACTTTGGACTTGCAATCTCAAAGCAAG  
ACAGGCATGGAGAAGCAGGGATTCTAAACCCAATTATGGAAAGAGAAGAAAAATTTGTTACAGATGAAGA  
ACTCCATGTCGTAACACAGGAGGTTGAAAAGCTGGAGGATGAATGGGAAAGAGAAGAACTTGATAAAAAAG  
CAAAAAGAAAAGGCATCTTTACTGGAAAGAGAAGAAAAGCCAGAGTGGGAACTGGAAGGTACCAACTAG  
GAATGTATCCAATTCAGAATATGGATAAATGCAAGGAGAGGAAGAAGAAAGACTGAAGAGAGAAATGCT  
ACTTGCTAACTGAATGAAATTGACAGAGAACTCCAAGATTCTCGAAATCTAAAAACCCTGTTTTGCCA  
TTGTTACCTGATTTTGAATCAAACACTACACTCCCAGAGAGAAGCCCAAAACATACAGGTTCTCTGAAT  
CCTCAGAGAGATTATTTAATGGGCATCTTTGCAAGACATCAGTTTCTCAACTCCAAAAGGAGAAGGTCA  
GAATTCAGGAAATGTTAGAAGTCCAGCCTCCCCTAATGAGTTCGCATTTGGTAGCTACGTGCCTTCGTTT  
GCAAAAACATCAGAGAGGTCAAATCCATTTAGTCAAAAAAGTAGTTTTTTGGATTTCAAAGAAACAGTA  
TGGAAAAACTTAGTAAAGATGGTGTAGATTTAATTACAAGAAAAGAGAAAAAGCTAATTTGATGGAACA  
GTTATTTGGTGCCAGTGGTAGCAGCACCATTTCTCCTCAAAGCAGTGACCCAAATTTCTGTGGCTTCCAGT  
AAAGGAGACATTGACCCTCTAAATTTTCTCCCTGGGAATAAAGGCAGCAGAGATCAAGAACATGATGAAG  
ATGAAGGCTTTTTCTCAGTGAAGGAAGAAGTTTTAATCCAAATAGGCACCGATTAACATGCAGACGA  
TAAACCAGCAGTAAAAGCAGCTGATTCTGTAGAAGATGAAATTGAAGAAGTAGCACTGAGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG226099 representing NM\_001122769  
 Red=Cloning site Green=Tags(s)


MGERAGSPGTDQERKAGKHYSYLSDFETPQSSGRSSLVSSSPASVRRKNPKRQTS DGQVHHQAPRKPS  
 KGLPNRKGVRVGFERSQSLNREPLRKDLDLVTKRILSARLLKINELQNEVSELQVKLAELLKENKSLKRLQ  
 YRQEKALNKFEDAENEISQLIFRHNNEITALKERLRKSQEKERATEKRVKDESELFRTKFSLQKLKEIS  
 EARHLPERDDLAKKLVSAELKDDTERRIKELSKNLELSTNSFQRQLLAERKRAYEAHDENKVLQKEVQR  
 LYHKLKEKERELDIKNIYSNRLPKSSPNKEKELALRKNAACQSDFADLCTKGVQTMEDFKPEEYPLTPET  
 IMCYENKWEPEGHLTDLQSQKQDRHGEAGILNPIMEREKVFVTDEELHVVKQEVEKLEDEWEREELDKK  
 QKEKASLLEREKPEWETGRYQLGMYPIQNMDKLQGEERLREMLLAKLNEIDRELQDSRNLKYPVLP  
 LLPDFESKLHSPERSPKTYRFSESSERLFGHHLQDISFSTPKGEGQNSGNVRSASPNEFAFGSYVPSF  
 AKTSERSNPFQKSSFLDFQRNSMEKLSKDGVDLITRKEKKANLMEQLFGASGSSTISSKSSDPNSVASS  
 KGDIDPLNFLPGNKGRDQEHDEDEGFFLSEGRSFPNRRHRLKHADKPAVKAADSVEDEIEEVALR


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**


**ACCN:** NM\_001122769

**ORF Size:** 2091 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001122769.3](#)

**RefSeq Size:** 4542 bp

RefSeq ORF: 2094 bp

Locus ID: 167691

UniProt ID: Q86VQ0

Cytogenetics: 6q14.1

**Gene Summary:** This gene encodes a protein that is thought to be involved in centrosomal or ciliary functions. Mutations in this gene cause Leber congenital amaurosis type V. Alternatively spliced transcript variants are described. [provided by RefSeq, Oct 2009]