

Product datasheet for **RG226103**

SLC5A9 (NM_001135181) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC5A9 (NM_001135181) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SLC5A9
Synonyms:	SGLT4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG226103 representing NM_001135181
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGAGCAAGGAGCTGGCAGCAATGGGGCCTGGAGCTTCAGGGGACGGGGTCAGGACTGAGACAGCTCCAC
ACATAGCACTGGACTCCAGAGTTGGTCTGCACGCCACGACATCAGCGTGGTGGTCATCTACTTTGTCTT
CGTCATTGCTGTGGGATCTGGTCGTCCATCCGTGCAAGTCGAGGGACCATTGGCGGCTATTTCTGGCC
GGGAGGTCATGAGCTGGTGGCCAATTGGAGCATCTCTGATGTCCAGCAATGTGGGCAGTGGCTTGTTC
TCGGCCTGGTGGGACAGGGGCTGCCGGAGGCCTTGCCGTAGGTGGCTTCGAGTGGAAATGAGGAAATC
AAGGTCTGGAGGAGACAGAGGGATCCATCCAAGGTCACACGGGAGGACTGGGGTCAGGTCCAGGCAACC
TGGCTGCTCTGGCCCTTGGCTGGGTCTTCGTCCCTGTGTACATCGCAGCAGGTGTGGTCACAATGCCGC
AGTATCTGAAGAAGCGATTTGGGGCCAGAGGATCCAGGTGTACATGTCTGTCTCTCATCTCTA
CATCTTACCAAGATCTCGACTGACATCTTCTCTGGAGCCCTTTCATCCAGATGGCATTGGGCTGGAAC
CTGTACCTCTCCACAGGATCCTGCTGGTGGTACTGCCGTCTACACCATTGCAGGTGGCCCTCATGGCCG
TGATCTACACAGATGCTCTGCAGACGGTATCATGGTAGGGGGAGCCCTGGTCTCATGTTTCTGGGCTT
TCAGGACGTGGGCTGGTACCCAGGCCTGGAGCAGCGGTACAGGCAGGCCATCCCTAATGTACAGTCCCC
AACACCACCTGTACCTCCCACGGCCCGATGCTTCCACATTCTTCGGGACCCTGTGAGCGGGGACATCC
CTTGGCCAGGTCTCATTTTCGGGCTCACAGTGTGGCCACCTGGTGTGGTGCACAGACCAGGTCAATTGT
GCAGCGGTCTCTCTCGGCCAAGAGTCTGTCTCATGCCAAGGGAGGCTCCGTGCTGGGGGCTACCTGAAG
ATCCTCCCATGTTCTTCATCGTCATGCCTGGCATGATCAGCCGGGCCCTGTTCCAGACGAGGTGGGT
GCGTGGACCCTGATGTCTGCCAAAGAATCTGTGGGGCCGAGTGGGATGTTCCAACATTGCCTACCCTAA
GTTGGTCAATGGCCCTCATGCCTGTTGGTCTGCGGGGGCTGATGATTGCCGTGATCATGGCCGCTCTCATG
AGCTCACTCACCTCCATCTTCAACAGCAGCAGCACCCCTGTTACCAATTGATGTGTGGCAGCGCTCCGCA
GGAAGTCAACAGAGCAGGAGCTGATGGTGGTGGGCAGAGTGTGGTGGTCTCTGGTGTGCATCAGCAT
CCTCTGGATCCCCATCATCAAAGCTCCAACAGTGGCAGCTCTTCGACTACATCCAGGCTGTACCCAGT
TACCTGGCCCCACCCATCACCGCTCTTCTCTGCTGGCCATCTTCTGCAAGAGGGTCACAGAGCCCCGAG
CTTTCTGGGGCTCGTGTGGCCTGGGAGTGGGGCTTCTGCGTATGATCCTGGAGTTCTCATACCCAGC
GCCAGCCTGTGGGAGGTGGACCGGAGGCCAGCAGTGTGAAGGACTTCCACTACCTGTACTTTGCAATC
CTCCTCTGCGGGCTCACTGCCATCGTCATTGTCTATTGTCAAGCCTCTGTACAACCTCCATCCCTGAGGAA
AGCTCACAGCCTCACATGGTGGACTCGGAACTGCCCTCTCTGAGCTGGAGAAGGAGGCCACGAGAG
CACACCGGAGATATCCGAGAGGCCAGCCGGGAGTGCCTGCAGGAGGTGGAGCGGCAGAGAAGTCCGAGC
CTGGGCCAGGAGCAGCCTGAAGCCCAAGCAGGTCTGGGGAAAGTTGCTCTGGAGCTGGTCTGTGGGC
TCTCTGGAACACCGGAGCAGGCCCTGAGCCAGCAGAGAAGGCTGCGCTAGAACAGAAGCTGACAAGCAT
TGAGGAGGAGCCACTCTGGAGACATGTCTGCAACATCAATGCTGTCTTTTGTGGCCATCAACATCTTC
CTCTGGGGCTATTTTGGC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG226103 representing NM_001135181
 Red=Cloning site Green=Tags(s)


MSKELAAMGPGASGDGVRTETAPHIALDSRVGLHAYDISVVVIYFVFVIAVGIWSSIRASRGTIGGYFLA
 GRMSWVPIGASLMSSNVGSLFIGLAGTGAAGGLAVGGFEWNRKSRSGGDRGIHPRSHGRTGVRSQAT
 WLLLALGWVFPVYIAAGVVTMPQYLKKRFGGQRIQVYMSVLSLILYIFTKISTDIFSGALFIQMALGWN
 LYLSTGILLVVTAVYTIAGGLMAVIYTDALQTVIMVGGALVLMFLGFQDVGWYPGLEQRYRQAIPNVTP
 NTTCHLPRPDAFHILRDPVSGDIPWPGLIFGLTVLATWCWCTDQVIVQRSLSAKSLSHAKGGSVLLGGYLK
 ILPMFFIVMPGMISRALFPDEVGCVDPVCQRICGARVGCNSIAYPKLVMALMPVGLRGLMIAVIMAALM
 SSSLTSIFNSSSTLFTIDVWQRFRRKSTEQELMVVGRVVFVFLVVISILWIPPIQSSNSGQLFDYIQAVTS
 YLAPPITALFLLAIFCKRVTEPGAFWLFGVGLLRMILEFSYPAPACGEVDRRPAVLKDFHYLYFAI
 LLCGLTAIVIVIVSLCTTPIPEEQLTRLTWTRNCPLSELEKEAHESTPEISERPAGECPAGGAAENSS
 LGQEPEAPSRSWGKLLWSWFCGLSGTPEQALSPA EKAAL EQKLT SIEEEPLWRHVCNINAVLLLA INIF
 LWGYFA

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001135181

ORF Size: 2118 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135181.2](#)

RefSeq Size: 3293 bp

RefSeq ORF: 2121 bp

Locus ID: 200010

UniProt ID: [Q2M3M2](#)

Cytogenetics: 1p33

Protein Families: Transmembrane

Gene Summary: Involved in sodium-dependent transport of D-mannose, D-glucose and D-fructose.
[UniProtKB/Swiss-Prot Function]