

Product datasheet for **RG226169**

PDE10A (NM_001130690) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDE10A (NM_001130690) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PDE10A
Synonyms:	ADSD2; HSPDE10A; IOLOD; LINC00473; PDE10A19
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG226169 representing NM_001130690
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAAGATGGACCTTCTAATAATGCGAGCTGCTTCCGAAGGCTGACCGAGTGCTTCTGAGCCCCAGTT
 TGACAGATGAAAAAGTGAAGGCATATCTTCTCTTCAACCCAGGTATTAGATGAATTTGTATCTGAAAG
 TGTTAGTGCAGAGACAGTAGAGAAATGGCTGAAGAGGAAGAACAACAATCAGAAGATGAATCAGCTCCT
 AAGGAAGTCAAGCAGTACCAAGATACGAATATGCAGGGAGTTGTATATGAACTAAACAGCTATATAGAAC
 AACGGTTGGACACAGGAGGAGACAACCAGCTACTCCTCTATGAACTGAGCAGCATCATAAAATAGCCAC
 AAAAGCCGATGGATTTGCACTGTATTTCTTGGAGAGTGAATAATAGCCTGTGTATATTCACGCCACCT
 GGGATAAAGGAAGGAAAACCCCGCTCATCCCTGCTGGGCCATCACTCAGGGCACCACCGTCTCTGCTT
 ATGTGGCCAAGTCCAGGAAAACACTGCTAGTAGAAGACATCCTTGGAGATGAACGATTTCCAAGAGGTAC
 TGGACTGGAATCAGGGACTCGTATCCAGTCTGTTCTTTGCTTACCAATTGCACTGCAATGGTGACTTG
 ATTGGTATTCTCGAGCTGTATCGGCACTGGGGCAAAGAAGCCTTCTGTCTTAGTACCAGGAGGTTGCAA
 CAGCAAATCTTGCTGGCTTCCAGTAGCAATACATCAGGTGCAGGTATGCAGAGGCTTGCCAAACAGAC
 AGAATTGAATGACTTCTACTCGACGTATCAAAAACATATTTTGATAACATAGTTGCAATAGATTCTCTA
 CTTGAACACATAATGATATATGCAAAAAACCTGGTGAATGCCGATCGTTGTGCGCTTTTCCAGGTGGACC
 ATAAGAACAAGGAGTTATATTCAGACCTTTTGTATTGGAGAGGAAAAGGAAGGAAAACCTGTCTTCAA
 GAAGACCAAGAGATAAGATTTCAATTGAGAAAGGAATTGCTGGCCAAGTAGCAAGAACAGGGGAAATC
 CTGAACATTCAGATGCCTATGCAGACCCACGCTTAAACAGAGAAGTAGACTTGTACACAGGCTACACCA
 CGCGGAACATCCTGTGCATGCCATCGTCAGCCGAGGCAGCGTGATAGGTGTGGTGCAGATGGTCAACAA
 AATCAGTGGCAGTGCCTTCTCTAAAACAGATGAAAACAACCTTCAAATGTTTGGCGCTTTTGTGCTTTA
 GCCTTACACTGTGCTAATATGTATCATAGAATTCGCCACTCAGAGTGCATTTACCGGTAACGATGGAAA
 AGCTGTCTACCATAGCATTTGTACTTCAGAAGAGTGGCAAGGTCTCATGCAATTCACCCTTCCCGTGCG
 TCTCTGCAAGAAATGAATTATCCACTTTGACATTGGTCTTTTGAACATGTGGCCTGGAATTTTT
 GTCTACATGGTTCATCGGTCCTGTGGGACATCCTGCTTTGAGCTTGAAGTTGTGTGTTTTATTATGT
 CTGTGAAGAAGAACTATCGGCGGTTCTTATCACAACGGAAGCATGCGGTCACTGTAGCACACTGCAT
 GTATGCCATACTCAGAACAATCACACGCTTTTACAGACCTTGAGCGCAAAGGACTGCTGATTGCGTGT
 CTGTGTATGACCTGGACCACAGGGCTTCAGTAACAGCTACCTGCAGAAGTTGACCACCCTTGCCCG
 CTCTACTCCACTTCCACCATGGAGCAGCACCCTTCTCCAGACTGTGTCCATCCTCCAGTTGGAAGG
 GCACAATATCTTCTCACTCTGAGCTCCAGTGAATATGAGCAGGTGCTTGAGATCATCCGAAAGCCATC
 ATTGCCACAGACCTTGCTTTATACTTTGGAAACAGGAAGCAGTTGGAAGAGATGTACCAGACCGGATCAC
 TAAACCTTAATAATCAATCACATAGAGACCGTGAATTTGGTTTGTATGACTGCCTGTGACCTTTGTTC
 TGTGACAAAACCTGTGGCCGTTACAAAATTGACGGCAAATGATATATATGCAGAATCTGGGCTGAGGGT
 GATGAAATGAAGAAATTTGGGAATACAGCCTATTCCTATGATGGACAGAGACAAGAAGGATGAAGTCCCC
 AAGGCCAGCTTGGGTTCTACAATGCCGTGGCCATTCCTGCTATACAACCTTACCCAGATCCTCCCTCC
 CACGGAGCCTCTTCTGAAAGCATGCAGGGATAATCTCAGTCAGTGGGAGAAGGTGATTCGAGGGGAGGAG
 ACTGCAACCTGGATTTTCATCCCATCCGTGGCTCAGAAGGCAGCTGCATCTGAAGAT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG226169 representing NM_001130690
 Red=Cloning site Green=Tags(s)


```

MEDGPSNNASCFRRLTECF LSPSLTDEKVKAYLSLHPQVLDEFVSESVSAETVEKWLKRKNNKSEDESAP
KEVSRYQDTNMQGVVYELNSYIEQRLDTGGDNQLLLYELSSIIKIATKADGFALYFLGECNNSLCIFTPP
GIKEGKPR LIPAGPITQGTTVSAYVAKSRKTL LVEDILGDERFPRGTGLESGTRI QSVLCLPIVTAIGDL
IGILELYRHWGKEAFCLSHQEVATANLAWASVAIHQVQVCRGLAKQTELNDFLLDVSKTYFDNIVAIDSL
LEHIMIYAKNLVNADRCALFQVDHKNKELYSDLFDIGEEKEGKPVFKKTK EIRFSIEKGIAGQVARTGEV
LNI PDAYADPRFNREVDLYTGYTTRNILCMPIVSRG SVIGVVQMVNKISGSAFSKTDENNFKMF AVFCAL
ALHCANMYHRIRHSECIYRVTMEKLSYHSICTSEEWQGLMQFTLPVRLCKEIELFHF DIGPFENMWPGIF
VYMVHRSCGTSCFELEKLCRFIMSVKKNYRRVPYHNWKHAVTVAHCMYAILQNNHTLFTDLERKGLLIAC
LCHDLDRHGFNSY LQKFDHPLAALYSTSTMEQHHSQT VSVILQLEGHNIFSTLSSEYEQVLEIIRKAI
IATDLALYFGNRKQLEEMYQTGSLNLNNQSHRDRVIGLMMTACDLC SVTKLWPVTKLTANDIYAEFWAEG
DEMKKLGIQPIPMMDRDKKDEVPQQLGFYNAVAIPCYTTLTQILPPTPELLKACRDNL SQWEKVI R GEE
TATWISSPSVAQKAAASED
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_001130690

ORF Size: 2367 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001130690.3
RefSeq Size:	8207 bp
RefSeq ORF:	2370 bp
Locus ID:	10846
UniProt ID:	Q9Y233
Cytogenetics:	6q27
Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism
Gene Summary:	<p>The protein encoded by this gene belongs to the cyclic nucleotide phosphodiesterase family. It plays a role in signal transduction by regulating the intracellular concentration of cyclic nucleotides. This protein can hydrolyze both cAMP and cGMP to the corresponding nucleoside 5' monophosphate, but has higher affinity for cAMP, and is more efficient with cAMP as substrate. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Dec 2011]</p>

Product images:



Circular map for RG226169