

Product datasheet for **RG226294**

ORP150 (HYOU1) (NM_001130991) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ORP150 (HYOU1) (NM_001130991) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	HYOU1
Synonyms:	GRP-170; Grp170; HSP12A; IMD59; ORP-150; ORP150
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)


ORF Nucleotide Sequence:

>RG226294 representing NM_001130991
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAGACAAAGTTAGGAGGCAGAGGCCGAGGAGCGAGTCTGTTGGCCTTGGTGGCTGTGCTCTTGG
 CAGACCTGTTGGCACTGAGTGATACACTGGCAGTGATGTCTGTGGACCTGGGCAGTGAGTCCATGAAGGT
 GGCCATTGTCAAACCTGGAGTGCCCATGAAAATTGCTTGAATAAGGAATCTCGGAGGAAAAACCCGGTG
 ATCGTGACCCTGAAAGAAAATGAAAGATTCTTTGGAGACAGTGCAGCAAGCATGGCGATTAAGAATCCAA
 AGGCTACGCTACGTTACTTCCAGCACCTCTGGGGAAGCAGGCAGATAACCCCATGTAGCTCTTTACCA
 GGCCCGCTTCCCGGAGCACGAGCTGACTTTCGACCACAGAGGCAGACTGTGCACCTTTCAGATCAGCTCG
 CAGCTGCAGTTCTCACCTGAGGAAGTGTGGGCATGGTCTCAATTATTCTGTTCTCTAGCTGAAGATT
 TTGCAGAGCAGCCATCAAGGATGCAGTGATCACCGTCCAGTCTTCTTCAACCAGGCCGAGCGCCGAGC
 TGTGCTGCAGGCTGCTCGTATGGCTGGCCTCAAAGTGCTGCAGCTCATCAATGACAACACCGCCACTGCC
 CTCAGCTATGGTGTCTTCCGCCGAAAGATATTAAACCACTGCCAGAATATCATGTTCTATGACATGG
 GCTCAGGCAGCACCGTATGCACCATTGTGACCTACCAGATGGTGAAGACTAAGGAAGCTGGGATGCAGCC
 ACAGCTGCAGATCCGGGGAGTAGGATTTGACCGTACCCTGGGGGGCTGGAGATGGAGCTCCGGCTTCGA
 GAACGCCTGGCTGGGCTTTTCAATGAGCAGCGCAAGGGTCAGAGAGCAAAGGATGTGCGGGAGAACCCTGC
 GTGCCATGGCCAAGCTGCTGCGTGAGGCTAATCGGCTCAAACCGTCTCAGTGCCAACCGTGACCACAT
 GGCACAGATTGAAGCCCTGATGGATGATGTGGACTTCAAGGCAAAAGTGAAGTGTGGAATTTGAGGAG
 TTGTGTGCAGACTTGTGGAGCGGTGCCTGGGCTGTACAGCAGGCCCTCCAGAGTCCGCAAGTGTGCT
 TGGATGAGATTGAGCAGGTGACTGCTGGTGGTGGGGCCACTCGGTCCCAGAGTTCAGAGGTGCTGTG
 GAAGGCCGTGGCAAGGAGGAGCTGGGGAAGAATCAATGCAGATGAAGCAGCCCATGGGGCAGTG
 TACCAGGCAGCTGCGCTCAGCAAAGCCTTTAAAGTGAAGCCATTTGTCGTCGAGATGCAGTGGTCTACC
 CCATCCTGGTGGAGTTCACGAGGGAGGTGGAGGAGGAGCCTGGGATTCACAGCCTGAAGCACAATAAACG
 GGTACTCTTCTCGGATGGGGCCCTACCCTCAACGCAAAGTATCACCTTTAACCGCTACAGCCATGAT
 TTCAACTTCCACATCACTACGGCGACCTGGGCTTCTGGGGCCTGAAGATCTTCGGGTATTTGGCTCCC
 AGAATCTGACCACAGTGAAGCTAAAAGGGTGGGTGACAGCTTCAAGAAGTATCTGACTACGAGTCCAA
 GGGCATCAAGGCTCACTTCAACCTGGATGAGAGTGGCGTGTGCTCAGTCTAGACAGGTTGGAGTCTGATTT
 GAGACTGTTGAGGACAGCGCAGAAGAGGAATCTACTCTACCAAACCTGGCAACACCATTTCCAGCC
 TGTTTGGAGGCGGTACCACACCAGATGCCAAGGAGAATGGTACTGATACTGTCCAGGAGGAAGAGGAGAG
 CCCTGCAGAGGGGAGCAAGGACGAGCCTGGGGAGCAGGTGGAGCTCAAGGAGGAAGCTGAGGCCCAAGT
 GAGGATGGCTCTCAGCCCCACCCCTGAACCTAAGGGAGATGCAACCCCTGAGGGAGAAAAGGCCACAG
 AAAAGAAAATGGGGACAAGTCTGAGGCCCAGAAACCAAGTGAGAAGGCAGAGGCAGGGCCTGAGGGCGT
 CGCTCCAGCCCCAGAGGGAGAGAAGAAGCAGAAGCCCGCCAGGAAGCGCGAATGGTAGAGGAGATCGGG
 GTGGAGCTGGTTGTTCTGGACCTGCCTGACTTGCCAGAGGATAAGCTGGCTCAGTCGGTGCAGAACTTC
 AGGACTTGACACTCCGAGACCTGGAGAAGCAGGAACGGGAAAAGCTGCCAACAGCTTGAAGCATTTCAT
 ATTTGAGACCCAGGACAAGCTGTACCAGCCGAGTACCAGGAAGTGTCCACAGAGGAGCAGCGTGAGGAG
 ATCTCTGGGAAGCTCAGCGCCGATCCACCTGGCTGGAGGATGAGGGTGTGGAGCCACACAGTGTGAT
 TGAAGGAGAAGCTGGCTGAGCTGAGGAAGCTGTGCCAAGGGCTGTTTTTTCGGGTAGAGGAGCGCAAGAA
 GTGGCCGAACGGCTGTCTGCCCTCGATAATCTCCTCAACCATTCCAGCATGTTCTCAAGGGGGCCCGG
 CTCATCCCAGAGATGGACCAGATCTTCACTGAGGTGGAGATGACAACGTTAGAGAAAAGTCAATGAGA
 CCTGGGCTGGAAGAATGCAACTCTGGCCGAGCAGGCTAAGCTGCCCGCCACAGAGAAGCCTGTGTTGCT
 CTCAAAGACATTGAAGCTAAGATGATGGCCCTGGACCGAGAGGTGCAGTATCTGCTCAATAAGGCCAAG
 TTTACCAAGCCCCGGCCCGCCTAAGGACAAGAATGGGACCCGGGCAGAGCCACCCCTCAATGCCAGTG
 CCAGTGACCAGGGGGAGAAGGTATCCCTCCAGCAGGCCAGACTGAAGATGCAGAGCCATTTCCAGAAC
 TGAGAAAAGTAGAGACTGGATCCGAGCCAGGAGACTGAGCCTTTGGAGTTAGGAGTCTGGAGCAGAA
 CCTGAACAGAAAACAATCGACAGGACAGAAGCGGCCTTTGAAGAACGACGAACTA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Plasmid Map:


ACCN: NM_001130991

ORF Size: 2997 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130991.3](#)

RefSeq Size: 4498 bp

RefSeq ORF: 3000 bp

Locus ID: 10525

UniProt ID: [Q9Y4L1](#)

Cytogenetics: 11q23.3

Protein Families: Transmembrane

Gene Summary: The protein encoded by this gene belongs to the heat shock protein 70 family. This gene uses alternative transcription start sites. A cis-acting segment found in the 5' UTR is involved in stress-dependent induction, resulting in the accumulation of this protein in the endoplasmic reticulum (ER) under hypoxic conditions. The protein encoded by this gene is thought to play an important role in protein folding and secretion in the ER. Since suppression of the protein is associated with accelerated apoptosis, it is also suggested to have an important cytoprotective role in hypoxia-induced cellular perturbation. This protein has been shown to be up-regulated in tumors, especially in breast tumors, and thus it is associated with tumor invasiveness. This gene also has an alternative translation initiation site, resulting in a protein that lacks the N-terminal signal peptide. This signal peptide-lacking protein, which is only 3 amino acids shorter than the mature protein in the ER, is thought to have a housekeeping function in the cytosol. In rat, this protein localizes to both the ER by a carboxy-terminal peptide sequence and to mitochondria by an amino-terminal targeting signal. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014]