

Product datasheet for **RG226532**

CXCR3 (NM_001142797) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CXCR3 (NM_001142797) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CXCR3
Synonyms:	CD182; CD183; CKR-L2; CMKAR3; GPR9; IP10-R; Mig-R; MigR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG226532 representing NM_001142797 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**

ATGGAGTTGAGGAAGTACGGCCCTGGAAGACTGGCGGGACAGTTATAGGAGGAGCTGCTCAGAGTAAATCACAGACTAAATCAGACTCAATCACAAAAGAGTTCCTGCCAGGCCTTTACACAGCCCTTCTCCCCGTTCCC GCCCTCACAGGTGAGTGACCACCAAGTCTAAATGACGCCGAGGTTGCCGCCCTCTGGAGAACTTCAGCTTCTCTATGACTATGGAGAAAACGAGAGTACTCGTGTACCTCCCCGCCCTGCCACAGGACTTCAGCCTGAACCTCGACCGGCCTTCTGCCAGCCCTCTACAGCCTCCTTTTCTGCTGGGGCTGCTGGGCAACGGCGCGGTGGCAGCCGTGCTGCTGAGCCGGCGACAGCCCTGAGCAGCACCACACCTTCTGCTCACCTAGCTGTAGCAGACACGCTGCTGGTGTGACTGCGCTCTGGGAGTGGAGCTGCCGTCCAGTGGGTCTTTGGCTCTGGCCTCTGCAAAGTGGCAGGTGCCCTCTTCAACATCAACTTCTACGCAGGAGCCCTCTGCTGGCCTGCATCAGCTTTGACCGCTACCTGAACATAGTTTCATGCCACCCAGCTTACCGCCGGGGCCCCGGCCCGGTGACCCTCACCTGCCTGGCTGTCTGGGGGCTGCGCTGCTTTTCCGCTCCCAGACTTCATCTTCTGTGCGCCACCACGACGAGCGCCTCAACGCCACCCACTGCCAATACAACTCCCACAGGTGGGCCGCACGGCTCTGCGGGTGTGACGCTGGTGGCTGGCTTTCTGCTGCCCTGCTGGTCATGGCCTACTGCTATGCCACATCCTGGCCGTGCTGCTGGTTCCAGGGGCCAGCGGCCCTGCCGGCCATGCGGCTGGTGGTGGTCTGGTGGCCTTTGCCCTCTGCTGGACCCCTATACCTGGTGGTGGTGGACATCCTCATGGACCTGGGCGCTTTGGCCGCAACTGTGGCCGAGAAAGCAGGGTAGACGTGGCCAAGTCGGTCACCTCAGGCCCTGGGCTACATGCACTGCTGCCTCAACCCGCTGCTCTATGCCTTTGTAGGGTCAAGTCCGGGAGCGGATGTGGATGCTGCTTTGCGCCTGGGCTGCCCAACCAGAGAGGGCTCCAGAGGCAGCCATCGCTTCCCCGCCGGATTATCCTGGTCTGAGACCTCAGAGGCCTCTACTCGGGCTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >RG226532 representing NM_001142797
Red=Cloning site Green=Tags(s)

MELRKYGPRLAGTVIGGAAQSKSQTKSDSITKEFLPGLYTAPSSPFPPSQVSDHQVLNDAEVAALLENF
 SSSYDYGENESDSCCTSPPCPQDFSLNFDRAFLPALYSLLFLLGLLGNGAVAALLSRRTALSSTDTFLL
 HLAVADTLLVLTLPWAVDAAVQWVFGSLCKVAGALFNINFYAGALLACISFDRLNIVHATQLYRRG
 PPARVTLTCLAVWGLCLLFFALPDFIFLSAHHDERLNATHCQYNFPQVGRALRVLQLVAGFLLPLLVMAY
 CYAHILAVLLVSRGQRRLRAMRLVVVVVAFALCWTPYHLVVLVDILMDLGALARNCGRESRVDVAKSVT
 SGLGYMHCCLNPLLYAFVGVKFRERMMWMLLLRLGCPNQRGLQRQPSSRRDSSWSETSEASYSGL

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja3357_a03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001142797


ORF Size: 1245 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001142797.2
RefSeq Size:	1914 bp
RefSeq ORF:	1248 bp
Locus ID:	2833
UniProt ID:	P49682
Cytogenetics:	Xq13.1
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Chemokine signaling pathway, Cytokine-cytokine receptor interaction
Gene Summary:	<p>This gene encodes a G protein-coupled receptor with selectivity for three chemokines, termed CXCL9/Mig (monokine induced by interferon-g), CXCL10/IP10 (interferon-g-inducible 10 kDa protein) and CXCL11/I-TAC (interferon-inducible T cell a-chemoattractant). Binding of chemokines to this protein induces cellular responses that are involved in leukocyte traffic, most notably integrin activation, cytoskeletal changes and chemotactic migration. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. One of the isoforms (CXCR3-B) shows high affinity binding to chemokine, CXCL4/PF4 (PMID:12782716). [provided by RefSeq, Jun 2011]</p>

Product images:



Circular map for RG226532