

## Product datasheet for **RG226552**

### **APLP2 (NM\_001142276) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	APLP2 (NM_001142276) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	APLP2
Synonyms:	APLP-2; APPH; APPL2; CDEBP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RG226552 representing NM\_001142276  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGCGGCCACCGGGACCGCGGCCCGCGACGCCACGGGCAGGCTCCTGCTTCTGCTGCTGGTGGGGCTCA  
CGGCGCCTGCCTTGGCGCTGGCCGGCTACATCGAGGCTCTTGCAGCCAATGCCGGAACAGGATTTGCTGT  
TGCTGAGCCTCAAATCGCAATGTTTTGTGGGAAGTTAAATATGCATGTGAACATTCAGACTGGGAAATGG  
GAACCTGATCCAACAGGCACCAAGAGCTGCTTTGAAACAAAAGAAGAAGTTCTTCAGTACTGTCAGGAGA  
TGTATCCAGAGCTACAGATCACAAATGTGATGGAGGCAAACCAGCGGGTATGATTGACAACCTGGTGCCG  
GAGGGACAAAAGCAATGCAAGAGTCGCTTTGTTACACCTTTCAAGTGTCTCGTGGGTGAATTTGTAAGT  
GATGTCCTGCTAGTTCAGAAAAGTGCCAGTTTTTCCACAAGAGCGGATGGAGGTGTGTGAGAATCACC  
AGCACTGGCACACGGTAGTCAAAGAGGCATGTCTGACTCAGGGAATGACCTTATATAGCTACGGCATGCT  
GCTCCCATGTGGGGTAGACCAGTTCATGGCACTGAATATGTGTGCTGCCCTCAGACAAAGATTATTGGA  
TCTGTGTCAAAGAAGAGGAAGAGGAAGATGAAGAGGAAGAGGAAGGAAGATGAAGAGGAAGACTATG  
ATGTTTATAAAAGTGAATTTCTACTGAAGCAGATCTGGAAGACTTCACAGAAGCAGCTGTGGATGAGGA  
TGATGAGGATGAGGAAGAAGGGGAGGAAGTGGTGGAGGACCGAGATTACTACTATGACACCTTCAAAGGA  
GATGACTACAATGAGGAGAATCCTACTGAACCCGGCAGCGACGGCACCATGTCAGACAAGGAAATTA  
ATGATGTCAAAGCTGTCTGCTCCAGGAGGCGATGACGGGGCCCTGCCGGGCCGTGATGCCTCGTTGGTA  
CTTCGACCTCTCCAAGGAAAGTGCCTGCGCTTATATATGGTGGCTGCGGGCGCAACAGGAACAATTTT  
GAGTCTGAGGATTATTGTATGGCTGTGTGTAAGCGATGATTCTCCAACCTCCTCGCAACCAATGATG  
TTGATGTGATTTTCGAGACCTCTGCAGATGATAATGAGCATGCTCGCTTCCAGAAGGCTAAGGAGCAGCT  
GGAGATTCGGCACCGCAACCGAATGGACAGGGTAAAGAAGGAATGGGAAGAGGCAGAGCTTCAAGCTAAG  
AACCTCCCAAAGCAGAGAGGCAGACTCTGATTAGCACTTCCAAGCCATGGTTAAAGCTTTAGAGAAGG  
AAGCAGCCAGTGAGAAGCAGCAGCTGGTGGAGACCCACCTGGCCCGAGTGAAGCTATGCTGAATGACCG  
CCGTGGATGGCTCTGGAGAACTACCTGGCTGCCTTGCAGTCTGACCCGCCACGGCCTCATCGCATTCTC  
CAGGCCTTACGGCGTTATGTCCGTGCTGAGAACAAGATCGCTTACATACCATCCGTATTACCAGCATG  
TGTTGGCTGTTGACCCAGAAAAGGCGGCCAGATGAAATCCCAGGTGATGACACATCTCCACGTGATTGA  
AGAAAGGAGGAACCAAAGCCTCTCTGCTCTACAAAGTACCTTATGTAGCCCAAGAAATCAAGAGGAA  
ATTGATGAGCTCCTTCAGGAGCAGCGTGCAGATATGGACCAGTTCCTGCCTCAATCTCAGAGACCCCTG  
TGGACGTCCGGGTGAGCTCTGAGGAGAGTGAGGAGATCCACCGTTCACCCCTTCCACCCCTTCCAGC  
CCTACCTGAGAACGAAGGATCTGGAGTGGGAGAGCAGGATGGGGGACTGATCGGTGCCGAAGAGAAAGT  
ATTAACAGTAAGAATAAAGTGGATGAAAACATGGTCATTGACGAGACTCTGGATGTTAAGGAAATGATTT  
TCAATGCCGAGAGAGTTGGAGGCCTCGAGGAAGAGCGGGAATCCGTGGGCCCACTGCGGGAGGACTTCAG  
TCTGAGTAGCAGTGTCTCATTGGCCTGCTGGTCATCGCAGTGGCCATTGCCACGGTTCATCGTCATCAGC  
CTGGTGTGCTGAGGAAGAGGCAGTATGGCACCATCAGCCACGGGATCGTGGAGGTTGATCCAATGCTCA  
CCCCAGAAGAGCGTCACCTGAACAAGATGCAGAACCATGGCTATGAGAACCCACCTACAATACCTGGA  
GCAGATGCAGATT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG226552 representing NM\_001142276  
Red=Cloning site Green=Tags(s)

```

MAATGTAAAAATGRLLLLLLVLGTAPALALAGYIEALAANAGTGFVAEAPQIAMFCGKLNMHVNIQTGKW
EPDPTGKTKSCFETKEEVLQYQCQEMYPELQITNVMEANQRVSIIDNWCRRDKKQCKSRFVTPFKCLVGEFVS
DVLLVPEKQCFHKERMEVCENHQHWHVTVVKEACLQGMTLYSYGMLLPCGVDQFHGTEYVCCPQTKIIG
SVSKEEEEEDEEEEEDEEEDYDVKSEFPTEADLEDFTEAAVDEDEDEEEEGEEVVEDRDYYYDTFKG
DDYNEENPTEPGSDGTMSDKEITHDVKAVCSQEAMTGPCRAVMPRWYFDLSK GKCVRFIYGCGGNRNNF
ESEDYCMVCKAMIPPTPLPTNDVDVYFETSADDNEHARFQKAKEQLEIRHRNRMDRVKKEWEEAELQAK
NLPKAERQTLIQHFQAMVKALEKAASEKQQLVETHLARVEAMLNDRRRMALENYLAALQSDPPRPHRIL
QALRRYVRAENKDRHTIRHYQHVLAVDPEKAAQMKSQVMTHLHVIEERRNQLSLLYKVPYVAQEIQEE
IDELLQEQRADMDQFTASISETPVDVRSSESESEIIPPFHPFPFPALPENEGSGVGEQDGLIGAEKVK
INSKNKVDENMVIDETLDVKEMIFNAERVGGLEEEERESVGPLREDFSLSSSALIGLLVIAVAIATVIVIS
LVMLRKRQYGTISHGIVEVDPMLTPEERHLNKMQNHYENPTYKYLEQMQUI
    
```

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



## Plasmid Map:



ACCN: NM\_001142276

ORF Size: 2253 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001142276.2](#)

**RefSeq Size:** 3791 bp

**RefSeq ORF:** 2256 bp

**Locus ID:** 334

**UniProt ID:** [Q06481](#)

**Cytogenetics:** 11q24.3

**Protein Families:** Druggable Genome, Transmembrane

**Gene Summary:** This gene encodes amyloid precursor- like protein 2 (APLP2), which is a member of the APP (amyloid precursor protein) family including APP, APLP1 and APLP2. This protein is ubiquitously expressed. It contains heparin-, copper- and zinc- binding domains at the N-terminus, BPTI/Kunitz inhibitor and E2 domains in the middle region, and transmembrane and intracellular domains at the C-terminus. This protein interacts with major histocompatibility complex (MHC) class I molecules. The synergy of this protein and the APP is required to mediate neuromuscular transmission, spatial learning and synaptic plasticity. This protein has been implicated in the pathogenesis of Alzheimer's disease. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2011]