

Product datasheet for **RG226585**

C22orf46 (NM_001142964) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	C22orf46
Synonyms:	FLJ23584
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG226585 representing NM_001142964
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTTGTTGTCACTGCTGGGTGCCTGTGCTGTGGTGGGGCCATTCCATGGCCCTGAGTGGGAGCCAGTGC
AGGGCCTGCTCTCCAGAATCACAGCTGCAGGGACCCCTCAGTGCTGTGGCAACCTGCTTGTCTCTGCCT
CTTTCTGGTCTGGCAGGTCCGGCACTGTTGGCACCAGGTACCAGGACCCGCTTCAGCACAAGGAATGTC
ATCAAGGTGCCACTGCAGAAGCGGGCAGTGCCCTCCATGAGGTGCGAACTGTCTTCAAGCTGACTCCTG
AATTCTTCAGTCTGGAAGTCCAGGGGCTAGATTCTCAACAATGCGCACAAGGCAGAGATGGGGATA
CCGGAGGAGCCTCCAGGAATCATGGGCCAGAACCTGCTCTCTCCACAGCACCCATGTCCAGGCCACCT
TCGGGTGTCCACCCACTCTGAGCCCATCTTTGTACCACCTCCATTTCAAACACCTGTTTACTGCCTC
AGAACAGTTCTGGAAGCATGGCAGGTGCCTTGGTGTCTCCATGATGGTCAGACTCGCCCTGCCTTGA
CATGTGTCAAGAGATGGAGCAGCTGCTGCTTCACTCACAGGAAAGGTTGGTGTCACTGGAGCCTGTCATC
AGTGTGAGGTCTCGCCCCACCTCCATGACCTTAACCACCTCTCTTCCAACTTACTTTCAGCTGAGAGGC
TGCAGTTCTGCCCCAGAGAGCTCCTGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG226585 representing NM_001142964
Red=Cloning site Green=Tags(s)

MLLSLLGACAVVGPFGPEWEPVQGLLSQNHSCRDPQCCGNLLVLCLFLVWQVRHCWHQVTRTRFSTRNV
IKVPLQKRAVPSMRCETVFKLTPEFFSPGKSRGLDSQQCAQRQWGYRRSLQESWAQNLLSPQHPCPGPP
SGVHTHSEPIFCTTISINTCLLPQNSSWKAWQVPWCLHDGQTRPALDMCQEMEQLLLHSQERLVSLEPVI
SVRSRPTSMTLTSLPNLLSAERLQFCPQAPAP

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001142964

ORF Size: 729 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001142964.1, NP_001136436.1


RefSeq Size: 4991 bp

RefSeq ORF: 732 bp

Locus ID: 79640

Cytogenetics: 22q13.2

Product images:



Circular map for RG226585