

Product datasheet for **RG226660**

ALS2 (NM_001135745) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALS2 (NM_001135745) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ALS2
Synonyms:	ALS2CR6; ALSJ; IAHSP; PLSJ
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG226660 representing NM_001135745 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACTCAAAGAAGAGAAGCTCAACAGAGGCAGAAGGATCCAAGGAAAGAGGCCTGGTCCATATCTGGC
AGGCAGGATCCTTTCCATAACACCAGAGAGATTGCCAGGCTGGGAGGAAAGACTGTTTTGCAGGCAGC
CCTCGGAGTGAACATGGAGTTCTTCTGACTGAAGATGGTGAGGTCTACAGCTTTGGGACTTCCCTGG
AGAAGTGGACCAGTGGAGATTTGTCCAAGTAGCCCATCTAGAAAATGCCCTGGTTGGCAATATGTTA
TTACTGTGGCAACAGGAAGCTTCCATAGTGGAGCAGTGACAGACAATGGTGTGCGGTACATGTGGGAGA
GAATTCGCTGGCCAGTGTGCAGTAGCCAACCAGCAGTATGTGCCGGAACCAAATCCTGTGAGCATTGCT
GATTCTGAGGCCAGCCCTTTGTTAGCAGTCAGGATTTTACAGTTGGCGTGTGGCGAGGAGCACACTCTGG
CATTGTCAATAAGCAGAGAGATTTGGGCATGGGGTACCGGTTGTGAGTTGGGTCTCATTACCACTGCCTT
CCCAGTGACAAAGCCGAAAAGGTAGAACATCTTGCTGGGCGAGTGGTGTCAAGTTGCCTGTGGTGTCT
TTCCACAGCTTAGCCCTGTACAATGCCTCCCTCCCAGGATCTGAAGCCAGTCCCAGAACGATGCAACC
AGTGCAGCCAGCTCTTGATTACTATGACTGACAAAAGAAGACCATGTGATTATATCAGACAGTCATTGTTG
CCCATTAGGTGTGACACTGACAGAATCTCAGGCAGAAAACCATGCCAGCACTGCTCTCAGCCCCTCCACT
GAAACCCCTTGACAGGCAGGAAGAAGTATTTGAGAACACTCTGTAGCAAATGATCAGTCTGTTGCTACTG
AACTGAATGCAGTAAGTGTCTCAGATCACAAAGCAGCGATGCCATGTCCTCTCAACAAAATGTCCATGGGAA
AACTGAAATTTCTCTGCCAGAAACATACCATCATACCCTGACACCCAAGCAGTCAATGAATACCTACGG
AAACTGTCAGATCATTAGTAAGAGAGGACTCAGAGCATGGTGAAAAGCCAGTGCATCTCAGGTACCTG
CTCAATTTTATAAAAATAAAGTGTGTCTAGAGTTGAACTGTATGGGCTTTTCTTTAGAACTTTAAAA

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG226660 representing NM_001135745
 Red=Cloning site Green=Tags(s)

MDSKRSSTEAEGSKERGLVHIWQAGSFPITPERLPGWGGKTVLQAALGVKHGVLLEDGEVYSFGTLPW
 RSGPVEICPSSPILENALVGQYVITVATGSFHSGAVTDNGVAYMWGENSAGQCAVANQQYVPEPNPVSIA
 DSEASPLLAVRILQLACGEEHTLALSISREIWAAGTGCCQLGLITTAFPVTKPKQKVEHLAQRVVLQVACGA
 FHSLALVQCLPSQDLKPVPERCNQCSQLLITMTDKEDHVIISDSHCCPLGVTLTESQAENHASTALSPST
 ETLDRQEEVFENTLVANDQSVATELNAVSAQITSSDAMSSQQNVMTTEISSARNIPSPYPTDQAVNEYLR
 KLSDSHSVREDSEHGKPVPSQVPAQFYKIKVCLELNCMGFSLETLK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001135745

ORF Size: 1188 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001135745.2
RefSeq Size:	2908 bp
RefSeq ORF:	1191 bp
Locus ID:	57679
UniProt ID:	Q96Q42
Cytogenetics:	2q33.1
Protein Families:	Druggable Genome
Protein Pathways:	Amyotrophic lateral sclerosis (ALS)
Gene Summary:	The protein encoded by this gene contains an ATS1/RCC1-like domain, a RhoGEF domain, and a vacuolar protein sorting 9 (VPS9) domain, all of which are guanine-nucleotide exchange factors that activate members of the Ras superfamily of GTPases. The protein functions as a guanine nucleotide exchange factor for the small GTPase RAB5. The protein localizes with RAB5 on early endosomal compartments, and functions as a modulator for endosomal dynamics. Mutations in this gene result in several forms of juvenile lateral sclerosis and infantile-onset ascending spastic paralysis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008]

Product images:



Circular map for RG226660