

Product datasheet for **RG226714**

NR2F2 (NM_001145156) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NR2F2 (NM_001145156) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NR2F2
Synonyms:	ARP-1; ARP1; CHTD4; COUPTF2; COUPTFB; COUPTFII; NF-E3; SRXX5; SVP40; TFCOUP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG226714 representing NM_001145156 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAATGGTAGTCAGCACGTGGCGGACCCCCAGGACGAGGTGCCCGGCTCACAGGGCAGCCAGGCCCT
CGCAGGGCGCCCGTGCCCGGCCCGCCCGGCCCGCCCGCACACGCCACAGACGCCCGGCAAGGGG
CCCAGCCAGCACGCCAGCCAGACGGCGCGCGGTGGCCAGGGCGGCCCTGGCGGCCCGGGTAGCGACAAG
CAGCAGCAGCAGCAACACATCGAGTGCCTGGTGTGCGGAGACAAGTCGAGCGCAAGCACTACGGCCAGT
TCACGTGCGAGGGCTGCAAGAGCTTCTTCAAGCGCAGCGTGGGAGGAACCTGAGCTACAGTGGCCGC
CAACCGAACTGTCCATCGACCAGCACCATCGCAACCAGTGCCAGTACTGCCGCCTCAAAAAGTGCCCTC
AAAGTGGGCATGAGACGGGAAGCGGTGCAGAGGGGAGGATGCCGCCGACCCAGCCGACCCACGGGAGT
TCGCGCTGACCAACGGGGATCCCCTCAACTGCCACTCGTACCTGTCCGGATATATTTCCCTGCTGTTGCG
CGCGGAGCCCTATCCCACGTGCGCTTCGGCAGCCAATGCATGCAGCCCAACAACATCATGGGTATCGAG
AACATTTGCGAACTGGCCGCGAGGATGCTCTTCAGCGCCGTCGAGTGGGCCCGGAACATCCCCTTCTCC
CCGACCTGCAGATCACGGACCAGGTGGCCCTGCTTCGCCCTCACCTGGAGCGAGCTGTTTGTGTTGAATGC
GGCGCAGTGCTCCATGCCCTCCACGTGCCCCGCTCCTGGCCCGCCCGGCTGCATGCTTCGCCATG
TCCGCCGACCGGGTGGTGCCTTTATGGACCACATACGGATCTTCAAGAGCAAGTGAGAAAGCTCAAGG
CGCTGCAGTTGACTCAGCCGAGTACAGCTGCCTCAAGGCCATAGTCTGTTACCTCAGATGCCTGTGG
TCTCTCTGATGTAGCCATGTGGAAGCTTGCAGGAAAAGTCTCAGTGTGCTTTGGAAGAATACGTTAGG
AGCCAGTACCCCAACCAGCCGACGAGATTCGGAAGCTTTTGCTTCGCTCCCTCCCTCCGACCGTCT
CCTCCTCAGTCATAGAGCAATTGTTTTTCGTCGTTTGGTAGGTAACCCCATCGAAACCCTCATCCG
GGATATGTTACTGTCGGCAGCAGTTTTAACTGGCCGTATATGGCAATTCAA

ACGGTACGGCGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >RG226714 representing NM_001145156
Red=Cloning site Green=Tags(s)

MAMVVSTWRDPQDEVPGSQGSQASQAPPVPGPPPAPHTPQTGQGGPASTPAQTAAAGQGGPGGPGSDK
 QQQQHHIECVVCGDKSSGKHYGQFTCEGCKSFFKRSVRRNL SYTCRANRNCPIDQHHRNQCYCRLKKCL
 KVGMRREAVQRGRMPPTQPTHGQFAL TNGDPLNCHSYLSGYISLLLRAPYPTSRFGSQCMQPNNIMGIE
 NICELAARMLFSAVEWARNIPFPDLQITDQVALLRLTWESEFVLNAAQCSMPLHVAPLLAAAGLHASP
 SADRVVAFMDHIRIFQEQQVEKLKALHVDSAEYSCLKAIVLFTSDACGLSDVAHVESLQEKSQCALEEYVR
 SQYPNQPTRFGKLLRLPSLRTVSSSVIEQLFFVRLVGKTPIETLIRDMLLSGSSFNWPYMAIQ


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001145156

ORF Size:	1245 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001145156.1 , NP_001138628.1
RefSeq Size:	3501 bp
RefSeq ORF:	786 bp
Locus ID:	7026
UniProt ID:	P24468
Cytogenetics:	15q26.2
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Gene Summary:	This gene encodes a member of the steroid thyroid hormone superfamily of nuclear receptors. The encoded protein is a ligand inducible transcription factor that is involved in the regulation of many different genes. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]