

Product datasheet for **RG226737**

TMEM106C (NM_001143842) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	TMEM106C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG226737 representing NM_001143842
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGGGTCTCAGCATTCCGCTGCTGCTCGCCCTCCTCCTGCAGGCGAAAGCAAGAAGATGACAGGGACG
GTTTGTGGCTGAACGAGAGCAGGAAGAAGCCATTGCTCAGTTCCTCATGTGGAATTCACCGGAGAGA
TAGCATCACCTGTCTCAGTGCCAGGGGACAGGCTACATTCCAACAGAGCAAGTAAATGAGTTGGTGGCT
TTGATCCACACAGTGATCAGAGATTGCGCCCTCAGCGAACTAAGCAATATGTCCTCCTGTCCATCCTGC
TTTGTCTCTGGCATCTGGTTGGTGGTTTTCTCCTGTTTCCGCATTAGTCCTTGTGGATGATGACGG
CATCAAAGTGGTGAAAGTCACATTTAATAAGCAAGACTCCCTTGTAAATCTCACCATCATGGCCACCCTG
AAAATCAGGAACCTCAACTTCTACACGGTGGCAGTGACAGCCTGTCCAGCCAGATTAGTACATGAACA
CAGTGGTCAGTACATATGTGACTACTAACGTCTCCCTTATTCCACCTCGGAGTGAGCAACTGGTGAATTT
TACCGGAAGGCCGAGATGGGAGGACCGTTTTCTATGTGTACTTCTTGCACGGTACCTGAGATCCTG
GTGCACAACATAGTGATCTTCATGCGAACTTCAGTGAAGATTTATACATTGGCCTCATGACCCAGAGCT
CCTTGGAGACATCACTATGTGGATTGTGGAGGAAATTCACAGCTATT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG226737 representing NM_001143842
Red=Cloning site Green=Tags(s)


MGSQHSAAARPSSCRKQEDDRDGLLAEREQEEAIAQFPYVEFTGRDSITCLTCQGTGYIPTQVNELVA
LIPHSDQRLRPQRTKQYVLLSILLCLLASGLVVFLLFPHSVLVDDDGIKVVKVTFNKQDSLVIITIMATL
KIRNSNFYTVAVTSLSSQIQYMNTVVSTVYTTNVSLIPRSEQLVNFTGKAEMGGPFSYVYFFCTVPEIL
VHNIVIFMRTSVKISYIGLMTQSSLETHHYVDCGNGSTAI

TRTRPLE - GFP Tag - V



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001143842

ORF Size: 750 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001143842.2
RefSeq Size: 1606 bp
RefSeq ORF: 753 bp
Locus ID: 79022
UniProt ID: Q9BVX2
Cytogenetics: 12q13.11
Protein Families: Transmembrane

Product images:



Circular map for RG226737