

Product datasheet for **RG226794**

CRISP2 (NM_001142407) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CRISP2 (NM_001142407) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CRISP2
Synonyms:	CRISP-2; CT36; GAPDL5; TPX1; TSP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG226794 representing NM_001142407 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTTTACTACCGGTGTTGTTTCTGGTACTGTGCTGCTTCCATCTTTACCTGCAGAAGGAAAGGATC
CCGCTTTTACTGCTTTGTTAACCACCCAGTTGCAAGTGCAAAGGGAGATTGTAATAAACACAATGAACT
AAGGAAAGCAGTCTCTCCACCTGCCAGTAACATGCTAAAGATGGAATGGAGCAGAGAGGTAACAACGAAT
GCCCAAAGGTGGCAAACAAGTGCCTTTACAACATAGTGATCCAGAGGACCGCAAACCAAGTACAAGT
GTGGTGAAGTCTCTATATGTCAAGTGACCCTACTTCCTGGTCTTCTGCAATCCAAGCTGGTATGACGA
GATCCTAGATTTTGTCTATGGTGTAGGACCAAAGAGTCCCAATGCAGTTGTTGGACATTATACTCAGCTT
GTTTGGTACTCGACTTACCAGGTAGGCTGTGGAATTGCCTACTGTCCCAATCAAGATAGTCTAAAATACT
ACTATGTTTGCCAATATTGTCCTGCTGGTAATAATATGAATAGAAAGAATACCCCGTACCAACAAGGAAC
ACCTTGTGCCGTTGCCCTGATGACTGTGACAAAGGACTATGCACCAATAGTTGCCAGTATCAAGATCTC
CTAAGTAACTGTGATTCTTGAAGAATACAGCTGGCTGTGAACATGAGTTACTCAAGGAAAAGTGAAGG
CTACTTGCCTATGTGAGAACAAAATTTAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG226794 representing NM_001142407
Red=Cloning site Green=Tags(s)

MALLPVLFLVTLLPSLPAEGKDP AFTALLTTQLQVQREIVNKHNELRKAVSPPASNMLKMEWSREVTTN
 AQRWANKCTLQHSDPEDRKTSTRCGENLYMSSDPTSWSSAIQSWYDEILDFVYGVGPKSPNAVVGHYTQL
 VWYSTYQVGGCIAYCPNQDSLKYYYYVCQYCPAGNNMNRKNTPYQQGTPCAGCPDDCDKGLCTNSCQYQDL
 LSNCDLKNNTAGCEHELLKEKCKATCLCENKIY

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001142407

ORF Size: 729 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001142407.3
RefSeq Size:	1365 bp
RefSeq ORF:	732 bp
Locus ID:	7180
UniProt ID:	P16562
Cytogenetics:	6p12.3
Protein Families:	Secreted Protein
Gene Summary:	May regulate some ion channels' activity and thereby regulate calcium fluxes during sperm capacitation.[UniProtKB/Swiss-Prot Function]