

Product datasheet for **RG226878**

CXorf49B (NM_001145139) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	CXorf49B
Synonyms:	CXorf49
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG226878 representing NM_001145139
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAACTACGACTCACTATAGGGCGGCCGGAATTCGTCTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGAGCTCCCCGACAAGGTGTCCGTTTGTGGGCGGTTTCGACCTGGAGGGTGGCAAAAAGGCTGGCT
 CCCGCACGGCCAGCCCCGAGCCCCAGGGGCCACAGCCACGGCCTCGATTGGGGGTGCCGGCAGCGG
 CGATGGCAAGAGCGAGAGCGGGTTACAGATCCAGAGGGCTTCAGCTTCGAGTCTGAGAGCGAATTGATA
 GAGCAAGGAAGGGTGGTCTCTGGGGCCGGGAAGGCCGCCAGGCACCCCGGTGGATGACCAAGGGGACG
 TTGTGGACTACTCATTCTACCTGGCTGACGAACAGCCGCCATCGTGCCGCCGCCAGCGTCCAGGGACA
 CCGCTTCCCAAGAAGGTGCCGTGCCGAAGGGTCCGCTGAGAATTGGGCAGATGCGGAGGTCCGTCCAGT
 GGGAGAGACGTGCTGGGCCACAGCCCTGGAATGGCAGCAGGCCTCTGCCGGCGTCTCCACCTCTGCG
 GTCTGGGCCAGTGCAGGCTGGAAGAACCCGGAAGGGGCTCGAAGAGCAGATGGAGCCTCCGCGTGGA
 TCCCCAGCAGCCCTCTGCGAAGGCCCCACCAGGCTGCCTACCCAGACTCTGATTCCGAGATGAGAGC
 AGCGACTTACCACTGATGAAGGTAGGCATTTCGCCGAACGAAGGAAGCCAGGCCAAGCCCGCAGCCCCA
 AGAAGCGAGCAGACATCCAGACAGGCAAGCTTCACTGCAAGGAGAGTTACCTGCCTGTGCCGGGCCG
 TTTCTGACCTCTGCTCCCGCGGACTCACTCCAGTCGAGAGAGGCGGCTGTGGGAGAGCTGGAGGAC
 TCTCCCAAGAAAATGCAGAGCAGGGCCTGGGAAAGGTGGAGGTGAGGCCAGCTGCTCAGGAGCTG
 CTGCTGACGGGGCCTGCCCGAGGCTTTTCGAGGAGGAAGATGGCCGGGGGGAAGAAGTCCCTAGGGGG
 TGCTCTCACTGGCCCTGGGAGAGGCTTTCTGCTGCGGAGAGAGACTCTCAGCCGCTCCCGGAG
 CCGGCCACCTTCCCGCATTCTCTGGTGTGCGGCCACAGGGGATGTCCAAGAAACCCAAAGCCTAAGC
 ACAGCAGCCCTGGGAAGAAACAGCAGGGAGGAAGACCAGGGAGTCCAGGCTGCGGCCAGAGAAGATAA
 TGACCCAAATAGAGATGAGGTCCCAAGAGCCCACTTCCACCCACAGGCCAGGACTGCCTCGCCTGTCT
 GTGCGTCTGGAGAATTCAGCAGTAGCGACCCCAACATCAGAGCTCCCAACTCCGGGAATTCAGAGC
 CCTCGGCCTACAGCCGGGAGGCTCTGTCAGAGCCATGCACCTCCGTAACAGCAGCCGCTGT
 CCATCCCCAAGACCGAAAGCAGCAGCAGCCCCGGGAGCCAGGGCTGTCTCGGTGCATCTGGCTG
 CAGAGGGAAATTGAGACCTTACACAGCAGCTAGCGGCCATGCAGTTCTCACTGACAAGTTCAGGACC
 TT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG226878 representing NM_001145139
 Red=Cloning site Green=Tags(s)

MSSPDKVSVCGAGFDLEGGKKAGSRTASPGAPGAHSHGLDLGVPGSGDGKSESFTDPEGFSFESESELI
 EQGRVVLWGREGRPVDDQGDVVDYSFYLADEPAAIVPPPSVQGHFPPEGAAAEGSAENWADAIEVGPS
 GRDVLGHSPGKWKQASAGRLHLCGPVRAWKNPERGSKSRWSLRVDPQQPSAKGPTRLPTHDSADAES
 SDLPLMKVIGICRNEGSQAKPGSPKKRADTSRQASFHCKESYLPVPGRFLTSAPRGLTPVAERPVGLEED
 SPQKKMQSRAWGKVEVRPSCSGAAAAGALPQGLSRRKMAGGKSLGGASQLALGRGFPACGERLSAAPPE
 PATFPPFSGVRPQGMKKPKPKHSSPGKKPAGRKTRESQAAARENDPNRDEVPRALPTHRPLRLS
 VRRGEFSSSDPNIRAPQLPGTSEPSAYSPGGLVPRRHAPSGNQPPVHPRPERQQQPPGAQGCPRCIWL
 QREIEDLTQQLAAMQFLTDKFQDL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001145139

ORF Size: 1542 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001145139.1](#), [NP_001138611.1](#)

RefSeq Size: 1851 bp


RefSeq ORF: 1545 bp

Locus ID: 100132994

UniProt ID: [A8MYA2](#)

Cytogenetics: Xq13.1

Product images:



Circular map for RG226878