

Product datasheet for **RG226918**

PASIC1 (LMNTD1) (NM_001145727) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	PASIC1
Synonyms:	IFLTD1; LMNARSI; PASIC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG226918 representing NM_001145727
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGCAACAGCCAGTGAGAGCAGGAGACATCGGTGTTTCCAGAGTCCAGACCCACCACAATGCCAGCCCT
TCTCAGTTCGAAGAAACAGGATGAATCACCATGATTGGGGATGGAGAAGATTATTTCTTTCTTTGTT
TGGTGATTCAAAGAACTTACAGCAGCTCAAACTACACTCAGAAAACTTAAATACTTTCTATGATT
CTTGAAGAAGTTGGCAATTTACCTCCAGTCTCTTGGAGATGTTGAAATAGCTGAAGTGAATGTCAAGG
GTTTGTTTCGTGAAGCTCATTAACTCTCCCTTGACAAAGAAATGGCAATTGGAGATCATATTCTCCAGCA
AAATGTGAATGGACAAACCATTTCTTTGTACCGATTCTTCCAAACATCGTAATGCAGGCAAATTCACA
GTAACAGTGTGGGCAGCAGCATCTGAAGCAAAGCATCAACCTCCATCAGATTTCTTTGGAAGGAACAAG
ACAAGTTTAGAGCAAGTCCTGATTGTATAACAATCCTGTGCAAACCGAACGGTCAAGCCATTGCGTGTA
CACCCCTATCCACTGGAAGCAAGCGTGGGAAAAATTAGATGCTGACGTTGAATTTAACAGATGTTCAAGTA
GTATCTCCAACATTCGAAAGCGTGTGTTTCAGTGGACAGCATCTACAGCTACAATAACTAAAGAAAAAC
AAGATCAACCTAAGAAAGATATCTCAAATTATCAGGTGGAACAAGCTCAAGTTCTTCTTAAGAGAGAGAA
GGAAATCCCACCAACCGTTTTCCCTAATCGCAGCCCTTGGTGCCAGAATCCCTATGTCTCTGCACATCCT
TACTGTCCTCTGATTGAACCACACAATACATCCACCGCTGGAGGCAGATTGGATAGACAGCCAGGACTC
GGTCAACCAGACCTAATCGAGCCTCAGGTCTAAGAAAAAGAAGACATCTGAGTCACAAAAGCAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG226918 representing NM_001145727
 Red=Cloning site Green=Tags(s)

MQQPVRAGDIGVSRVQTHHNAPFSVPKKQDESPMIGDGEDYFLSLFGDSKKLTAHSNYTQKTLKYFSMI
 LEEVGQFTSSSLGDVEIAEVNVKGLFVKLINSSLDKEMAIGDHILQQNVNGQTISLYRFLPNIVMQANST
 VTVWAAASEAKHQPPSDFLWKEQDKFRASPDCTILCKPNGQAIWYTPIHWKQAWEKLDADVEFNRCVS
 VSPTFRKRVPQWTASTATITKEKQDQPKKDISNYQVEQAQVLLKREKEIPPTVFPNRSPPWCQNPVYSAHP
 YCPLIEPHNTSTAGGRLDQRPRTRSTRPNRASGSKKKKTSQSQQ

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



CTATAGGGCGGCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGATCGCCGGCGCCAGATCT

EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI

CAAGCTTAAGTACTAGTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC --- --- ---

HindIII NheI RsrII MluI NotI XhoI GFP Tag

T R T R P L E M E S D - - -

--- --- GAA GAA AGA GTT TAA ACGGCCGCGCGGAGCT

- - E E R V Stop

PmeI FseI

ACCN: NM_001145727

ORF Size: 975 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001145727.3](#)

RefSeq Size: 1285 bp


RefSeq ORF: 978 bp

Locus ID: 160492

UniProt ID: [Q8N9Z9](#)

Cytogenetics: 12p12.1

Product images:



Circular map for RG226918