

Product datasheet for **RG226928**

PASIC1 (LMNTD1) (NM_001145728) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	PASIC1
Synonyms:	IFLTD1; LMNARSI; PASIC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG226928 representing NM_001145728
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGAAAGATACACAAGACATTCAGGAAGCTTCGAAGGCAATGCAGAATAAAGTCCATGAGCAGGAAGATA
AGAATGAGAAACAAAAACAAAGAGAAGACAACTTGGAGTATATCTTTAGTACATTTTCCCAAAGAT
GTTGGGTTCAGTTGCCACAACACTGCCATTGTCAATCTTCAATCCAGTGGATGCCTCTTGGTTACTAT
CTGTCTAGTCCTCAGATTAGTAGAGTAATATCAACAACCTGGACAATTGACTTCTAAAGCTACTGTAG
GTAGCTGTTCCAGAGTGGAAAACAGCTTGGATGCCAGCCCCTTCTCAGTTCGGAAGAAACAGGATGAATC
ACCCATGATTGGGGATGGAGAAGTATTTCTTTCTTTGTTGGTGATTCAAAGAACTTACAGCACAC
TCAAACACTCAGAAAACCTTAAATACTTTCTATGATTCTTGAAGAAGTTGGCCAATTTACCTCCA
GTTCTCTTGGAGATGTTGAAATAGCTGAAGTGAATGTCAAGGGTTTGTTCGTGAAGCTCATTAACCTTC
CCTTGACAAAGAAATGGCAATTGGAGATCATATTCTCCAGCAAAATGTGAATGGACAAACATTTCTTTG
TACCGATTCTTCCAAACATCGTAATGCAGGCAAATTCACAGTAACAGTGTGGGCAGCAGCATCTGAAG
CAAAGCATCAACCTCCATCAGATTTCTTTGGAAGGAACAAGACAAGTTTAGAGCAAGTCTGATTGTAT
AACAACTCTGTGCAAAACCGAAGCTCAAGCCATTGCGTGGTACACCCCTATCCACTGGAAGCAAGCGTGG
GAAAAATTAGATGCTGACGTTGAATTTAACAGATGTTCAAGTAGTATCTCCAACATTCCGAAAGCGTGTGT
TTCAGTGGACAGCATCTACAGCTACAATACTAAAGAAAAACAAGATCAACCTAAGAAAGATATCTCAAA
TTATCAGGTGGAACAAGCTCAAGTTCTTCTTAAGAGAGAGAAGGAAATCCACCAACCGTTTTCCCTAAT
CGCAGCCCTTGGTGCCAGAATCCCTATGTCTCTGCACATCCTTACTGTCTCTGATTGAACACACAATA
CATCCACCGCTGGAGGCAGATTGGATAGACAGCCAGGACTCGGTCAACCAACCTAATCGAGCCTCAGG
GTCTAAGAAAAAGAAGACATCTGAGTCACAAAAGCAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG226928 representing NM_001145728
 Red=Cloning site Green=Tags(s)

MKDTQDIQEASKAMQNKVHEQEDKNEKQKQREDKLGVSLSLVHFSKMLGSVATTLPLSSSNSSGMPLGYY
 LSSPQISRVTISTTGLTSKATVVGSCSRVENSLEASPFVSPKKQDESPMIGDGEDYFLSLFGDSKKLTAH
 SNYTQKTLKYFSMILEEVGQFTSSSLGDVEIAEVNVKGLFVKLINSSLDKEMAIGDHILQQNVNGQTISL
 YRFLPNIVMQANSTVTVWAAASEAKHQPPSDFLWKEQDKFRASPCITILCKPNGQAIAWYTPIHVKQAW
 EKLDADVEFNRCVSVPTFRKRVFQWTASTATITKEKQDQPKKDISNYQVEQAQVLLKREKEIPPTVFPN
 RSPWCQNPVYSAHPYCLIEPHNTSTAGGRDRQPRTRSTRPNRASGSKKKKTSSESQKQ

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001145728

ORF Size: 1227 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001145728.2](#), [NP_001139200.1](#)

RefSeq Size: 1891 bp


RefSeq ORF: 1230 bp

Locus ID: 160492

UniProt ID: [Q8N9Z9](#)

Cytogenetics: 12p12.1

Product images:



Circular map for RG226928