

Product datasheet for **RG226952**

C17orf49 (NM_001142798) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C17orf49 (NM_001142798) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	C17orf49
Synonyms:	BAP18; HEPIS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG226952 representing NM_001142798 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACGTCAGCGTCCACAAAGGTCGGAGAGATCTTCTCGCGGCCGCGCCGCCTTCACGAAGCTCGGGG
AGCTGACGATGCAGCTGCATCCCGTGGCCGACTCTTCTCCTCGGGTGCGAAGTGACGGAGACGGAAAT
AGAGATGCTGAGGGCTGCTGTGAAGCGATTGGGGACGATCTTAATCACATCAGCTGTGCATCAAGGAA
CGGACAGTGGCCAGATAAAGGCCACTGTGAAACGCAAGGTATATGAAGATTCTGGCATCCCCCTCCAG
CTGAGTACCCCAAGAAAGGGCCCAAGAAGGTGGCATCTGGTGTCTGTACCTCCTCCAGCTGCCCTCC
TCCAGCAGCTCCAGTGTCCCTGAGGCCGGGGTCCCCCATAAAGAAACAGAAGGCTGATGTGACACTC
AGTGCTCTGAACGACTCCGATGCCAACAGTGACGTGGTGGATATTGAAGGCTAGGAGAACTCCTCCAG
CTAAGAACTCAACTTCGACCAGGACAGCCTGACCCTGGATTCTGGCCTTCTCATGACCTGCTGATCC
TCCTCTCTCTCTGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:	>RG226952 representing NM_001142798 Red=Cloning site Green=Tags(s)
-------------------	---

MTSASTKVGEIFSAAGAAFTKLGLTLMQLHPVADSSPAGAKWTETIEMLRAAVKRFGDDLNHISCVIKE
RTVAQIKATVKRKVYEDSGIPLPAESPKKGPKKVASGVLSPPPAAPPSSSSVPEAGGPPIKKQKADVTL
SALNDSANDVDIEGLGETPPAKKLNFDQDSLTLDSGLLMTSADPPLLSC

TRTRPLE - GFP Tag - V

Restriction Sites:	Sgfl-MluI
--------------------	-----------



[View online »](#)

Cloning Scheme:


ACCN: NM_001142798

ORF Size: 576 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142798.3](#)

RefSeq Size: 913 bp

RefSeq ORF: 579 bp


Locus ID: 124944

UniProt ID: [Q8IXM2](#)

Cytogenetics: 17p13.1

Gene Summary: Component of chromatin complexes such as the MLL1/MLL and NURF complexes.
[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG226952