

Product datasheet for **RG227129**

FXYD3 (NM_001136008) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: FXYD3 (NM_001136008) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: FXYD3
Synonyms: MAT8; PLML
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG227129 representing NM_001136008
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGAAGGTGACCCTGGCCTGCTTGTGTTCTGGCAGGCTTTCCTGTCCTGGACGCCAATGACCTAG
AAGATAAAACAGTCCTTCTACTATGACTGGCACAGCCTCCAGGTTGGCGGGCTCATCTGCGCTGGGGT
TCTGTGCGCCATGGGCATCATCATCGTCATGAGTGAGTGGAGGAGCTCGGGGAGCAGGCCGGCCGGGGC
TGGGGCTCCCCTCCCCTGACCACTCAGCTCTCCCAACAGGTGCAAATGCAAATGCAAGTTTGGCCAGA
AGTCCGGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG227129 representing NM_001136008
Red=Cloning site Green=Tags(s)
MQKVTLGLLVFLAGFPVLDANDLEDKNSPFYYDWHSLQVGLICAGVLCAMGIIIVMSEWRSSGEQAGRG
WGSPLTTQLSPTGAKCKCKFGQKSG


TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001136008

ORF Size: 288 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001136008.1</u> , <u>NP_001129480.1</u>
RefSeq Size:	778 bp
RefSeq ORF:	291 bp
Locus ID:	5349
UniProt ID:	<u>Q14802</u>
Cytogenetics:	19q13.12
Protein Families:	Ion Channels: Other, Transmembrane
Gene Summary:	This gene belongs to a small family of FXYD-domain containing regulators of Na ⁺ /K ⁺ ATPases which share a 35-amino acid signature sequence domain, beginning with the sequence PFXYD, and containing 7 invariant and 6 highly conserved amino acids. This gene encodes a cell membrane protein that may regulate the function of ion-pumps and ion-channels. This gene may also play a role in tumor progression. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, Oct 2008]