

Product datasheet for **RG227169**

SPART (NM_001142296) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPART (NM_001142296) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SPART
Synonyms:	SPG20; TAHCCP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG227169 representing NM_001142296
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGCAAGAGCCACAAAATGGAGAACCTGCTGAAATTAAGATCATCAGAGAAGCATATAAGAAGGCCT
TTTTATTTGTTAAACAAAGGTCTGAATACAGATGAATTAGGTCAGAAGGAAGAAGCAAAGAAGCTACTATAA
GCAAGGAATAGGACACCTGCTCAGAGGGATCAGCATTTTCATCAAAAGAGTCTGAACACACAGGACTGCG
TGGGAATCTGCTAGACAGATGCAACAGAAAATGAAAGAACTCTACAGAAATGTACGCACCAGGCTGGAAA
TTCTAGAGAAGGGTCTTGCCACTTCTCTGCAGAAATGATCTTCAGGAGGTGCCAAGTTATATCCAGAATT
TCCACCTAAAGACATGTGTGAAAAATTACCAGAGCCTCAGTCTTTTAGTTCAGCTCCTCAGCATGCTGAA
GTAATGGAAACACCTCAACTCCAAGTGCAGGGCAGTTGCTGCACCTGCTTCTGTCTTTACCATCAC
AAAGTTGTCCAGCAGAAGCTCCTCTGCTTATACTCCTCAAGCTGCTGAAGGTCCTACTACTGTATCCTA
TGGAACAGATTCTGGGGAGTTTTTCATCAGTTGGAGAGGAGTTTTATAGGAATCATTCTCAGCCACCGCT
CTTGAGACCTTAGGGCTGGATGCAGATGAATTGATTTTGATACCAAATGGAGTACAGATTTTTTTTGTA
ATCCTGCAGGGGAGGTTAGTGCACCTTCGTATCCTGGGTACCTTCGAATTGTGAGTTTTTTGGATAATTC
TCTCGATACGGTTCTAAACCGTCTCCCGGTTTCTTCAGGTTTGTGACTGGTTATATCCTCTAGTTCCT
GATAGATCTCCGTTCTGAAATGTACTGCGGGAGCCTACATGTTTCTGATACAATGCTACAAGCAGCAG
GATGCTTTGTGGGGTCTGCTCCTCTGAGTTACCAGAGGATGATAGAGAGCTCTTTGAGGATCTGTT
AAGGCAATGTCTGACCTTCGGCTCCAGGCAACTGGAACAGAGCAGAAGAAGAAAATGAATCCAAATC
CCTGGAAGAAGTAGACCCTCCTCTGACCAACTAAAAGAAGCCTCTGGCACTGATGTGAAACAGTTGGACC
AAGGCAATAAGGATGTACGTCATAAAGGAAAAGTGGAAAAAGGGCTAAAGATACTTCAAGTGAAGAAGT
TAACCTGAGTCACATTGTACCATGTGAGCCAGTTCAGAAAGAAAAGCCAAAAGAATTACATGAATGGAGT
GAAAAAGTGGCTCACAACATTTTGTGAGGTGCTTCTGGGTGAGTTGGGGTTTAGTCAAAGGTGCTGAGA
TTACTGTTAAGGCAATCCAGAAAGGTGCTTCTAAACTCCGAGAGCGGATTCAACCAGAAGAAAAACCCGT
GGAAGTTAGTCCAGCTGTACCAAGGGACTTTATATAGCGAAGCAAGCTACAGGAGGAGCAGAAAAGTC
AGTCAGTTCCTGGTTGATGGAGTTTGCCTGTAGCAAATGCGTTGAAAAAGAACTAGCTCCACATGTCA
AGAAGCATGGAAGCAAATTTGTTCCAGAATCTCTAAAAAAGACAAGATGGGAAATCTCCTCTGGATGG
TGCTATGGTTGTAGCAGCGAGTAGTGTCAAGGATTTTCAACTGTCTGGCAAGGATTGGAATGTGCAGCT
AAATGCATCGTTAACAATGTTTCAGCAGAACTGTACAACTGTCAGATACAAATACGGATATAATGCAG
GAGAAGCTACCCACCATGCGGTGGATTCTGCGGTCAATGTTGGCGTAACTGCCTACAATATAACAACAT
TGGTATCAAAGCAATGGTGAAGAAAAGTCAACACAAAACAGGACACACTCTCCTTGAGGACTATCAGATA
GTTGATAATTCTCAGAGGGAAAATCAAGAAGGAGCAGCAAATGTCAACGTGAGAGGGGAGAAGGATGAGC
AGACGAAGGAAGTAAAGGAGGCAAGAAGAAAGATAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG227169 representing NM_001142296
 Red=Cloning site Green=Tags(s)


MEQEPQNGEPAEIKIIREAYKKAFLFVNKGLNTDELGQKEEAKNYYKQGIGHLLRGISISSKESEHTGTG
 WESARQMQQMKETLQNVTRLEILEKGLATSLQNDLQEVPKLYPEFPPKDMCEKLEPQSFSSAPQHA
 VNGNTSTPSAGAVAAPASLSLPSQSCPAEAPPAYTPQAAEGHYTVSYGTDSGEFSSVGEFYNHSQPPP
 LETLGLDADELILIPNGVQIFFVNPAGEVSAPSYPGYLRIVRFLDNLDTVLNRPPGFLQVCDWLYPLVP
 DRSPVLKCTAGAYMFPDMLQAAGCFVGVLSSELPEDDRELFEDLLRQMSDLRLQANWNRAEEENEFI
 PGRTRPSSDQLKEASGTDVKQLDQGNKDVHRHKGKRAKDTSSSEEVNLSHIVPCEPVPEEKPKELHEWS
 EKVAHNILSGASWVSWGLVKGAIEITGKAIQKGASKLRERIQPEEKPEVSPAVTKGLYIAKQATGGA
 SQFLVDGVCTVANCVGKELAPHVKKHSGSKLVPELKKDKDGKSPLDGAMVVAASSVQGFSTVWQGLECA
 KCIVNNVSAETVQTVRYKYGYNAGEATHHAVDSAVNVGVTAYNINNIGIKAMVKKATQTGHTLLEDYQI
 VDNSQRENQEGAANVNRGEKDEQTKVEKAKKKDK

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001142296

ORF Size: 1998 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142296.1](#), [NP_001135768.1](#)

RefSeq Size: 5014 bp

RefSeq ORF: 2001 bp

Locus ID: 23111

UniProt ID: [Q8N0X7](#)

Cytogenetics: 13q13.3

Gene Summary: This gene encodes a protein containing a MIT (Microtubule Interacting and Trafficking molecule) domain, and is implicated in regulating endosomal trafficking and mitochondria function. The protein localizes to mitochondria and partially co-localizes with microtubules. Stimulation with epidermal growth factor (EGF) results in protein translocation to the plasma membrane, and the protein functions in the degradation and intracellular trafficking of EGF receptor. Multiple alternatively spliced variants, encoding the same protein, have been identified. Mutations associated with this gene cause autosomal recessive spastic paraplegia 20 (Troyer syndrome). [provided by RefSeq, Nov 2008]