

Product datasheet for **RG227208**

NBPF4 (NM_001143989) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NBPF4 (NM_001143989) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NBPF4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG227208 representing NM_001143989
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGTGGTATCTGCCGACCTTTGTCCAGCGAGAGGGCAGAGATGAACATCCTAGAAATCAACCAGGAAT
 TGGCTCGCAGCTGGCAGAGAGCAATCAGCAGTTCCGAGACCTCAAAGAGAAATTCCTTATAACTCAAGC
 TACTGCCTACTCCCTGGCCAACCACTGAGAAATACAAGTGTGAAGAGTACAAAGACATCATAGACTCT
 GTGCTGAGGGATGAAGTGCAGTCCATGGAGAAGCTGGCAGAGAAGCTCAGGCAAGCTGAGGAGCTCAGGC
 AGTATAAAGCCCTGGTTCACTCTCAGGCAAAAGAGCTGACCCAGTTACGGGAGAAGTTACGGGAAGGGAG
 AGATGCCTCCCGTGGCTGAACAAGCATCTGAAAACCTCCTCACTCTGATGACCCTGACAAGTCCAG
 GGTGAGGACCTCCGAGAGCAGCTGGCTGAGGGGCACAGGCTGGCAGAGCACCTTTGTCACAAGCTGAGCC
 CAGAAAATGATGAAGATGAAGATGAGGATGAAGACGACAAAGACGAGGAGTTGAGAAAGTACAGGAATC
 ACCTGCCCCCAGAGAGGTGCAGAAGACTGAAGAAAAGGAAGTCCCTCAGGACTCACTGGAGGAATGTGCT
 GTCACCTGTTCAAATAGTCACAACCTTCTAACTCCAACCAAGCCTCACAGGAGCACCAAAATCACATTTA
 AGGAACACGAAGTCGACTCTGCTCTGGTTGTAGAGAGTGAACACCTCATGATGAAGAGGAGGAAGCTCT
 AAACATTCCTCCAGAAAATCAAATGACCATGAGGAGGAGGAGGGAAAGCGCCAGTGCCCCCAGACAC
 CATGACAAGTCCAACCTTACCGGCATCGTGAAGTCTCTTTCTTGGCATTGGATGAACAGAAAGTTTGCT
 CCGCTCAGGATGTTGCCAGGGATTACTCCAATCCCAATGGGATGAAACCTCACTTGCTTCTCGAAAA
 GCAAAGTGATCTTGAAGAGGTGAAAGGACAAGAAACAGTTGCTCCAGGCTCAGCAGGGGACCGCTGAGA
 GTGGACAAGCATGAAATCCCCAGGAGTCACTGGATGGATGTTGCTTGACTCCTTCCATCCTTCTGACC
 TGACTCCCTCCTACCACCTTATTGGAGCACTTTGTAAGTCTTTTGAAGACAAGCAAGTCAAGTGGCTCT
 TGTAGACAAAATTAAAAAGGATCAAGAGGAGATAGAAGACCAAAAGCCACCATGCCCCAGGCTCAGCCAG
 GAGCTGCCAGAGGTGAAGGAGCAGGAAGTCCCAGAGGACTCTGTGAATGAAGTTTACTTGACTCCCTCAG
 TTCACCATGACGTGTCTGACTGCCACGACCTTATAGCAGCACCTTGCTCATTGGAGGATCAGCTTGC
 CTGCTCTGCTCTGGATGTAGCTCCCCACCGAGGCGGCTGTCCCAAGGGACTTGAGTGAGAGCTTG
 AGCCACCACCACTCAGAGGTGCAAGTTTACAGGCACAGCTGGAACCAAGCACCTTGGTCCCAAGTTGTC
 TGCGACTACAGCTGGATCAAGGTTTCACTGTGGGAACGCTTGCCCAAGCGGGCTTCTCCACCAC
 CTGCAGCTTCTCAGCAATGCTGATTCTGGGAACCAATGGCCCTTCAAGAGCTGGTTTTAGAGCCCTCT
 CTGGGGATGAAGAACCCTCCCCAGCTGGAAGATGATGCACTTGAAGGCTCAGCAAGCAACACACAAGGGC
 GTCAAGTCACTGGCCGATTTCGTGCCTCCCTTGTCTGATACTGAAGACCATCAGAAGAAGACTCCCGTT
 CAGCAAGTGGAGACTGGCATTGATTCGCTGGCCCGCATGCTGAGAGCGCAGAGATACCAATACTGCT
 GGAAGGACGCAAAGGATGGCAGGA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:


>RG227208 representing NM_001143989
 Red=Cloning site Green=Tags(s)

MVVSADPLSSERAEMNILEINQELRSQLAESNQFRDLKEKFLITQATAYSLANQLKKYKCEEYKDIIDS
 VLRDELQSMELAEKLRAEELRQYKALVHSQAKELTQLREKLREGDASRWLNKHLKTLTPDDPKSQ
 GQDLREQLAEGHRLAEHLVHKLSPENDEDEDEDDKDEEVEKVQESPAPREVQKTEEKEVPQDSLEECA
 VTCNSNHNPSNSNQPHRSTKITFKEHEVDSALVVESEHPHDEEEALNIPPENQNDHEEEGKAPVPPRH
 HDKSNSYRHREVSFLALDEQKVCSAQDVARDYSNPKWDETSLGFLEKQSDLEEVEKQETVAPRLSRGPLR
 VDKHEIPQESLDGCLTPSILPDLTSPYHPYWSTLYSFEDKQVSLALVDKIKKDQEEIEDQSPPCPRLSQ
 ELPEVKEQEVPEDSVNEVYLTPSVHHDVSDCHQPYSSSTLSLEDQLACSALDVASPTAACPQGTWSGDL
 SHHQSEVQVSQAQLEPSTLVPSCRLQLDQGFHCGNLAQRGLSSTTCSFSANADSGNQWPFQELVLEPS
 LGMKNPPQLEDDALEGSANTQGRQVTGRIRASLVILKTI RRRLPF SKWRLAFRFAGPHAESAEPNTA
 GRTQRMAG

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001143989

ORF Size: 1914 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001143989.2](#), [NP_001137461.1](#)

RefSeq Size: 3375 bp

RefSeq ORF: 1917 bp


Locus ID: 148545

UniProt ID: [Q96M43](#)

Cytogenetics: 1p13.3

Gene Summary: This gene is a member of the neuroblastoma breakpoint family (NBPF) which consists of dozens of recently duplicated genes primarily located in segmental duplications on human chromosome 1. This gene family has experienced its greatest expansion within the human lineage and has expanded, to a lesser extent, among primates in general. Members of this gene family are characterized by tandemly repeated copies of DUF1220 protein domains. Gene copy number variations in the human chromosomal region 1q21.1, where most DUF1220 domains are located, have been implicated in a number of developmental and neurogenetic diseases such as microcephaly, macrocephaly, autism, schizophrenia, cognitive disability, congenital heart disease, neuroblastoma, and congenital kidney and urinary tract anomalies. Altered expression of some gene family members is associated with several types of cancer. This gene family contains numerous pseudogenes. [provided by RefSeq, Mar 2013]

Product images:



Circular map for RG227208