

Product datasheet for **RG227260**

PSD93 (DLG2) (NM_001142700) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PSD93 (DLG2) (NM_001142700) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	DLG2
Synonyms:	chapsyn-110; PPP1R58; PSD-93; PSD93
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG227260 representing NM_001142700
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCAACGGCCAAGTGTTCCTCGGGCTGAGAATTATCAGCTTTTGTGGGATACTATTGCTTCCTTAAAC
 AATGTGAACAAGCTATGCAACATGCATTTATCCGGTCAATGGGACAGAAATTGAATATGAATTTGAAGA
 AATTACACTGGAGAGGGGAATTCTGGCCTGGGATTAGTATTGCTGGGGGACAGATAATCCCCACATT
 GGAGATGACCCTGGCATAATTTATACGAAGATTATACCAGGAGGTGCTGCAGCAGAGGATGGCAGACTCA
 GGGTCAATGATTGTATCTTGCGGTGAATGAGGTTGATGTGCAGAGTTTCCACAGTAAAGCGGTGGA
 AGCCCTGAAGGAAGCAGGTCTATCGTTCGGCTGTATGTGCGTAGAAGACGACCTATTTGGAGACCGTT
 GTGAAATCAAAGTGTCAAAGGCCCTAAAGGTTTAGGCTTCAGTATTGCAGGAGGTGTGGGAACCAAC
 ACATTCCTGGAGACAACAGCATTATGTAATAAAATATAGATGGAGGAGCTGCACAAAAGATGAAG
 GTTGCAAGTAGGAGATAGACTACTAATGGTAAACAACACAGTTTGAAGAAGTAACACACGAAGAGGCA
 GTAGCAATATTAAGAACACATCAGAGGTAGTTTATTTAAAAGTTGGCAAACCCACTACCATTTATATGA
 CTGATCCTTATGGTCCACCTGATTAATACTACTCTTATTCCTCCACCAATGGAAAACCATCTACTCTCTGG
 CAACAATGGCACTTTAGAATATAAAACCTCCCTGCCACCCATCTCTCCAGGAAGTACTACCAATTCCA
 AAGCACATGCTTGTGACGACGACTACACCAGTCATCCCAACATAGCACCGCAACTCGTCAGCCTTCAA
 TGACTCTCCAACGGGCCGTCTCCCTGGAAGGAGAGCCTCGCAAGGTAGTCTGCACAAAGGCTCCACTGG
 CCTGGGCTTCAACATTGTGCGTGGGAAGATGGAGAAGTATTTTGTGTCCTTCATTCTGGCTGGTGA
 CCAGCAGACCTAAGTGGGAGCTCCAGAGAGGAGACCAGATCCTATCGGTGAATGGCATTGACCTCCGTG
 GTGCATCCACGAGCAGGCAGCTGCTGCACTAAAGGGGGCTGGACAGACAGTGACGATTATAGCACATA
 TCAACCTGAAGATTACGCTCGATTTGAGGCCAAAATCCATGACCTACGAGAGCAGATGATGAACACAGC
 ATGAGCTCCGGTCCGGATCCCTGCGAACCAATCAGAAACGCTCCCTCTACGTCAGAGCCATGTTGACT
 ACGACAAGAGCAAGGACAGTGGGCTGCCAAGTCAAGGACTTAGTTTTAAATATGGAGATATTCTCCACGT
 TATCAATGCCTCTGATGATGAGTGGTGGCAAGCCAGGAGAGTCATGCTGGAGGGAGACAGTGAGGAGATG
 GGGTCATCCCCAGCAAAGGAGGGTGGAAAGAAAGAACGTGCCGATTGAAGACAGTGAAGTTTAAATG
 CCAAACCTGGAGTGATTGATTGAAAGGGGACATCCCCGATTAGGTGACGACGGTTATGGAACAAAGAC
 TCTGAGAGGACAAGAAGACCTCATTCTTCTATGAGCCTGTTACAAGGCAGGAAATAAACTACACCCGG
 CCGGTGATTATCCTGGGCCCATGAAGGATCGGATCAATGACGACTTGATATCTGAATCCCTGATAAAT
 TTGGCTCCTGTGTGCCTATACTACGAGGCCAAAGCGAGACTACGAGGTGGATGGCAGAGACTATCACTT
 TGTCAATTCAGAGAACAATGGAGAAAGATATCCAAGAGCACAAGTTTATAGAAGCCGGCCAGTACAAT
 GACAATTTATATGGAACCAAGTGTGCAGTCTGTGAGATTTGTAGCAGAAAGAGGCAAAACACTGTATACTTG
 ATGTATCAGGAAATGCTATCAAGCGTTACAAGTTGCCAGCTCTATCCCATTGCCATCTTCATAAAACC
 CAGGTCTCTGGAACCTCTTATGGAGATGAATAAGCGTCTAACAGAGGAACAAGCCAAGAAAACCTATGAT
 CGAGCAATTAAGCTAGAACAAGAATTTGGAGAATTTTACAGCTATTGTCCAAGGAGATACTTTAGAAG
 ATATATAACCAATGCAAGCTTGTATTGAAGAGCAATCTGGCCCTTTCATCTGGATCCCTCAAAGGA
 AAAGTTA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG227260 representing NM_001142700
 Red=Cloning site Green=Tags(s)


MQRPSVSRAENYQLLWDTIASLKQCEQAMQHAFIPVNGTEIEYEFEEITLERGNSGLGFSIAGGTDNPHI
 GDDPGIFITKIIPGAAAEDGRLRVNDCILRVNEVDVSEVSHSKAVEALKEAGSIVRLYVRRRPILETV
 VEIKLFKGPKGLGFSIAGGVGNQHIPGDNISYVTKIIDGGAAQKDGRLQVGDRLLMVNNYSLEEVTHEEA
 VAILKNTSEVVYLKVGKPTTIYMTDPYGPDPDITHSYSPPMENHLLSGNNGTLEYKTSLPPISPGRYSPIP
 KHMLVDDDYTSHSQHSTATRQPSMTLQRAVSLEGEPRKVVHLHGSTGLGFNIVGGEDGEGIFVSFILAGG
 PADLSGELQRGDQILSVNGIDLRGASHEQAAAALKGAGQVTVIIAQYQPEDYARFEAKIHDLREQMMNHS
 MSSGSGSLRTNQKRSLYVRAMFDYDKSKDGLPSQGLSFKYGDILHVINASDDEWWQARRVMLEGDSEEM
 GVIPSKRRVERKERARLKTVKFNAKPGVIDSKGDIPGLGDDGYGKTLRGQEDLILSYEVPTRQEIYNYTR
 PVIILGPMKDRINDDLISEFPDKFGSCVPHTTRPKRDYVDGRDYHFVISREQMEKDIQEHKFIEAGQYN
 DNLYGTSVQSVRFVAERGGKHCILDVSGNAIKRLQVAQLYPIAIFIKPRSLEPLMEMNKRLEEQAKKTYD
 RAIKLEQEFGEYFTAIVQGDITLEIYNQCKLVIEEQSGPFIWIPSKEKL

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001142700

ORF Size: 2247 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142700.2](#)

RefSeq Size: 7096 bp

RefSeq ORF: 2250 bp

Locus ID: 1740

UniProt ID: [Q15700](#)

Cytogenetics: 11q14.1

Protein Families: Druggable Genome

Gene Summary: This gene encodes a member of the membrane-associated guanylate kinase (MAGUK) family. The encoded protein forms a heterodimer with a related family member that may interact at postsynaptic sites to form a multimeric scaffold for the clustering of receptors, ion channels, and associated signaling proteins. Multiple transcript variants encoding different isoforms have been found for this gene. Additional transcript variants have been described, but their full-length nature is not known. [provided by RefSeq, Dec 2008]