

Product datasheet for **RG227281**

Periostin (POSTN) (NM_001135934) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Periostin (POSTN) (NM_001135934) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Periostin
Synonyms:	OSF-2; OSF2; PDLPOSTN; PN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG227281 representing NM_001135934
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATTCCCTTTTACCCATGTTTTCTCTACTATTGCTGCTTATTGTTAACCCATAAACGCCAACAAATC
 ATTATGACAAGATCTTGGCTCATAGTCGTATCAGGGTTCGGGACCAAGGCCCAAATGTCTGTGCCCTTCA
 ACAGATTTTGGGCACAAAAAGAAATACTTCAGCACTTGTAAAGAACTGGTATAAAAAAGTCCATCTGTGGA
 CAGAAAACGACTGTGTTATATGAATGTTGCCCTGGTTATATGAGAATGGAAGGAATGAAAGGCTGCCAG
 CAGTTTTGCCATTGACCATGTTTATGGCACTCTGGGCATCGTGGGAGCCACCACAACGCAGCGCTATTC
 TGACGCCTCAAACTGAGGGAGGAGATCGAGGGAAAGGGATCCTTCACTTACTTTGCACCGAGTAATGAG
 GCTTGGGACAACCTGGATTCTGATATCCGTAGAGGTTTGGAGAGCAACGTGAATGTTGAATTACTGAATG
 CTTTACATAGTCACATGATTAATAAGAGAATGTTGACCAAGGACTTAAAAATGGCATGATTATTCCTTC
 AATGTATAACAATTTGGGGCTTTTCATTAACCATTATCCTAATGGGGTGTCACTGTTAATTGTGCTCGA
 ATCATCCATGGGAACAGATTGCAACAAATGGTGTGTCCATGTCAATTGACCGTGTGCTTACACAAATTG
 GTACCTCAATTCAGACTTCATTGAAGCAGAAGATGACCTTTCATCTTTAGAGCAGCTGCCATCACATC
 GGACATATTGGAGGCCCTTGAAGAGACGGTCACTTCACACTCTTTGCTCCCACCAATGAGGCTTTTGAG
 AAACCTCCACGAGGTGCTCTAGAAAGGATCATGGGAGACAAAGTGGCTTCCGAAGCTCTTATGAAGTACC
 ACATCTTAAACTCTCCAGTGTCTGAGTCTATTATGGGAGGAGCAGTCTTTGAGACGCTGGAAGGAAA
 TACAATTGAGATAGGATGTGACGGTGACAGTATAACAGTAAATGGAATCAAAATGGTGAACAAAAAGGAT
 ATTGTGACAAATAATGGTGTGATCCATTTGATTGATCAGGTCCTAATTCCTGATTCTGCCAAACAAGTTA
 TTGAGCTGGCTGGAAAACAGCAAACCACCTTCACGGATCTGTGGCCCAATTAGGCTTGGCATCTGCTCT
 GAGGCCAGATGGAGAATACACTTTGCTGGCACCTGTGAATAATGCATTTTCTGATGATACTCTCAGCATG
 GATCAGCGCTCCTTAAATTAATTCTGCAGAATCACATATTGAAAGTAAAAGTTGGCCTTATGAGCTTT
 ACAACGGGCAAATACTGGAAACCATCGGAGGCAACAGCTCAGAGTCTTCGTATATCGTACAGCTGTCTG
 CATTGAAAATTCATGCATGGAGAAAGGAGTAAGCAAGGGAGAAACGGTGCATTACATATTCGCGGAG
 ATCATCAAGCCAGCAGAGAAATCCCTCCATGAAAAGTTAAAACAAGATAAGCGCTTAGCACCTTCTCA
 GCCTACTTGAAGCTGCAGACTTGAAGAGCTCCTGACACAACCTGGAGACTGGACATTATTTGTGCCAAC
 CAATGATGCTTTTAAAGGAATGACTAGTGAAGAAAAAGAAATTCGTACGGGACAAAAATGCTCTTCAA
 AACATCATTCTTTATCACCTGACACCAGGATTTTCATTGAAAAGGATTTGAACCTGGTGTACTAACA
 TTTTAAAGACCACACAAGGAAGCAAAATCTTTCTGAAAGAAGTAAATGATACACTTCTGGTGAATGAATT
 GAAATCAAAAGAACTGACATCATGACAACAAATGGTGAATTCATGTTGTAGATAAACTCCTCTATCCA
 GCAGACACACCTGTTGGAATGATCAACTGCTGGAAATACTTAATAAATTAATCAAAATACATCCAAATTA
 AGTTTGTTCGTGGTAGCACCTTCAAAGAAATCCCCGTGACTGTCTATAAGCCAATTATTAATAAATACAC
 CAAAATCATTGATGGAGTGCCTGTGAAATAACTGAAAAAGAGACACGAGAAGAACGAATCATTACAGGT
 CCTGAAATAAATACTAGGATTTCTACTGGAGGTGGAGAAACAGAAGAACTCTGAAGAAATGTTAC
 AAGAAGAGGTACCAAGGTACCAAAATTCATTGAAGGTGGTATGGTCATTTATTTGAAGATGAAGAAAT
 TAAAAGACTGCTTCAGGGAGACACACCCGTGAGGAAGTTGCAAGCCAACAAAAAGTTCAAGGATCTAGA
 AGACGATTAAGGGAAGGTCGTTCTCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG227281 representing NM_001135934
 Red=Cloning site Green=Tags(s)

MIPFLPMFSLLLLLIVNPINANNHYDKILAHSRIRGRDQGNVICALQQILGTTKKYFSTCKNWKYSICG
 QKTTVL YECCPGYMRMEGMKGC PAVLPIDHVVYGLGIVGATTTQRYSDASKLREEIEGKGSFTYFAPSNE
 AWDNLDSDIRRGLESNVNVELLNALHSHMINKRMLTKDLKNGMIIPSMYNNLGLFINHYPNGVVTVNCAR
 I IHGNQIATNGVVHVIDRVL TQIGTSIQDFIEAEDDLSSFRAAAITSDILEALGRDGHFTLFAPTNEAFE
 KLPRGVLERIMGDKVASEALMKYHILNLTQCSESIMGGAVFETLEGNTIEIGCDGDSITVNGIKMVNKKD
 IVTNGV IHLIDQVLIPDSAKQVIELAGKQQTTFDLVAQLGLASALRPDGEYTL LAPVNNAFSDDTL SM
 DQRL LKLILQNHILKVKVGLNELYNGQILETIGGKQLRVFYRTAVCIENSCMEKGSQGRNGAIHIFRE
 I IKPAEKSLHEK LKQDKRFSTFLSLEADLKELLTQPGDWTLFVPTNDAFKGMTSEEKEILIRDKNALQ
 NIILYHLTPGVFIGKGFEPGVN ILKTTQGSKIFLKEVNDLLVNELKSKESDIMTNGVIHVVDKLLYP
 ADTPVGN DQ LLEILNKL IKYIQIKFVRG STFKEIPVTYKPI IKKYTKI IDGVPVEITEKETREERIITG
 PEIKYTRISTGGGETEETLKKLLQEEVTKVTKFIEGGDGHLEFEDEEIKRLLQGDTPVRKLQANKKVQGS
 RRLREGRSQ

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_001135934

ORF Size: 2337 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001135934.2
RefSeq Size:	3219 bp
RefSeq ORF:	2340 bp
Locus ID:	10631
UniProt ID:	Q15063
Cytogenetics:	13q13.3
Protein Families:	Druggable Genome, Secreted Protein
Gene Summary:	<p>This gene encodes a secreted extracellular matrix protein that functions in tissue development and regeneration, including wound healing, and ventricular remodeling following myocardial infarction. The encoded protein binds to integrins to support adhesion and migration of epithelial cells. This protein plays a role in cancer stem cell maintenance and metastasis. Mice lacking this gene exhibit cardiac valve disease, and skeletal and dental defects. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]</p>

Product images:



Circular map for RG227281