

Product datasheet for **RG227290**

SORBS2 (NM_001145672) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SORBS2 (NM_001145672) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SORBS2
Synonyms:	ARGBP2; PRO0618
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG227290 representing NM_001145672
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAACACAGGGCGTGATTCTCAGTCACCAGACTCAGCTTGGAGATCTTACAATGATGGCAATCAGGAGA
CACTGAACGGAGATGCTACATATTCCTCTCTTGCAGCAAAGGTTTTAGAAGCGTTCGACCAAACCTACA
AGATAAAAGATCACCAACTCAGAGCCAGATAACAGTGAATGGAACTCAGGAGGTGCCGTGAGTCCCATG
AGTTACTATCAGAGGCCGTTTTCCCTCGGCATATTCTCTCCAGCCTCACTCACTCCAGCATTGTCA
TGCAGCACGGCACATCCCTCGATTCCACAGACACATATCCCAGCATGCGCAGTCTCTGGATGGCACCAC
CAGCAGCTCTATCCCCTGTACCGATCCTCAGAGGAAGAGAAGAGAGTGACAGTCATCAAAGCCCCGCAT
TACCCAGGGATCGGGCCCGTGGATGAATCCGGAATCCCCACAGCAATTAGAACGACAGTCGACCGGCCCA
AGGACTGGTACAAGACGATGTTAAGCAAATTCACATGGTGCACAAGCCGGATGATGACACAGACATGTA
TAATACTCCTTATACATACAATGCAGGTCTGTACAACCCACCCTACAGTGCTCAGTCACACCCTGTGCA
AAGACCCAAACCTACAGACCTCTTCCAAAAGCCACTCCGACAACAGCCCAATGCCTTTAAGGATGCGT
CCTCCCCAGTGCCTCCCCACATGTTCCACCTCCAGTCCCGCCGCTTCGACCAAGAGATCGGTCTTCAAC
AGAAAAGCATGACTGGGATCCTCCAGACAGAAAAGTGGACACAAGAAAATTTTCGGTCTGAGCCAAGGAGT
ATTTTTGAATAAGAACCTGGCAAGTCATCAATCTTTCAGCATGAAAGACCACCTCCTAAAAAGCCTCTGG
ACTATGTTCAAGATCATTCTTCTGGTGTTCCTCAATGAGGCCTCCTGTATCAGTCCTCTATAGACAGAAG
CCTGGAAGACCCATGAGTTCGCAAGCATGGCCAGTACTTCAGGAAGCGGAGGAAGAGCGAGCCTGCA
GTGGTCCACCACGGGGCTTGGGAGATCAAAGTGCAGCAGGACTAGCCCAGGCCGAGTGGACCTCCAG
GATCAAGCACCACCTTTACAAGTCTTTCAGTCTTCTCCTTCTTCCCCATCAAGAGCAAAAAGACCG
TGAGTCCCCTAGAAGTTACTCATCCACTTTGACTGACATGGGGAGAAGTGCACCAAGGGAAAGAAGAGGA
ACTCCAGAAAAAGAGAAATTGCCTGCAAAAAGCTGTTTATGATTTTAAGGCTCAGACATCTAAGGAGTTGT
CATTTAAGAAAGGAGATACTGTCTACATCCTCAGGAAAATTGATCAAAATTGGTATGAGGGAGAACACCA
CGGGAGAGTGGGCATCTTCCCGATCTCATACGTAGAGAAACTCACACCTCCTGAGAAAGCACAGCCTGCA
AGACCACCTCCGCCAGCCAGCCCGGAGAAATCGGAGAAGCTATAGCCAAATACTCAACTTCAACGCAGACA
CAAATGTGGAGCTGTCACTGAGAAAGGGAGATAGAGTTATTCTTCTTAAAAGAGTTGATCAAACTGGTA
TGAAGGTAATAATCCAGGAACCAACAGACAAGGCATCTTCCCTGTTTCTATGTGGAGGTCGTCAAGAAG
AACACAAAAGGTGCTGAGGACTACCTGACCCTCCAATACCCACAGCTATTCTAGTGATAGGATTCACA
GCTTGAGCTCAAATAAGCCACAGCGTCTGTGTTACTCATGAAAATATTCAAGGTGGGGGGGAACCGTT
TCAGGCTCTGATAACTATACTCCAGGAATGAAGATGAGCTGGAGCTCAGAGAAAGTGTGTCATTGAT
GTCATGGAAAAGTGTGATGACGGCTGGTTTGTGGGGACCTCAAGAAGAACCAAAATCTTTGGTACTTTCC
CCGGAACTACGTCAAGAGGCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG227290 representing NM_001145672
 Red=Cloning site Green=Tags(s)


MNTGRDSQSPDSAWRSYNDGNQETLNGDATYSSLAAGFRSVRPNLQDKRSPTQSQITVNGNSGGAVSPM
 SYYQRPFSAYSLSASLNSSIVMQHGTSLDSTDTYPQHAQSLDGTSSSIPLYRSSEEEKRVTVIKAPH
 YPGIGPVDESGIPTAIRTTVDRPKDWYKTMFKQIHMVHKPDDDDTMYNTPYTYNAGLYNPPYSAQSHPA
 KTQTYRPLSKSHSDNSPNAFKDASSPVPPPVPVPLRPRDRSSTEKHDWDDPPDRKVDTRKFRSEPRS
 IFEYEPGKSSILQHERPPPKPLDYVQDHSSGVFNEASLYQSSIDRSLERPMSSASMASDFRKRKSEPA
 VGPPRGLGDQASRTSPGRVDLPGSSTLTKSFTSSSPSPSRAKDRESPRSSTLTDMGRSAPRERRG
 TPEKEKLPKAVYDFKAQTSKELSFKKGDTVYILRKIDQNWYEGEHHRVGIFFPISYVEKLTPEKAQPA
 RPPPPAQGEIGEIAKYNFNADTNVELSLRKGDRIILLKRVQDNWYEGKIPGTNRQGIFFVSYVEVVKK
 NTKGAEDYDPPPIPHSYSSDRIHSLSSNKQRPVFTHEINIQGGGEPFQALYNYTPRNEDELELRESDDVID
 VMEKDDGWVGTSRRTKFFGTFFGNYVKRL

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001145672

ORF Size: 1983 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145672.1](#), [NP_001139144.1](#)

RefSeq Size: 4356 bp

RefSeq ORF: 1986 bp

Locus ID: 8470

UniProt ID: [O94875](#)

Cytogenetics: 4q35.1

Gene Summary: Arg and c-Abl represent the mammalian members of the Abelson family of non-receptor protein-tyrosine kinases. They interact with the Arg/Abl binding proteins via the SH3 domains present in the carboxy end of the latter group of proteins. This gene encodes the sorbin and SH3 domain containing 2 protein. It has three C-terminal SH3 domains and an N-terminal sorbin homology (SoHo) domain that interacts with lipid raft proteins. The subcellular localization of this protein in epithelial and cardiac muscle cells suggests that it functions as an adapter protein to assemble signaling complexes in stress fibers, and that it is a potential link between Abl family kinases and the actin cytoskeleton. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]