

Product datasheet for **RG227325**

SORBS2 (NM_001145675) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SORBS2 (NM_001145675) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SORBS2
Synonyms:	ARGBP2; PRO0618
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG227325 representing NM_001145675
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAAGCAACAACACCTTTGCAGACAGTCGACCGGCCAAGGACTGGTACAAGACGATGTTTAAGCAAA
 TTCACATGGTGCACAAGCCGGATGATGACACAGACATGTATAACTCCTTATACATACAATGCAGGTCT
 GTACAACCACCCTACAGTGCTCAGTCACACCCTGCTGCAAAGACCCAAACCTACAGACCTCTTTCCAAA
 AGCCACTCCGACAACAGCCCAATGCCTTTAAGGATGCGTCTCCCCAGTGCCTCCCCACATGTTCCAC
 CTCCAGTCCCAGCGCTTCGACCAAGAGATCGGTCTTCAACAGAAAAGCATGACTGGGATCTCCAGACAG
 AAAAGTGGACACAAGAAAATTCGGTCTGAGCCAAGGAGTATTTTTGAATATGAACCTGGCAAGTCATCA
 ATTCTTCAGCATGAAAGACCAGCCTCTTGTATCAGTCTCTATAGACAGAAGCCTGAAAGACCCATGA
 GTTCTGCAAGCATGGCCAGTGACTTCAGGAAGCGGAGGAAGAGCGAGCCTGCAGTGGTCCACCACGGGG
 CTTGGGAGATCAAAGTGCAGCAGGACTAGCCAGGCCGAGTGGACCTCCAGGATCAAGCACCACTCTT
 ACAAAAGTCTTCACTAGCTCTTCTCCTTCTCCCATCAAGAGCAAAAGACCGTGAGTCCCTTAGAAGTT
 ACTCATCCACTTTGACTGACATGGGGAGAAGTGCACCAAGGGAAAAGAGGAACTCCAGAAAAAGAGAA
 ATTGCTGCAAAAAGCTGTTTATGATTTTAAAGGCTCAGACATCTAAGGAGTTGTCATTTAAGAAAAGGAGAT
 ACTGTCTACATCCTCAGGAAAATGATCAAAAATGGTATGAGGGAGAACCACCGGGAGAGTGGGCATCT
 TCCCGATCTACGTAGAGAACTCACACCTCCTGAGAAAGCACAGCCTGCAAGACCACCTCCGCCAGC
 CCAGCCCGGAGAAATCGGAGAAGCTATAGCCAAATACAACCTCAACGCAGACACAATGTGGAGCTGTCA
 CTGAGAAAAGGAGATAGAGTTATTCTTCTTAAAAGAGTTGATCAAACTGGTATGAAGGTAATCCAG
 GAACCAACAGACAAGGCATCTTCCCTGTTTCTATGTGGAGTCTGCAAGAAGAACAACAAAAGTGTCTGA
 GGACTACCCTGACCTCCAATACCCACAGCTATTCTAGTGATAGGATTACAGCTTGAGCTCAAATAAG
 CCACAGCGTCTGTGTTTACTCATGAAAATATTCAAGGTGGGGGGAACCGTTTCAGGCTCTGTAACT
 ATACTCCAGGAATGAAGATGAGCTGGAGCTCAGAGAAAGTGATGTCATTGATGTCATGGAAAAGTGTGA
 TGACGGCTGGTTTGTGGGACCTCAAGAAGAACCAAAATCTTTGGTACTTTCCCGGAACTACGTCAAG
 AGGCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG227325 representing NM_001145675
 Red=Cloning site Green=Tags(s)

MKATTPLQTVDRPKDWYKTMFKQIHMVHKPDDDDTDMYNTPTYNAGLYNPPYSAQSHPAAKTQTYRPLSK
 SHSDNSPNAFKDASSPVPPPHVPPVPLRPRDRSSTEKHDWPPDRKVDTRKFRSEPRSIFFEYEPGKSS
 ILQHERPASLYQSSIDRSLERPMSASMASDFRKRKSEPAVGPPRGLGDQASRTSPGRVDLPGSSTTL
 TKSFTSSPSSPSRAKDRESRYSSTLDMGRSAPRERRGTPEKEKLPKAVYDFKAQTSKELSFKKGD
 TVYILRKIDQNWYEGEHHGRVGFIPISYVEKLTPEKAQPARPPPAQPGEIGEIAKYNFNADTNVELS
 LRKGDRVILLKRVDQNWYEGKIPGTNRQGFIPVSYVEVVKKNTKGAEDYDPPPIPHSYSSDRIHSLSSNK
 PQRPFVTHENIQGGGEPFQALYNYTPRNEDELELRESDVIDVMEKCDDGWFVGTSRRTKFFGTFPGNYVK
 RL



TRTRPLE - GFP Tag - V


Restriction Sites:

Sgfl-Mlul

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_001145675

ORF Size: 1476 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145675.2](#)

RefSeq Size: 3993 bp

RefSeq ORF: 1479 bp

Locus ID: 8470

Cytogenetics: 4q35.1

Gene Summary: Arg and c-Abl represent the mammalian members of the Abelson family of non-receptor protein-tyrosine kinases. They interact with the Arg/Abl binding proteins via the SH3 domains present in the carboxy end of the latter group of proteins. This gene encodes the sorbin and SH3 domain containing 2 protein. It has three C-terminal SH3 domains and an N-terminal sorbin homology (SoHo) domain that interacts with lipid raft proteins. The subcellular localization of this protein in epithelial and cardiac muscle cells suggests that it functions as an adapter protein to assemble signaling complexes in stress fibers, and that it is a potential link between Abl family kinases and the actin cytoskeleton. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]