

## Product datasheet for **RG227328**

### Annexin A2 (ANXA2) (NM\_001136015) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Annexin A2 (ANXA2) (NM_001136015) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ANXA2
Synonyms:	ANX2; ANX2L4; CAL1H; HEL-S-270; LIP2; LPC2; LPC2D; P36; PAP-IV
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG227328 representing NM_001136015 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCTACTGTTACGAAATCCTGTGCAAGCTCAGCTTGGAGGGTGATCACTCTACACCCCAAGTGCAT  
ATGGGTCTGTCAAAGCCTATACTAACTTTGATGCTGAGCGGGATGCTTTGAACATTGAAACAGCCATCAA  
GACCAAAGGTGTGGATGAGGTCACCATTGTCAACATTTTGACCAACCGCAGCAATGCACAGAGACAGGAT  
ATTGCCTTCGCCTACCAGAGAAGGACAAAAAGGAACTTGCATCAGCACTGAAGTCAGCCTTATCTGGCC  
ACCTGGAGACGGTGATTTTGGCCTATTGAAGACACCTGCTCAGTATGACGCTTCTGAGCTAAAAGCTTC  
CATGAAGGGGCTGGAAACCGACGAGGACTCTCTATTGAGATCATCTGCTCCAGAACCAACCAGGAGCTG  
CAGGAAATTAACAGAGTCTACAAGGAAATGTACAAGACTGATCTGGAGAAGGACATTATTTCCGGACACAT  
CTGGTGACTTCCGCAAGCTGATGGTTGCCTGGCAAAGGGTAGAAGAGCAGAGGATGGCTCTGTCATTGA  
TTATGAACTGATTGACCAAGATGCTCGGGATCTCTATGACGCTGGAGTGAAGAGGAAAGGAACTGATGTT  
CCCAAGTGGATCAGCATCATGACCGAGCGGAGCGTGCCCCACCTCCAGAAAGTATTTGATAGGTACAAGA  
GTTACAGCCCTTATGACATGTTGGAAAGCATCAGGAAAGAGGTTAAAGGAGACCTGGAAAATGCTTTCT  
GAACCTGGTTCAGTGCATTGAGAAAGCCCTGTATTTGCTGATCGGCTGTATGACTCCATGAAGGGC  
AAGGGGACGCGAGATAAGGTCCTGATCAGAATCATGGTCTCCCGCAGTGAAGTGGACATGTTGAAAATTA  
GGTCTGAATTCAGAGAAAGTACGGCAAGTCCCTGTACTATTATATCCAGCAAGACACTAAGGGCGACTA  
CCAGAAAGCGCTGCTGTACCTGTGTGGTGGAGATGAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG227328 representing NM\_001136015  
Red=Cloning site Green=Tags(s)


MSTVHEILCKLSLEGDHSTPPSAYGVSVKAYTNFDAERDALNIETAIKTKGVDEVTIVNIL TNRSNAQRQD  
 IAFAYQRRTKKELASALKSALSGHLETVILGLLKTPAQYDASELKASMKGLGTDEDSLIEIICSRTNQEL  
 QEINRVYKEMYKTDLEKDIISDTSGDFRKL MVALAKGRR AEDGSVIDYELIDQDARDLYDAGVKRKGTDV  
 PKWISIMTERSVPHLQKVFDRYKSYSPYDML ESIRKEVKGDLENAFLNLVQCIQNKPLYFADRLYDSMKG  
 KGTRDKVLIIRIMVSRSEVDMLKIRSEFKRKYGKSLYYYIQQDTKGDYQKALLYLCGGDD

TRTRPLE - GFP Tag - V


**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001136015

**ORF Size:** 1017 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001136015.3</a>
<b>RefSeq Size:</b>	1665 bp
<b>RefSeq ORF:</b>	1020 bp
<b>Locus ID:</b>	302
<b>UniProt ID:</b>	<a href="#">P07355</a>
<b>Cytogenetics:</b>	15q22.2
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Stem cell - Pluripotency
<b>Gene Summary:</b>	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions as an autocrine factor which heightens osteoclast formation and bone resorption. This gene has three pseudogenes located on chromosomes 4, 9 and 10, respectively. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. Annexin A2 expression has been found to correlate with resistance to treatment against various cancer forms. [provided by RefSeq, Dec 2019]