

Product datasheet for **RG227337**

SQSTM1 (NM_001142298) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SQSTM1 (NM_001142298) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: SQSTM1
Synonyms: A170; DMRV; FTDALS3; NADGP; OSIL; p60; p62; p62B; PDB3; ZIP3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG227337 representing NM_001142298
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGTCGCTCACCGTGAAGGCCTACCTTCTGGGCAAGGAGGACGCGCGCGGAGATTCGCCGCTTCA
GCTTCTGCTGCAGCCCCGAGCCTGAGGCGGAAGCCGAGGCTGCGGCGGGTCCGGGACCCTGCGAGCGGCT
GCTGAGCCGGGTGGCCGCCCTGTCCCCGCGTGC GGCTGGCGGCTTCCAGGCGCACTACCGCGATGAG
GACGGGGACTTGGTTGCCTTTCCAGTGACGAGGAATTGACAATGGCCATGTCCTACGTGAAGGATGACA
TCTTCCGAATCTACATTAAGAGAAAAAGAGTGCCGGCGGGACCACCGCCACCGTGTGCTCAGGAGGC
GCCCGCAACATGGTGCACCCCAATGTGATCTGCGATGGCTGCAATGGGCTGTGGTAGGAACCCGCTAC
AAGTGCAGCGTCTGCCAGACTACGACTTGTGTAGCGTCTGCGAGGAAAAGGGCTTGACCCGGGGCACA
CCAAGCTCGCATTCCCCAGCCCCCTCGGGCACCTGTCTGAGGGCTTCTCGCACAGCCGCTGGCTCCGGAA
GGTGAAACACGGACACTTCGGGTGGCCAGGATGGGAAATGGGTCCACCAGGAAACTGGAGCCACGTCTCT
CCTCGTGCAGGGGAGGCCCGCCCTGGCCCCACGGCAGAATCAGCTTCTGGTCCATCGGAGGATCCGAGTG
TGAATTTCTGAAGAAGTGGGGAGAGTGTGGCAGCTGCCCTTAGCCCTCTGGGCATTGAAGTTGATAT
CGATGTGGAGCACGGAGGAAAAGAAGCCGCTGACCCCCGTCTCTCCAGAGAGTCCAGCACAGAGGAG
AAGAGCAGCTCACAGCCAAGCAGCTGCTGCTGATCCCAGCAAGCCGGGTGGGAATGTTGAGGGCCGCA
CGCAGTCTTGGCGGAGCAGATGAGAAAGATCGCCTTGGAGTCCGAGGGGCGCCCTGAGGAACAGATGGA
GTCGGATAAAGTTCAGGAGGAGATGATGACTGGACCCATCTGTCTTCAAAGAAGTGGACCCGCTTACA
GGTGAAGTCCAGTCCCTACAGATGCCAGAATCCGAAGGGCCAAGCTCTCTGGACCCTCCAGGAGGGAC
CCACAGGGCTGAAGGAAGCTGCCTTGTACCCACATCTCCCGCCAGAGGCTGACCCGCGGCTGATTGAGTC
CCTCTCCAGATGCTGTCCATGGGCTTCTCTGATGAAGGCGGCTGGCTACCAGGCTCCTGCAGACCAAG
AACTATGACATCGGAGCGGCTCTGGACACCATCCAGTATTCAAAGCATCCCCCGCGTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >RG227337 representing NM_001142298
 Red=Cloning site Green=Tags(s)

MASLTVKAYLLGKEDAAREIRRFSCSPEPEAEAEAAAGGPCERLLSRVAALFPALRPGGFQAHYRDE
 DGDLVAFSSDEELTMAMSYVKDDIFRIYIKEKKECRRDHRPPCAQEAPRNMVHPNVICDGCNGPVVGT
 RYKCSVCPDYDLCSVCEGKGLHRGHTKLAFFSPFGHLSEGFSHSRWLKVKHGHFGWPGWEMGPPGNWSPR
 PRAGEARPGPTAESASGSPEDPSVNFLKNVGESVAAAL SPLGIEVDIDVEHGGKRSRLTPVSPESSSTEE
 KSSSQPSSCCSDPSKPGNVEGATQSLAEQMRKIALESEGRPEEQMESDNCSSGGDDWTHLSSKEVDPST
 GELQSLQMPSEGPSLDPSEQEGPTGLKEAALYPHLPPEADPRLIESLSQMLSMGFSDGGWLTRLLQTK
 NYDIGAALDTIQYSKHPPL

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001142298

ORF Size:	1068 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001142298.1 , NP_001135770.1
RefSeq Size:	2931 bp
RefSeq ORF:	1071 bp
Locus ID:	8878
UniProt ID:	Q13501
Cytogenetics:	5q35.3
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	This gene encodes a multifunctional protein that binds ubiquitin and regulates activation of the nuclear factor kappa-B (NF-kB) signaling pathway. The protein functions as a scaffolding/adaptor protein in concert with TNF receptor-associated factor 6 to mediate activation of NF-kB in response to upstream signals. Alternatively spliced transcript variants encoding either the same or different isoforms have been identified for this gene. Mutations in this gene result in sporadic and familial Paget disease of bone. [provided by RefSeq, Mar 2009]