

Product datasheet for **RG227388**

ZNF566 (NM_001145343) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF566 (NM_001145343) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZNF566
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG227388 representing NM_001145343 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCAGGAGTCAGTGATGTTTCAGTGATGTGTCCGTAGACTTCTCTCAGGAGGAGTGGGAATGCCTGA
ATGATGATCAGAGAGATTTATACAGAGATGTGATGTTGGAGAATTACAGCAACCTGGTTTCAATGGGGCA
TTCTATTTCTAAACCAATGTGATCTCCTACTTGGAGCAAGGGAAGGAGCCCTGGTTGGCTGACAGAGAG
CTAACAAGAGGCCAGTGGCCAGTCTGGAATCAAGATGTGAGACCAAGAAATTTCTGAAGAAAGAAA
TTTATGAAATAGAATCAACCCAGTGGGAAATAATGGAAAACTACAAGACGTGATTTTCAGTGCTCCAG
TTTCAGAGATGATTGGGAATGTAATCGGCAGTTTAAGAAAGAACTCGGCTCTCAGGGGGACATTTCAAT
CAATTGGTATTCATCATGAAGATCTGCCCACTTTGAGTCACCATCCATCCTTCACATTACAGCAAATCA
TTAACAGTAAAAAGAAATTTCTGTGCATCTAAAGAATATAGGAAAACCTTTAGACATGGCTCACAGTTTGC
TACACATGAGATAATTCATACCATTTGAGAAGCCTTATGAATGTAAGGAATGTGGAAAGTCCCTTTAGACAT
CCCTCAAGACTCACTCATCATCAGAAAATTCATACTGGCAAGAAACCTTTGAATGTAAGGAATGTGGAA
AAACCTTTATTTGTGGCTCAGACCTTACTCGACATCACAGAATTCACACTGGTGAGAAACCTATGAATG
TAAGGAATGTGGAAAGCCTTTAGTAGTGGTTCAAACCTCACTCGACATCAGAGAATTCACACAGGTGAG
AAGCCTTATGAATGCAAAGAATGCGGGAAGGCCTTTAGTAGTGGCTCAAACCTTACTCAACATCAGAGAA
TTCACTGGGGAAAAACCTATGAATGTAAGGAATGTGGCAATGCCTTTAGTCAGAGCTACAACCTTAT
TAAACATCAAAGAATCCATACAGGTGAGAAACCTTACGAATGTAAGGAATGTGAAAAGGCTTTTCGTTCT
GGCTCAGACCTTACTAGACATCAGAGAATTCATACTGGGGAGAAAACCTATGAATGTAAGATATGTGGGA
AGGCTTATTTCTCAGAGTTCACAGCTTATTAGTCATCATAGAATTCATACTAGTGAGAAAACCTATGAATA
TAGGGAATGTGGAAAGAACTTTAATTATGACCCACAACCTATTTCAGCATCAAATTTGTACTGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG227388 representing NM_001145343
Red=Cloning site Green=Tags(s)


MAQESVMFSDVSVDFSQEEWECLNDDQRDLYRDVMLENYSNLVSMGHSISKPNVISYLEQGKEPWLADRE
 LTRGQWPVLESRCETKKLFLKKEIYEIESTQWEIMEKL TRRDFQCSSFRRDDWECNRQFKKELGSQGGHFN
 QLVFTHEDLPTLSHHPSTLQQIINSKKKFCASKEYRKTFRHGSQFATHEIIHTIEKPYECKEKGKFSFRH
 PSRLTHHQIHTGKKPFECKEKGKTFICGSDLTRHHRIHTGEKPYECKEKGKAFSSGSNFTRHQRIHTGE
 KPYECKEKGKAFSSGSNFQHQRIHTGEKPYECKEKGNAFSQSSQLIKHQRIHTGEKPYECKEKEKAFRS
 GSDLTRHQRIHTGEKPYECKGKAYSQSSQLISHHRIHTSEKPYEYRECGKNFNYPQLIQHQNLVW

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001145343

ORF Size:	1257 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001145343.1 , NP_001138815.1
RefSeq Size:	5253 bp
RefSeq ORF:	1260 bp
Locus ID:	84924
UniProt ID:	Q969W8
Cytogenetics:	19q13.12
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]