

## Product datasheet for **RG227468**

### ADTRP (NM\_001143948) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ADTRP (NM\_001143948) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** ADTRP  
**Synonyms:** AIG1L; C6orf105; dj413H6.1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG227468 representing NM\_001143948  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACGAAGACTTCTACATGCATATACCACTTCCTTGTCTGAGCTGGTATACTTTCCTCAATTATTACA  
TCTCACAGGAAGGAAAAGACGAGGTGAAACCCAAAATCTGGCAAATGGTGAAGGTGAAATATATGAC  
GCTGCTTAATCTGCTCAAGAACAGGACTGCTGGGTTTGACATCTACCAGCCAGGAAGCTTAGGCAGCTC  
TTGCAGACCATTTCTACGGGTCACCTGCCTGGATGATGTGCTGAAAAGAACCAAGGGGAAAAGACA  
TTAAGTTCCTAACTGCCTCAGAGACCTGCTTTTCACCACTCTGGCTTTTCCTGTATCCAGTTTTGTATT  
TTTGGCATTCTGGATCCTTTTCTACAATCGAGATCTCATTTACCCCAAGGTCTAGATACTGTATC  
CCCGTGTGGCTGAATCATGCAATGCACACTTTTCATATTTCCCATCACATTGGCTGAAGTCGTCTCAGGC  
CTCACTCCTATCCATCAAGAAGACAGGACTCACCTTGCTGGCTGCTGCCAGCATTGCTTACATCAGCCG  
CATCCTATGGCTCTACTTTGAGACGGGTACCTGGGTGATCCTGTGTTTGCCAACTCAGCCTCTTGGGT  
CTAGCAGCTTTCTCTCTCAGCTACGTCTTCATCGCCAGCATCTACCTACTTGGAGAGAAGCTCAACC  
ACTGAAATGGGGTGACATGAGGCAGCCACGGAAGAAGAGGAAG

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG227468 representing NM\_001143948  
 Red=Cloning site Green=Tags(s)

MTKTSTCIYHFLVLSWYTFLNYYISQEGKDEVKPKILANGARWKYMTLLNLLKNRTAGFDIYQPGSFRQL  
 LQTIIFYGVTCLDDVLKRTKGGKDIKFLTAFRDLLFTLLAFPVSTFVFLAFWILFLYNRDLIYPKVLDTVI  
 PVWLNHAMHTFIFPITLAEVVLRPHSYPSKKTGLTLAAASIAYISRILWLYFETGTWVYPVFAKLSLLG  
 LAAFFLSYVFIASIIYLLGEKLNHWKWDMDRQPRKRRK

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM\_001143948

ORF Size: 744 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001143948.1](#), [NP\\_001137420.1](#)

**RefSeq Size:** 1866 bp

**RefSeq ORF:** 747 bp

**Locus ID:** 84830

**UniProt ID:** [Q96IZ2](#)

**Cytogenetics:** 6p24.1

**Protein Families:** Transmembrane

**Gene Summary:** Regulates the expression and the cell-associated anticoagulant activity of the inhibitor TFPI in endothelial cells (in vitro).[UniProtKB/Swiss-Prot Function]