

Product datasheet for **RG227473**

Leukotriene B4 Receptor (LTB4R) (NM_001143919) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Leukotriene B4 Receptor (LTB4R) (NM_001143919) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LTB4R
Synonyms:	BLT1; BLTR; CMKRL1; GPR16; LTB4R1; LTBR1; P2RY7; P2Y7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG227473 representing NM_001143919 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACACTACATCTTCTGCAGCACCCCTCACTAGGTGTAGAGTTCATCTCTCTGCTGGCTATCATCC
TGCTGTCAGTGGCGCTGGCTGTGGGGCTTCCCGCAACAGCTTTGTGGTGTGGAGTATCTGAAAAGGAT
GCAGAAGCGCTCTGCTACTGCCCTGATGGTGTGAACCTGGCCCTGGCCGACCTGGCCGATTGCTCACT
GCTCCCTTTTTCCTTCACTTCTGGCCCAAGGCACCTGGAGTTTTGGACTGGCTGGTGGCCGCTGTGTC
ACTATGTCTGCGGAGTCAGCATGTACGCCAGCGTCCTGCTTATCACGGCCATGAGTCTAGACCGCTCACT
GGCGGTGGCCCGCCCTTTGTGTCCAGAAAGCTACGCACCAAGGCGATGGCCCGGCGGTGCTGGCAGGC
ATCTGGGTGTTGTCTTTCTGCTGGCCACACCCGTCCTCGCTACCGCACAGTAGTGCCTGGAAAACGA
ACATGAGCCTGTGCTTCCCGCGGTACCCAGCGAAGGGCACCGGGCCTTCCATCTAATCTTCGAGGCTGT
CACGGGCTTCTGCTGCCCTTCTGGCTGTGGTGGCCAGCTACTCGGACATAGGGCGTGGCTACAGGCC
CGGCGCTTCCGCCGAGCCGCCGACCGGCCGCTGGTGGTGTGCTCATCATCTGACCTTCGCCGCTTCT
GGCTGCCCTACCACGTGGTGAACCTGGCTGAGGCGGGCCGCGCTGGCCGGCCAGGCCCGCGGTTAGG
GCTCGTGGGAAGCGGCTGAGCCTGGCCGCAACGTGCTCATCGCACTCGCTTCTGAGCAGCAGCGTG
AACCCCGTGTACGCGTGCCTGGCCGCGGCGCTGCTGCGCTCGGCGGGCGTGGGCTTCGTCGCCAAGC
TGCTGGAGGGCACGGGCTCCGAGGCGTCCAGCACGCGCCGCGGGGAGCCTGGGCCAGACCGCTAGGAG
CGGCCCGCGCTTTGGAGCCCGCCCTTCCGAGAGCCTCACTGCCTCCAGCCCTTCAAGTTAAACGAA
CTGAAC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG227473 representing NM_001143919
 Red=Cloning site Green=Tags(s)

MNTTSSAAPPSLGVFISLLAIILL SVALAVGLPGNSFVWVSILKRMQKRSVTALMVLNLALADLAVLLT
 APFFLHFLAQGTWSFGLAGCRLCHYVCGVSMYASVLLITAMSLDRSLAVARPFVSQLRRTKAMARRVLAG
 IWVLSFLLATPVLAYRTVVPWKTNMSLCFPRYPSEGHAFHLIFEAVTGFLLPFLAVVASYSDIGRRLQA
 RRFRRSRRTGRLVVLIIILTFAAFWLPYHVVNLAEAGRALAGQAAGLGLVGKRLSLARNVLIALAFLLSSSV
 NPVLYACAGGGLLR SAGVGFVAKLLEGTGSEASSTRRGGSLGQTARSGPAALEPGPSESLTASSPFKLINE
 LN


SGPTRRRLE - GFP Tag - V

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001143919

ORF Size:	1056 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001143919.1 , NP_001137391.1
RefSeq Size:	2525 bp
RefSeq ORF:	1059 bp
Locus ID:	1241
UniProt ID:	Q15722
Cytogenetics:	14q12
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction
Gene Summary:	Receptor for extracellular ATP > UTP and ADP. The activity of this receptor is mediated by G proteins which activate a phosphatidylinositol-calcium second messenger system. May be the cardiac P2Y receptor involved in the regulation of cardiac muscle contraction through modulation of L-type calcium currents. Is a receptor for leukotriene B4, a potent chemoattractant involved in inflammation and immune response.[UniProtKB/Swiss-Prot Function]