

## Product datasheet for **RG227565**

### NSMAF (NM\_001144772) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	NSMAF (NM_001144772) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NSMAF
Synonyms:	FAN; GRAMD5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG227565 representing NM\_001144772  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGACGGGCCAGGCCAGATTCCACCGCGGAGCCGAGGCATGCCGAGCTATCCGGCAGTGCAGGTCCCGGC  
 GGCCGCGGGCAGCTCACCCGGAAGGCCGCGCGGGTCCGGGCGTGGATCCCGGCACGCGAGCCCGGGCGT  
 TCGCAGGGGAGGATTTTCTTGTGCTGCTTAACCTGGAGGAGTACTACTTTGAACAGCATAGAGCCAAT  
 CACATTTTGCACAAGGGCAGTCACCATGAAAGGAAAATCAGAGGCTCCTTAAAAATATGTTCAAAAATCGG  
 TGATTTTTGAACCAGATTCAATATCCCAGCCCATCATCAAGATTCCTTTGAGAGACTGTAAAAATAGG  
 AAAGCATGGAGAAAATGGAGCCAATAGACACTTCACAAAGGCAAAAATCTGGGGTATTTCACTCATTTTC  
 AGTCAGGTATATTTCAATAAAGAACATAATGTTGTTGCACCATATAAAATAGAAAGGGGCAAAAATGGAAT  
 ATGTTTTTGAATTGGATGTTCCCGGAAAGTGAAGATGTTGTGGAGACGTTGCTTCAGCTTCACAGAGC  
 ATCCTGCCTTGACAAATGGGTGACCAACCGCCATGATAACAGCTATTTTGCAGTCTGTTTAGCTAGA  
 ACATCATTTGACAAAAACAGGTTCCAAAACATTTCTGAAAAGCTGCACATGGAATGCAAAGCAGAAATGG  
 TGACGCCTCTGGTACTAATCCTGGACAGTGTGCATCACGGACACAAAACCTGTATTTTACGCCCCTCAA  
 CGGCTACCCGAAACCTGTGGTCCAGATAACACTCCAAGATGTCCGCCCATCTACAAAAGGAGGCACGGC  
 CTCATGCCTCTGGGCTTGGAAAGTATTTTGCACAGAAGATGATCTGTGTTCCGACATCTACCTAAAGTTCT  
 ATGAACCTCAAGATAGAGATGATCTCTATTTTTACATTGCCACATACCTAGAGCACCATGTGGCGGAGCA  
 CACTGTGAGAGCTACATGCTGCAGTGGCAGCGTGGACACCTTCCAACATCAGTACCTCCTTCACCTC  
 AACAACTGGCCGACCGCAGCTGCAACGACCTCTCCAGTACCCTGTGTTCCATGGATAATACATGATT  
 ATTCAGCTCAGAAGTATGTTTCAAAATCCAGGAACCTTCCGGGATCTCAGTAAGCCAGTGGGACCT  
 AAATAAGGAACGGCTGGAGAGACTACTGACACGCTACCAGGAAATGCCTGAACCAAAGTTATGTATGGG  
 AGTCACTACTCTTCCCGGGTTATGTACTTTTTATCTTGTAGGATTGCACCAGAGTATATGCTGTGCC  
 TGCAGAATGGAAGATTTGATAATGCAGATAGAATGTTCAACAGTATTGCAGAAAATTTGAAAAACTGTCT  
 GGATGGTGAACGGATTTTAAAGAGTTAATCCAGAATCTATGGTGATGATGTGAGCTTCTAGTCAAT  
 AGCCTGAAGTTGGATTTGGGAAAGAGACAAGGAGGACAGATGGTTGACGACGTGGAGCTTCCCCCTTGGG  
 CTTCCAGTCCCGAGGACTTCTCCAGAAGAGCAAAGATGCATTGGAAAGCAATTATGTGCTGAACACCT  
 TCACGAGTGGATTGATCTAATATTTGGCTACAAAACAAAAGGGAGTATGATGAGTGGGGCCATAATGTA  
 TTTTATCCCTGACCTATGAAGGAGGTGTAGACTTGAACAGCATCCAGGATCCTGATGAGAAGGTAGCCA  
 TGCTTACGCAAAATCTTGAATTTGGGCAGACACCAAAAACAACTATTTGTGACACCACATCCTCGAAGGAT  
 CACCCCAAAGTTTAAAAGTTTGTCCAGACCTCCAGTTATAATGCTTCTATGGCAGATTTCCAGGTGAA  
 GAGTCTTTTGAAGACCTGACCGAAGAAAGCAAAACACTGGCCTGGAATAACATCACCAAACCTGCAGTTAC  
 ACGAGCACTATAAAATCCACAAAGAAGCAGTACTGGAATCACGGTCTCTCGCAATGGATCTTCAGTATT  
 CACAACATCCCAAGATCCACCTTGAAGATGTTTTCTAAAGAATCAAAAATGCTACAAAGAAGTATATCA  
 TTTTCAAATATGGCTTTATCGTCTGTTACTTTTACCAGGAGATGCCACTGTCAAACTTCTTCATGGG  
 ATAATAATGTCTATTTTTATTCCATAGCATTGGAAGACGCCAGGACACGTTAATGGGACATGATGATGC  
 TGTTAGTAAGATCTGTTGGCATGACAACAGGCTATATTCTGCATCGTGGGACTTACAGTGAAGGTGTGG  
 TCTGGTGTCTCTGCAGAGATGCCAGGCACCAAAAGACACCCTTTGACTTGTGGCCGAGCTGGAACATG  
 ATGTCAGTGTAGATACAATCAGTTTAAATGCTGCAAGCACAAGTGTAGTTTCCGGCACCAAGAAGGCAC  
 AGTGAATATTTGGGACCTCACAACGGCCACCTTAATGCACCAGATTCCATGCCATTACAGGGATTGTATGT  
 GACTGCTTTTACCCAGATAGTCGCCATGTCCTCAGCACAGGAACAGATGGCTGTCTTAATGTCATTG  
 ATGTGCAGACAGGAATGCTCATCTCCTCCATGACATCAGATGAGCCCCAGAGGTGCTTTGTCTGGGATGG  
 AAATCCGTTTTATCTGGCAGTCACTGTTGAACTGCTCGTTTGGGACCTCCTTGGAGCAAAAATCAGT  
 GAGAGAATACAGGGCCACACAGGTGCTGTGACATGTATGGATGAATGAACAGTGTAGCAGTATCATCA  
 CAGGAGGGGAAGACAGACAAATTATATTCTGGAATTGCAGTAT

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

**Protein Sequence:** >RG227565 representing NM\_001144772  
 Red=Cloning site Green=Tags(s)

```

MTGQARFHRGAEACRAIRQCRSRPRAAHPEGRAGSGRGRSHASPGVRRGGFSLLLLNL E E Y F E Q H R A N
H I L H K G S H H E R K I R G S L K I C S K S V I F E P D S I S Q P I I K I P L R D C I K I G K H G E N G A N R H F T K A K S G G I S L I F
S Q V Y F I K E H N V V A P Y K I E R G K M E Y V F E L D V P G K V E D V V E T L L Q L H R A S C L D K L G D Q T A M I T A I L Q S R L A R
T S F D K N R F Q N I S E K L H M E C K A E M V T P L V T N P G H V C I T D T N L Y F Q P L N G Y P K P V V Q I T L Q D V R R I Y K R R H G
L M P L G L E V F C T E D D L C S D I Y L K F Y E P Q D R D D L Y F Y I A T Y L E H H V A E H T A E S Y M L Q W Q R G H L S N Y Q Y L L H L
N N L A D R S C N D L S Q Y P V F P W I I H D Y S S E L D L S N P G T F R D L S K P V G A L N K E R L E R L L T R Y Q E M P E P K F M Y G
S H Y S S P G Y V L F Y L V R I A P E Y M L C L Q N G R F D N A D R M F N S I A E T W K N C L D G A T D F K E L I P E F Y G D D V S F L V N
S L K L D L G K R Q G G Q M V D D V E L P P W A S S P E D F L Q K S K D A L E S N Y V S E H L H E W I D L I F G Y K Q K G S D A V G A H N V
F H P L T Y E G G V D L N S I Q D P D E K V A M L T Q I L E F G Q T P K Q L F V T P H P R R I T P K F K S L S Q T S S Y N A S M A D S P G E
E S F E D L T E E S K T L A W N N I T K L Q L H E H Y K I H K E A V T G I T V S R N G S S V F T T S Q D S T L K M F S K E S K M L Q R S I S
F S N M A L S S C L L L P G D A T V I T S S W D N N V Y F Y S I A F G R R Q D T L M G H D D A V S K I C W H D N R L Y S A S W D S T V K W W
S G V P A E M P G T K R H H F D L L A E L E H D V S V D T I S L N A A S T L L V S G T K E G T V N I W D L T T A T L M H Q I P C H S G I V C
D T A F S P D S R H V L S T G T D G C L N V I D V Q T G M L I S S M T S D E P Q R C F V W D G N S V L S G S Q S G E L L V W D L L G A K I S
E R I Q G H T G A V T C I W M N E Q C S S I I T G G E D R Q I I F W K L Q Y
  
```

TRTRPLE – GFP Tag – V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**




**ACCN:** NM\_001144772

**ORF Size:** 2844 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001144772.1</a> , <a href="#">NP_001138244.1</a>
<b>RefSeq Size:</b>	3522 bp
<b>RefSeq ORF:</b>	2847 bp
<b>Locus ID:</b>	8439
<b>UniProt ID:</b>	<a href="#">Q92636</a>
<b>Cytogenetics:</b>	8q12.1
<b>Protein Families:</b>	Druggable Genome
<b>Gene Summary:</b>	This gene encodes a WD-repeat protein that binds the cytoplasmic sphingomyelinase activation domain of the 55kD tumor necrosis factor receptor. This protein is required for TNF-mediated activation of neutral sphingomyelinase and may play a role in regulating TNF-induced cellular responses such as inflammation. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Jan 2009]

Product images:



Circular map for RG227565