

## Product datasheet for **RG227650**

### **BDNF (NM\_001143805) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	BDNF (NM_001143805) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	BDNF
Synonyms:	ANON2; BULN2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG227650 representing NM_001143805 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACCATCCTTTTCCTTACTATGGTTATTTCACTTTGGTTGCATGAAGGCTGCCCCATGAAAGAAG  
CAAACATCCGAGGACAAGGTGGCTTGGCTACCCAGGTGTGCGGACCCATGGGACTCTGGAGAGCGTGAA  
TGGGCCAAGGCAGGTTCAAGAGGCTTGACATCATTGGCTGACACTTTGAAACACATGATAGAAGAGCTG  
TTGGATGAGGACCAGAAAGTTCGGCCCAATGAAGAAAACAATAAGGACGCAGACTTGTACACGTCCAGGG  
TGATGCTCAGTAGTCAAGTGCCTTTGGAGCCTCCTCTTTCTTTCTGCTGGAGGAATACAAAATTACCT  
AGACGCTGCAAACATGTCCATGAGGGTCCGGCGCCACTCTGACCCGCCCCGAGGGGAGCTGAGCGTG  
TGTGACAGTATTAGTGAGTGGTAACGGCGGCAGACAAAAGACTGCAGTGGACATGTGGGGCGGGACGG  
TCACAGTCCTTGAAGAGGTCCCTGTATCAAAAGGCCAACTGAAGCAATACTTCTACGAGACCAAGTGCAA  
TCCCATGGGTTACACAAAAGAAGGCTGCAGGGCATAGACAAAAGGCATTGGAATCCCAGTGCCGAACT  
ACCCAGTCGTACGTGCGGGCCCTTACCATGGATAGCAAAAAGAGAATTGGCTGGCGATTACATAAGGATAG  
ACACTTCTGTGTATGTACATTGACCATTAAGGGGAAGA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


**Protein Sequence:** >RG227650 representing NM\_001143805  
 Red=Cloning site Green=Tags(s)

MTILFLTMVISYFGCMKAAPMKEANIRGQGGLAYPGVTRHTGTLESVNGPKAGSRGLTSLADTFEHMIEEL  
 LDEDQKVRPNEENNKDADLYTSRVMLSSQVPLEPPLLFLLEEYKNYLDAAANMSMRVRRHSDPARRGELSV  
 CDSISEWVTAADKKTAVDMSSGGT VTVLEKVPVSKGQLKQYFYETKCNPMGYTKEGCRGIDKRHWNSQCRT  
 TQSYVRALTMDSKKRIGWRFIRIDTSCVCTLTIKRGR

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001143805

**ORF Size:** 741 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001143805.1</a> , <a href="#">NP_001137277.1</a>
<b>RefSeq Size:</b>	3827 bp
<b>RefSeq ORF:</b>	744 bp
<b>Locus ID:</b>	627
<b>UniProt ID:</b>	<a href="#">P23560</a>
<b>Cytogenetics:</b>	11p14.1
<b>Protein Families:</b>	Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Secreted Protein, Transmembrane
<b>Protein Pathways:</b>	Huntington's disease, MAPK signaling pathway, Neurotrophin signaling pathway
<b>Gene Summary:</b>	This gene encodes a member of the nerve growth factor family of proteins. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate the mature protein. Binding of this protein to its cognate receptor promotes neuronal survival in the adult brain. Expression of this gene is reduced in Alzheimer's, Parkinson's, and Huntington's disease patients. This gene may play a role in the regulation of the stress response and in the biology of mood disorders. [provided by RefSeq, Nov 2015]

## Product images:



Circular map for RG227650