

Product datasheet for **RG227887**

SH2B1 (NM_001145796) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SH2B1 (NM_001145796) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SH2B1
Synonyms:	PSM; SH2B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG227887 representing NM_001145796
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAATGGTGCCCTTCCCAGAGGACGGGGCCTCCCCCTCGTCTCCCCCGCTGCCCCACCCCGCCCC
 CTAGTTGGCGGGAGTTCTGTGAGTCCCACGCCGGCTGCGGCTCTGGACTTTGCCCGCCGTTTTCGCCT
 CTACCTGGCCTCCCACCCCAATATGCGGGGCCGGGGCCGAGGCTGCCTTCTCCCGCCGTTTTGTGAG
 CTCTTCTGCAGCACTTTGAAGCCGAGGTGGCCCGGGCCTCTGGCTCCCTGTGCGCACCCATCTGGCTC
 CCCTGAGCCCTGGTGGGAGATTTCGCCACATGACCTGTCCCTTGAGAGCTGCAGGTGGTGGGCCCT
 GGCTGTGTGGGCCCTTCTCGATCATCTGAGGACCTGGCCGGCCCCCTCCCTTCTCAGTCTTCTCTCC
 TCTACAACCTCTCAAGCCGAAGCTCAAGAAGCGTTTTCCCTGCGTTCAGTGGTTCGCTGTCCGAG
 GCTCAGTCCGTGGCATCTCGAGTGGCGGGGACCGTTGACCCTCCCTCCGCTGGGCCCTGGAGAC
 CTGCTCAGGCCCCAGTCTTAGGTGAAACAGCAACTCCAACCTCTGGCGGGCTGGGACCGTTGGT
 AGGGGACTGGTCACTGATGGAACGTCCCCTGGGAAAGATGGACTCACCGTTTTGAGAGGCTGAGACTCA
 GTCGGGGAGGGGGCGCCTTGAAGGATGGAGCAGGGATGGTGCAGAGGGAAGAGCTGCTGAGTTTCATGGG
 GGCTGAGGAGGCAGCCCCTGACCCAGCCGAGTGGGCCGGGGAGGAGGGGTGGCTGGGCCTCCTTACGGG
 GGAGGAGGGCAGCCCTCAGTGGCAGAAGTGTGCGCTGCTGCTTGAAGTGAAGGAGAAGGAGGAGGAGAA
 GTCGCCTGGAGTTCTTTGTACCACCCAAGGCCTCTCGGCCCGACTCAGCATCCCCTGCTTCTATCAC
 AGAGCTCCGACAACCACAGCCCTGGAGATGCCTGACCGGGAGAACACGTTTGTGGTTAAGGTGGAAGGT
 CCATCCGAGTATCATGGAGACAGTGGATGCCAGCATGTGAAGGCCTGGTGTCTGACATCCAAGAAT
 GCCTGAGCCCAGGACCCCTGCCCTGCTACCAGTCCCCGCCCATGACCCCTCCCTCTGGCCCTGGGACCTC
 ATTCTTTACAAGGGAGAACACAGACAGCCTGGAGCTGCCTGCCTGAATCACTCGGAGAGTCTACCCAGC
 CAGGACCTGCTGCTTGGACCCAGCGAGCAATGACCGCCTGTCGAGGGGGCATATGGGGCCTCTCAG
 ACCGCCCTCGGCATCCATCTCCCCAGCTCTGCCTCCATTGCCGCTCCCATTTTACTCGATGGAACCT
 GCTTCCCCCAGAGTTGCCCCCCGCATCCCCATTGAAGAGGGACCCCAACAGGGACAGTTCATCCCCTC
 TCAGCCCCCTACCCTCCCTTGGACTCCGAAACAGCCACAGGGTCTTCTGTCCAGGGGGAGCCAG
 AGGGCGGTGAGGGGGACAGCCCCTCTCAGGGTATCCTTGGTTCCACGGGATGCTCTCTCGGCTCAAGGC
 TGCACAGTTGGTGTGACTGGCGGCACTGGCTCCACGGTGTCTTCTGGTGCGCCAGAGTGAACAAGG
 CGGGGTGAATACGTCTCACCTTCAACTCCAGGGCAAGGCAAGCACCTGCGTTTGTGCTGAACGAGG
 AGGGTCAGTGCCGGTCCAGCACCTGTGGTTCCAGTCCATTTTCGATATGCTCGAGCACTCCGGGTGCA
 CCCCATCCCTTTGGAGTCGGGAGGCTCCAGTGTGTTGCTTGTGCTGCTGCTGCTGCTGCTGCTGCTGCT
 CAGCAGGGCCGGGAGCAGGCTGGGAGCCATGCGGGGGTGTGCGAGGGAGATGGATGCCACCCCGATGCC
 CCTGCACCCTCATGCCCTTCGGAGCGAGTACTGTGTAACCGACCACCTCCCA

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG227887 representing NM_001145796
 Red=Cloning site Green=Tags(s)

```
MNGAPSPEDGASPSSPPLPPPPPPSWREFCESHARAAALDFARRFRLYLASHPQYAGPGAEEAFSRRFAE
LFLQHFEAEVARASGSLSPILAPLSPGAEISPHDLSLESCRVGGLAVLGPSSRSEDLAGLPPSSVSSS
STTSSKPKLKKRFSLRVSVGRSVRGSVRGILQWRGTVDPPSSAGPLETSSGPPVLGGNSNSNSGGGAGTVG
RGLVSDGTPSPGERWTHRFRERLRLSRGGALKDGAGMVQREELL SFGMAEEAAPDPAGVGRGGVAGPPSG
GGGQPQWQKCRLLLRSEGEgggsrLEFFVPPKASRPRLSIPCSSITDVRTTTALEMPDRENTFVVKVEG
PSEYIMETVDAQHVKAWSIDIQECLSPGPCPATSPRPMTLPLAPGTSFLTRENTDSLELSCLNHSESLPS
QDLLLLGPSESNDRLSQGAYGGLSDRPSASISPPSSASIAASHFDSMELLPELPPRIPIIEEGPPTGTVHPL
SAPYPLDTPETATGSFLFQGEPEGGEDQPLSGYPWFHGMRLKAAQLVLTGGTGSHGVLVVRQSETR
RGEYVLTNFQGKAKHLRLSLNEEGQCRVQHLWFQSFIDMLEHFRVHPIPLESGGSSDVVLVSYVPPSSQR
QQGREQAGSHAGVCEGDGCHPDASCTLMFPGASDCVTDHLP
```

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001145796

ORF Size: 2013 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145796.1](#), [NP_001139268.1](#)

RefSeq Size: 3046 bp

RefSeq ORF: 2016 bp

Locus ID: 25970


UniProt ID: [Q9NRF2](#)

Cytogenetics: 16p11.2

Protein Pathways: Neurotrophin signaling pathway

Gene Summary: This gene encodes a member of the SH2-domain containing mediators family. The encoded protein mediates activation of various kinases and may function in cytokine and growth factor receptor signaling and cellular transformation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009]

Product images:



Circular map for RG227887