

Product datasheet for **RG228063**

DAOA (NM_001161812) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: DAOA (NM_001161812) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: DAOA
Synonyms: LG72; SG72
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG228063 representing NM_001161812
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGACGGCTATTTGGAAATGGCACAGAGGCATTTACAGAGATCATTATGTCCTTGGGTCTCTTACCT
TCCTCAGCCCTATGCAGAGAGCTGGGACTTCTAACCAATGGAACATGGGAAAGGTGATGACACTTGACTC
CGGTGATGAGGTTACCTTTCACAAGACCTGGCCAACCTGAGACCGGATCTCCTTACAGGCTGAAGAAGTA
AGCAGCCATGTTGAAAAGTCTTCATGGCAAGAACTATGAGTTCCTTGCCTATGAGGCCTTAAGGACC
GCAGGCAGCCTCTAGAACGAATGTGGACCTGCAACTACAACCAGCAAAAAGACCAGTCATGCAACCACAA
GGAAATAACTTCTACCAAGCTGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG228063 representing NM_001161812
Red=Cloning site Green=Tags(s)
MRTAIWKWHRGIYRDHYVLGSLTFLSPMQRAGTSNQWNMGKVMTLDSGDEVTFHKTWPTETGSPYRLEEV
SSHVGKVFMARYEFLAYEASKDRRQPLERMWTCNYNQQKDQSCNHKEITSTKAE

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001161812.1 , NP_001155284.1
RefSeq Size:	924 bp
RefSeq ORF:	378 bp
Locus ID:	267012
UniProt ID:	P59103
Cytogenetics:	13q34
Protein Families:	Druggable Genome
Gene Summary:	This gene encodes a protein that may function as an activator of D-amino acid oxidase, which degrades the gliotransmitter D-serine, a potent activator of N-methyl-D-aspartate (NMDA) type glutamate receptors. Studies also suggest that one encoded isoform may play a role in mitochondrial function and dendritic arborization. Polymorphisms in this gene have been implicated in susceptibility to schizophrenia and bipolar affective disorder. Alternatively spliced transcript variants encoding different isoforms have been identified.[provided by RefSeq, Mar 2011]