

Product datasheet for **RG228064**

zinc finger protein 138 (ZNF138) (NM_001160183) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: zinc finger protein 138 (ZNF138) (NM_001160183) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ZNF138
Synonyms: pHZ-32
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG228064 representing NM_001160183
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGACCACTGACATTTATGGATGTGGCCATAGAGTTCTTTGGAGGAGTGCCAGTGCCTGGACACTG
CACAGCGGAATGTATATAGGCATGTGATGTTAGAGAACTACAGAAACCTGGTTTTCTGGACCTGATTAC
CTGTCTGGAACAAGGAAAAGAGCCCTGGAATATGAAGAGACATGAGATGGTGGTAGCCAAACATTCAGCT
GTAGCATTTACCTTTGATCTCAGTGGACTCAACATTTTGTTCATTCAAAGTATCCTGCCATTTGTTTCAG
CACTTTATGTCATGGGAAACAGAAACCAGTGTCTGCAAAAAGCACCTAGAAGCCAGAAATAAAGATCTCTG
TGTTCTCGTTTTGCCCAAGACCTTTGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG228064 representing NM_001160183
Red=Cloning site Green=Tags(s)

MGPLTFMDVAIEFSLEEWQCLDTAQRNVYRHMVLENYRNLVFLDLITCLEQKPEWMMKRHEMVAKHSA
VAFTFDLSGLNLSFSQILPFVSYALYVMGNRNQCLQKHLEARNKDLVLPKTFG

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001160183

ORF Size: 378 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001160183.3
RefSeq Size:	2719 bp
RefSeq ORF:	381 bp
Locus ID:	7697
UniProt ID:	P52744
Cytogenetics:	7q11.21
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation as a repressor.[UniProtKB/Swiss-Prot Function]