

Product datasheet for **RG228083**

ZDHHC15 (NM_001146257) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZDHHC15 (NM_001146257) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ZDHHC15
Synonyms: DHHC15; MRX91
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG228083 representing NM_001146257
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCGCGAGGCTGGAAGATGGCTCTGTCTGGGGGCTGCGGTGCTGCCCGGGTACTGTCCTGGGTGC
 CAGTGCTCGTTATTGTCCTCGTCGTCTGGTCTACTATGCCTACGCTTTGAACTCTGCCTGGTTAT
 TTACCTCATACTCTACCATGCCATCTTTGTGTTCTTTACCTGGACCTACTGGAAGTCTATCTTTACTC
 CCACAGCAGCAAACCAGAAGTCCACTTGTCTACACAGACAAGGAGCGCTATGAAAATGAAGAAAGAC
 CTGAGGTCCAGAAGCAGATGCTTGTGATATGGCCAAAAGCTACCGGTTACACAAGAAGTGAAGTGG
 AGGTCAAGTTCATCAAAGGCAGCTAGAGAGGCAGCTCAGCAAGTATCTCAGAAAAGCTAAGTCATATATG
 TTCTCAAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG228083 representing NM_001146257
 Red=Cloning site Green=Tags(s)

MRRGWKMALSGGLRCCRRVLSWVPVLVIVLVLWSYYAYVFELCLVIYLILYHAIFVFFTWYWKSIFTL
 PQQPQKQFHLSYTDKERYENEERPEVQKQMLVDMAKKLPVYTRTGSQQFIQRQLERQLSKYLKAKSYM
 FSN

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



Cloning Scheme:

Plasmid Map:


ACCN: NM_001146257

ORF Size: 429 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001146257.2</u>
RefSeq Size:	1312 bp
RefSeq ORF:	432 bp
Locus ID:	158866
UniProt ID:	<u>Q96MV8</u>
Cytogenetics:	Xq13.3
Protein Families:	Transmembrane
Gene Summary:	The protein encoded by this gene belongs to the DHHC palmitoyltransferase family. Mutations in this gene are associated with mental retardatio X-linked type 91 (MRX91). Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]