

Product datasheet for **RG228109**

RPAIN (NM_001160244) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPAIN (NM_001160244) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: RPAIN
Synonyms: HRIP; RIP
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG228109 representing NM_001160244
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGAGTCGTTGAGGTCTCCGCGCGCTCCCTGTACAACTGGTGGGCTCGCCGCTTGAAAGAGG
CTTCCGGCAGAGATGCCTGGAGAGAATGAGAAACAGCCGGGACAGGCTCCTAACAGGTACCGCCAGGC
TGAAGCAGTGGCCAGGAATTCTCAGAACAGCTTTCTAGTTCAAGAGGTGATGGAAGAAGAGTGAAT
GCTTTCAGTCAGTGGAGAATTGTCCAGAAGACTTGGCTCAGCTGGAGGAGCTGATAGACATGGCTGTGC
TGGAGGAAATTCAACAGGAGCTGATCAACCAAGAGCAGTCCATCATCAGCGAGTATGAGAAGAGCTTGCA
GTTTGATGAAAAGTGTCTCAGCATCATGCTGGCTGAGTGGGAGGCAAACCACTCATCTGTCTGTATGT
ACAAAGTACAACCTGAGAATCACAAGCGGTGTGGTGTGTGTCAGTGTGGCCTGTCCATCCCATCTCATG
CCTGTGATACTTGGGCTGTGATCCTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG228109 representing NM_001160244
Red=Cloning site Green=Tags(s)

MAESLRSPRRSLYKLVGSPPWKEAFRQRCLERMNSRDRLNRYRQAGSSGPGNSQNSFLVQEVMEEEWN
ALQSVENCPEDLAQLEELIDMAVLEEIQQELINQEQSIISEYEKSLQFDEKCLSIMLAEWEANPLICPVC
TKYNLRITSGVVVCQGLSIPSHACDTWAVIL

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001160244

ORF Size: 516 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001160244.1, NP_001153716.1</u>
RefSeq Size:	1538 bp
RefSeq ORF:	519 bp
Locus ID:	84268
UniProt ID:	<u>Q86UA6</u>
Cytogenetics:	17p13.2
Gene Summary:	Mediates the import of RPA complex into the nucleus, possibly via some interaction with importin beta. Isoform 2 is sumoylated and mediates the localization of RPA complex into the PML body of the nucleus, thereby participating in RPA function in DNA metabolism. [UniProtKB/Swiss-Prot Function]