

Product datasheet for **RG228173**

Ataxin 3 (ATXN3) (NM_001164779) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ataxin 3 (ATXN3) (NM_001164779) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ATXN3
Synonyms:	AT3; ATX3; JOS; MJD; MJD1; SCA3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG228173 representing NM_001164779 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTCCATCTCCACGAGAAATGGTTAACTTGAATTCTCTCTTGACGGGTCCAGAATTAATATCAG
ATACATATCTTGCACTTTTCTGGCTCAATTACAACAGGAAGTTATTCTATATTTGTCGTTAAGGGTGA
TCTGCCAGATTGCGAAGCTGACCAACTCCTGCAGATGATTAGGGTCCAACAGATGCATCGACCAAACTT
ATTGGAGAAGAATTAGCACAACTAAAAGAGCAAAGAGTCCATAAAACAGACCTGGAACGAGTGTTAGAAG
CAAATGATGGCTCAGGAATGTTAGACGAAGATGAGGAGGATTTGCAGAGGGCTCTGGCACTAAGTCGCCA
AGAAATTGACATGGAAGATGAGGAAGCAGATCTCCGAGGGCTATTGAGCTAAGTATGCAAGGTAGTTCC
AGAAACATATCTCAAGATATGACACAGACATCAGGTACAATCTTACTTCAGAAGAGCTTCGGAAGAGAC
GAGAAGCCTACTTTGAAAAACAGCAGCAAAAGCAGCAACAGCAGCAGCAGCAGCAGCAGCAGCGGGGACCT
ATCAGGACAGAGTTCACATCCATGTGAAAGGCCAGCCACCAGTTCAGGAGCACTTGGGAGTGATCTAGGT
GATGCTATGAGTGAAGAAGACATGCTTCAGGCAGCTGTGACCATGTCTTTAGAACTGTCAGAAATGATT
TGAAAACAGAAGGAAAAAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG228173 representing NM_001164779
Red=Cloning site Green=Tags(s)

MESIFHEKWFNLSLLTGPELISDTYLALFLAQLQQEGYSIFVVKGDLPDCEADQLQMIRVQQMHRPKL
 IGEELAQLKEQRVHKTDLERVLEANDGSGMLDEDEEDLQRALALSRQEIDMEDEEADLRRAIQLSMQGSS
 RNISQDMTQTSNTLSEELRKRREAYFEKQQKQQQQQQQQQDLSGQSSHP CERPATSSGALGSDLG
 DAMSEEDMLQAAVTMSLETVRNDLKTEGKK

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001164779

ORF Size: 720 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001164779.2
RefSeq Size:	6560 bp
RefSeq ORF:	723 bp
Locus ID:	4287
Cytogenetics:	14q32.12
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	Machado-Joseph disease, also known as spinocerebellar ataxia-3, is an autosomal dominant neurologic disorder. The protein encoded by this gene contains (CAG) <i>n</i> repeats in the coding region, and the expansion of these repeats from the normal 12-44 to 52-86 is one cause of Machado-Joseph disease. There is a negative correlation between the age of onset and CAG repeat numbers. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2016]