

Product datasheet for **RG228198**

RDM1 (NM_001163130) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RDM1 (NM_001163130) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	RDM1
Synonyms:	RAD52B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG228198 representing NM_001163130 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGCATTTACTTGTCCACCCCCACAGCATTCTCTGTTACAGCATTCTCTCAGTTTGGCCTTCTGTATT
 CAGTCCGGGTCTTCCCAAATGCTGCAGTGGCCCATCTGGTTTCTATGCCGTCATTAAGTTTATTCTGC
 AAGGGCTGCCACAGAGCCCAAAGGCATGCGACCGGAAGCAGCTTTTTCAGAAATCTCCAGTCAAGGTT
 CGTCTTGGCACCAGACATAAGGCAGTTCAACATCAAGCCCTTGCCCTGAACAGTTCCAAATGCCAAGAAC
 TGGCGAATTACTACTTTGGTTTCAATGGGTGTTCCAAAGGATCATCAAGCTTCAGGAGCTTTCTGACCT
 TGAAGAAAGGGAAATGAAGATAGCATGGTGCCACTTCCGAAGCAAAGCCTGAAGTTCTTCTGTGCTTTA
 GAAGTGGTGTGCCATCCTGTGATTGCAGGAGTCTGGCATTGGCTTGGTGGAGGAGCCTATGGATAAGG
 TGGAGGAAGGACCATTATCATTCTTATGAAAAGGAAGACAGCCAGAAGCTTGCTATTCAGAAGGCTTT
 GTCAGATGCATTCCAGAACTGTTGATTGTTGTTCTAGAAAGTGGTAAATAGCTGTGGAGTACAGACCC
 AGTGAAGACATCGTAGGTGTCAGATGCGAAGAAGAACTACACGGTTTAATTCAAGTCCCTTGCTCTCCCT
 GGAAGCAGTATGGCCAAGAGGAGGAAGGGTATCTCTCGGATTTACGCTTGGAGGAGGAAGAGTTCAAGCT
 GCCAGAACTTGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


Protein Sequence: >RG228198 representing NM_001163130
 Red=Cloning site Green=Tags(s)

MHLLVPPPHSLFTAFSQFGLLYSVRVFPNAVAHPGFYAVIKFYSAARAHRAQKACDRKQLFQKSPVKV
 RLGTRHKAVQHQAALNSSKQELANYFNGFNGCSKRIKLQELSDLEERENEDSMVPLPKQSLKFFCAL
 EVVLPSCDCRSPGIGLVEEPMDKVEEGPLSFLMKRKTAKLAIQKALSDAFQKLLIVVLESGKIAVEYRP
 SEDIVGVRCEEELHGLIQVPCSPWKQYQGEEGYLSDFSLEEEEFRLPELD

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001163130

ORF Size: 783 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163130.1](#), [NP_001156602.1](#)

RefSeq Size: 1065 bp

RefSeq ORF: 786 bp


Locus ID: 201299

UniProt ID: [Q8NG50](#)

Cytogenetics: 17q12

Gene Summary: This gene encodes a protein involved in the cellular response to cisplatin, a drug commonly used in chemotherapy. The protein encoded by this gene contains two motifs: a motif found in RAD52, a protein that functions in DNA double-strand breaks and homologous recombination, and an RNA recognition motif (RRM) that is not found in RAD52. The RAD52 motif region in RAD52 is important for protein function and may be involved in DNA binding or oligomerization. Alternatively spliced transcript variants encoding different isoforms have been reported. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG228198