

Product datasheet for **RG228207**

HSDL1 (NM_001146051) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	HSDL1
Synonyms:	SDRI2C3
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG228207 representing NM_001146051
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCTGCTGTTGACAGTTTCTACCTCTTGACAGGAAATCGCCAGGTCTTGCAATTGCTATATGGAAG
CTCTAGCTTTGGTTGGAGCCTGGTATACGGCCAGAAAAAGCATCACTGTCATCTGTGACTTTTACAGCCT
GATCAGGCTGCATTTTATCCCCGCCTGGGAGCAGAGCAGACTTGATCAAGCAGTATGGAAGATGGGCC
GTTGTCAGCGGTGCAACAGATGGGATTGAAAAGCCTACGCTGAAGAGTTAGCAAGCCGAGGTCTCAATA
TAATCCTGATTAGTCGGAACGAGGAGAAGTTGCAGTATTTCACTCAGCTGTCCGAGGACAAGCTCTGGGA
CATCATAAATGTGAACATTGCCGCCGCTAGTTTGTGTTGTTACCGGGAATGGTGGAGAGA
AAGAAAGGTGCCATCGTCACGATCTTTCTGGCTCCTGCTGCAAACCCACTCCTCAGCTGGCTGCATTTT
CTGCTTCTAAGGCTTATTTAGACCACTTCAGCAGAGCCTTGCAATATGAATATGCCTCTAAAGGAATCTT
TGTACAGAGTCTAATCCCTTTCTATGTAGCCACCAGCATGACAGCACCAGCACTTTCTGCACAGGTGC
TCGTGTTGGTGCCTTCGCCAAAAGTCTATGCACATCATGCTGTTTCTACTCTTGGGATTTCCAAAAGGA
CCACAGGATATTGGTCCATTCTATTCAGTTTCTTTTGCACAGTATATGCCTGAATGGCTCTGGGTGTG
GGGAGCAAATATTCTCAACCGTTCACTACGTAAGGAAGCCTTATCCTGCACAGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG228207 representing NM_001146051
Red=Cloning site Green=Tags(s)

MAAVDSFYLLYREIARSCNCYMEALALVGAWYTARKSITVICDFYSLIRLHFIPRLGSRADLIKQYGRWA
VYSGATDGIGKAYAEELASRLNIILISRNEEKLYFTQLSEDKLWDIINVNIAAASLMVHVLPGMVER
KKGAIVTISGSCCKPTQLAAFSAKAYLDHFSRALQYEQYASKGIFVQSLIPFYVATSMTPSNFLHRC
SWLVPSPKVYAHHAVSTLGISKRTTGYWSHSIQFLFAQYMPPEWLWWGANILNRSRLRKEALSCTA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001146051

ORF Size: 825 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001146051.1, NP_001139523.1</u>
RefSeq Size:	3503 bp
RefSeq ORF:	828 bp
Locus ID:	83693
UniProt ID:	<u>Q3SXM5</u>
Cytogenetics:	16q24.1
Protein Families:	Druggable Genome

Product images:



Circular map for RG228207