

Product datasheet for **RG228237**

LEF1 (NM_001166119) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: LEF1 (NM_001166119) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: LEF1
Synonyms: LEF-1; TCF1ALPHA; TCF7L3; TCF10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG228237 representing NM_001166119
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCAGAGGTGGCCAGACAAGCACAAACCTCTCAGGAGCCCTACCACGACAAGGCCAGAGAACACCCCG
 ATGACGGAAAGCATCCAGATGGAGGCCTCTACAACAAGGGACCCTCCTACTCGAGTTATTCGGGTACAT
 AATGATGCCAAATATGAATAACGACCCATACATGTCAAATGGATCTCTTTCTCCACCCATCCCGAGAACA
 TCAAATAAAGTGCCCGTGGTGCAGCCATCCCATGCGGTCCATCCTCTACCCCCCTCATCACTTACAGTG
 ACGAGCACTTTTCTCCAGGATCACACCCGTACACATCCCATCAGATGTCAACTCCAACAAGGCATGTC
 CAGACATCCTCCAGCTCCTGATATCCCTACTTTTTATCCCTTGCTCCGGGTGGTGTGGACAGATCACC
 CCACCTCTGGCTGGTTTTCCCATCATATGATTCCTGGTCCCTCCTGGTCCCCACACAACCTGGCATCCCTC
 ATCCAGCTATTGTAACACCTCAGGTCAAACAGGAACATCCCCACACTGACAGTGACCTAATGCACGTGAA
 GCCTCAGCATGAACAGAGAAAGGAGCAGGAGCCAAAAGACCTCACATTAAGAAGCCTCTGAATGCTTTT
 ATGTTATACATGAAAGAAATGAGAGCGAATGTCGTTGCTGAGTGTACTCTAAAAGAAAGTGCAGCTATCA
 ACCAGATTCTGGCAGAAGGTGGCATGCCCTCTCCCGTGAAGAGCAGGCTAAATATTATGAATTAGCACG
 GAAAGAAAAGACAGCTACATATGCAGCTTTATCCAGGCTGGTCTGCAAGAGACAATTATGGTAAGAAAAAG
 AAGAGGAAGAGAGAGAACTACAGGAATCTGCATCAGGTACAGGTCCAAGAATGACAGCTGCCTACATC


ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001166119.1 , NP_001159591.1
RefSeq Size:	2660 bp
RefSeq ORF:	912 bp
Locus ID:	51176
UniProt ID:	Q9UJU2
Cytogenetics:	4q25
Protein Families:	Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors
Protein Pathways:	Acute myeloid leukemia, Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Basal cell carcinoma, Colorectal cancer, Endometrial cancer, Melanogenesis, Pathways in cancer, Prostate cancer, Thyroid cancer, Wnt signaling pathway
Gene Summary:	This gene encodes a transcription factor belonging to a family of proteins that share homology with the high mobility group protein-1. The protein encoded by this gene can bind to a functionally important site in the T-cell receptor-alpha enhancer, thereby conferring maximal enhancer activity. This transcription factor is involved in the Wnt signaling pathway, and it may function in hair cell differentiation and follicle morphogenesis. Mutations in this gene have been found in somatic sebaceous tumors. This gene has also been linked to other cancers, including androgen-independent prostate cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]

Product images:



Circular map for RG228237