

Product datasheet for **RG228246**

OR56A5 (NM_001146033) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	OR56A5
Synonyms:	OR56A5P
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG228246 representing NM_001146033
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACATTACCCAGCAACAACCTCCACTTCCCCAGTCTTTGAATTCTTCCTCATTTGTTTCCCAGTTTCC
AGAGCTGGCAGCACTGGCTGTCTCTGCCCCCTCAGCCTCCTCTTCTCCTGGCCATGGGGGCAATGCCAC
CCTTCTGATCACCATCTATCTGGAAGCCTCTCTGCACCAGCCCCGTACTACCTGCTCAGCCTCCTCTCC
CTGCTGGACATCGTACTCTGCCTCACCCTCATCCCCAAGTCTGGCCATCTTCTGGTTTGACCTCAGAT
CAATCAGTTCCCTGCCTGCTTCTTCCAGGTGTTTCATCATGAACAGTTTTCTGACTATGGAGTCCCTGCAC
ATTCATGATCATGGCCTATGACCGCTATGTGGCCATCTGCAAGCCCCCTACAGTACTCATCCATCATCACT
GATCAATTTGTGCTAGGGCTGCCATCTTTGTTGTGGCCAGGAATGGCCTTCTTACTATGCCTATCCCCA
TACTTTCTTCTCGACTCAGATACTGTGCAGGACACATCATCAAGAACTGCATCTGTACTAACGTGTCTGT
GTCTAAACTCTTGTGATGACATCACCTTGAATCAGAGCTACCAGTTTGTATAGGTTGGACCTGCTG
GGCTCTGACCTCATCCTTATTGTTCTCTTACTTTTTATCTTGAAGAACTGTGCTAAGGATTAAGGGTG
AGGGAGATATGGCCAAAGCTCTAGGTACTTGTGGTCCCACTTCATCCTCATCCTCTTCTTACCACAGT
CCTGCTGGTTCTGGTCATCACTAACCTGGCCAGGAAGAGAATTCCTCCGGATGTCCCATCCTGCTCAAC
ATCCTGCACCACCTTATCCCCCAGCTCTGAACCCATTGTTTATGGTGTGAGAACCAAGGAGATCAAGC
AGGGAATCCAGAACCTGCTGAGGAGGTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG228246 representing NM_001146033
Red=Cloning site Green=Tags(s)

MTLPSNNSTSPVFEFFLICFPSFQSWQHWLSLPLSLLFLLAMGANATLLITIYLEASLHQPLYLLSLLS
LLDIVLCLTVIPKVLAIWFDLRSISFPACFLQVFIMNSFLTMESCTFMIMAYDRYVAICKPLQYSSIIT
DQFVARAAIFVVARNGLLTMPILSSRLRYCAGHIKNCICTNVSVSKLSCDDITLNQSYQFVIGWTL
GSDLILIVLSYFFILKTVLRIKGEQMAKALGTCGSHFILILFFTTLVLLVITNLARKRIPDPVILLN
ILHHLIPPALNPVIYGVRTKEIKQGIQNLLRRL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGATCGCCGGCGCCAGATCT

EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI

CAAGCTTAAGCTAGTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC ---

HindIII NheI RsrII MluI NotI XhoI GFP Tag

T R T R P L E M E S D - - -

--- GAA GAA AGA GTT TAA ACGGCCGCCGCGGAGCT

- - E E R V Stop

ACCN: NM_001146033

ORF Size: 939 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001146033.1</u> , <u>NP_001139505.1</u>
RefSeq Size:	942 bp
RefSeq ORF:	942 bp
Locus ID:	390084
UniProt ID:	<u>P0C7T3</u>
Cytogenetics:	11p15.4
Protein Pathways:	Olfactory transduction
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Mar 2009]

Product images:



Circular map for RG228246