

Product datasheet for **RG228265**

ZBT24 (ZBTB24) (NM_001164313) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZBT24
Synonyms:	BIF1; ICF2; PATZ2; ZNF450
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG228265 representing NM_001164313
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGGCAGAAACATCGCCAGAGCCTTCTGGGCAGCTTGTGTACACTCAGACGCTCACAGTGACACTGTGC
TGGCCAGTTTTGAGGATCAGAGGAAGAAAGGCTTCCTCTGTGACATTACTTTAATCGTGGAGAATGTACA
TTTCCGGGCCACAAAGCCTTACTTGCTGCCAGTAGTGAATACTTCTCAATGATGTTTGAGAAGAGGGG
GAAATCGGCCAATCCATTTATATGCTGGAAGGCATGGTGCAGACACCTTGGTATCCTGCTGGAATTTA
TCTACACAGGTTATCTCCATGCCAGTGAGAAAAGTACAGAACAATCCTGGCTACTGCTCAGTTCTTAAA
AGTCTATGACCTGGTAAAGGCTTACACAGACTTCCAAAATAATCATAGCTCCCCAAAGCCAACAACCTTTG
AACACTGCTGGTGGCCAGTGGTTGTTATCTCTAATAAGAAAAACGATCCTCCAAAGCGGAAACGGGGAA
GACCAAAAAAGTCAATACATTGCAGGAGGAGAAATCAGAACTGGCTGCAGAGGAAGAAATACAGTTAAG
AGTGAACAATTCAGTTCAGAAATAGACAAAACCTTTGTGGTTAAAGGAGACAGTGGTGTACTGAATGAGCAA
ATTGCAGCAAAAGAAAAGGAAGATCGGAGCCTACTTGTGAGCCAAGTAGAGAGGAGGAAATGCCAGTTG
AAAAAGATGAGAACTATGATCCCAAGACCGAGGATGGCCAGGCAAGCCAGAGTCGATACAGCAAGCGGAG
GATTTGGAGATCCGTCAAACCTTAAAGATTACAACTTGTGGGGATCAAGAGGACCATGGTTCAGCCAAG
AGGATCTGTGGAAGGAGAAAGCGCCCTGGAGGCCCTGAGGCCCGCTGTAAGAACTGTGGCAAGGTCTTTA
AGTACAATCACTTTTATAGCAATCCACCAGAGGAGCCACACAGGTAATGATGTTTCAAAGCTGATTGCAG
TGTGTTGCAGAATTGGGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG228265 representing NM_001164313

Red=Cloning site Green=Tags(s)

MAETSPESGQLVHSDAHS DTVLASFEDQRKKGFLCDITLIVENVHFRAHKALLAASSEYFSMMFAEEG
EIGQSIYMLEGMVADTFGILLEFIYTGYLHASEKSTEQILATAQFLKVYDLVKAYTDFQNNHSSPKPTTL
NTAGAPVVVISNKKNDPPKRKRGRPKKVNLTQEEKSELAAAAEIIQLRVNNSVQNRQNFVVKGD SGVLNEQ
IAAKEKEESEPTCEPSREEEMPVEKDENDPKTEDGQASQSRYSKRRIWRSVKLKDYLKLVGDQEDHGS AK
RICGRRKRPGGPEARCKDCGVFKYNHFLAIHQRSHTGNDVFKADCSVLQNWE


TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGATCGCCGGCGCCAGATCT

EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI

CAAGCTTAAGCTAGCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC --- ---

HindIII NheI RsrII MluI NotI XhoI GFP Tag

T R T R P L E M E S D - - -

--- GAA GAA AGA GTT TAA ACGGCCGCCGCGGAGCT

- - E E R V Stop

ACCN:

NM_001164313

ORF Size:

999 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001164313.1](#), [NP_001157785.1](#)

RefSeq Size: 1260 bp

RefSeq ORF: 1002 bp

Locus ID: 9841


UniProt ID: [O43167](#)

Cytogenetics: 6q21

Protein Families: Transcription Factors

Gene Summary: This gene encodes a protein similar to a protein in rodents which is induced by bone morphogenic protein 2 in vitro. [provided by RefSeq, Aug 2011]

Product images:



Circular map for RG228265