

Product datasheet for **RG228310**

DISC1 (NM_001164553) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: DISC1 (NM_001164553) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: DISC1
Synonyms: C1orf136; SCZD9
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG228310 representing NM_001164553
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCAGGCGGGGTCTCAGGGCGCCCCAGCCGCCCGCGGGCGGGCGGTGAGCCACCGCAGGCA
GCCGGATTGCTTACCACCTGCAGCGTCTTTTCGGAGGCGGGTGGCACGGAGCCGGGTACATGAG
AAGCTCGACAGGGCTGGGATCGGGTCTTTCCCGAGCAGTGGGCACACTGTTCCGGTCCAGGAGGG
GTGTCTGGCGAGGAGTCCCACTCGGAGTCCAGGGCCAGACAGTGTGGCCTTGACTCGAGAGGCCTCT
TGGTCCGAGCCCTGTTCCAAGAGTGCAGCAGCCCTACTGTGACCTCTGTGAGAGGAACCTCGGCGCA
CTTGGGATTCAGCTCAGAGGTGGCACCAGATTGCCTGACAGGCTTAGCTGGCCGTGTGGCCTGGGAGT
GCTGGGTGGCAGCAAGAGTTGCAGCCATGGATAGTTCTGAGACCCTGGACGCCAGCTGGGAGGCAGCCT
GCAGCGATGGAGCAAGGCGTGTCCGGGCAGCAGGCTCTCTGCCATCAGCAGAGTTGAGTAGCAACAGCTG
CAGCCCTGGCTGTGGCCCTGAGGTCCCCCAACCCTCTGGCTCTCACAGTGCCTTTACCTCAAGCTTT
AGCTTTATTCGGCTCTCGCTTGGCTCTGCCGGGAACGTGGAGAAGCAGAAGGCTGCCACCATCCAGAG
AGGCTGAGTCCCATTGCCAGAGCCCCAGGAGATGGGAGCCAAAGCTGCCAGCTTGGACGGGCTCACGA
GGACCCGCGATGTCTCTCGGCCCTTCAGTCTCTTGGCTACACGGGTCTCTGCAGACTTGGCCAGGCC
GCAAGGAACAGCTCCAGGCCAGAGCGTGACATGCATTCTTTACCAGACATGGACCCTGGCTCCTCCAGTT
CTCTGGATCCCTCACTGGCTGGCTGTGGTGGTGTGAGGAGCAGCGGCTCAGGGGATGCCACTCTTGGGA
CACCTGCTCAGGAAATGGGAGCCAGTGTGCGGGACTGCTGCTGAGAAACCGGAGGCAGATGGAGGTA
ATATCCTTAAGATTAACCTTCAGAACTTCAGGAAGTGCAGTTGAGAAATGATGATTATGATAAAGCTG
GAACCAATTGCTTTGGATCCACCATGGAAGCCTCGACATCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG228310 representing NM_001164553
Red=Cloning site Green=Tags(s)

MPGGGPQGAPAAAGGGGVSHRAGSRDCLPPAACFRRRRLARRPGYMRSSSTGPGIGFLSPAVGTLFRFPGG
 VSGEESHSESRARQCGLDSRGLLVRSVPVSKSAAAPT VTSVRGTSAHFGIQLRGGTRL PDRLSWPCGPGS
 AGWQQEF AAMDSSSETLDASWEAACSDGARRVRAAGSLPSAELSSNSCSPGCGPEVPPTPPGSHSAFTSSF
 SFIRLSLGSAGERGEAEGCPPSREAE SHCQSPQEMGAKAASLDGPHEDPRCLSRPF SLLATRV SADLAQA
 ARNSSRPERDMHSLPMDPGSSSLDPSLAGCGDGS SGGDAHSWDTLLRKWEPVLRDCLLRNRRQMEV
 ISLRKLKQLQEDAVENDDYDKAGTNCFGSTMEASTS


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001164553

ORF Size:	1161 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001164553.2
RefSeq Size:	2706 bp
RefSeq ORF:	1164 bp
Locus ID:	27185
Cytogenetics:	1q42.2
Gene Summary:	This gene encodes a protein with multiple coiled coil motifs which is located in the nucleus, cytoplasm and mitochondria. The protein is involved in neurite outgrowth and cortical development through its interaction with other proteins. This gene is disrupted in a t(1;11) (q42.1;q14.3) translocation which segregates with schizophrenia and related psychiatric disorders in a large Scottish family. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]