

Product datasheet for **RG228336**

PODNL1 (NM_001146255) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	PODNL1
Synonyms:	SLRR5B
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG228336 representing NM_001146255
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGGGAGGCCACCCAGTGGCCGAGCCTGCTGCTGCTCCTGCTGTTGCCGGGGCCCCGCCCGTCGCCG
GCTTGAAGACGCTGCCTTCCCCACCTGGGGGAGAGCTTGACGCCCTGCCCGGGCCTGTCCCCTGCC
CTGCTCCTGCCCGGAGTCGACACTGTGGACTGTGATGGCTTGGACCTTCGAGTGTTCCCGACAACATC
ACCAGAGCCGCTCAGCACCTCTCCCTGCAGAACACCAGCTCCAGGAATCCCTACAATGAGCTGTCCC
GCCTCAGTGGCCTGCGAACCTCAACCTCCACAACAACCTCATCTCTCCGAAGGCCTGCCTGACGAGGC
CTTCGAGTCCCTCACCCAGCTGCAGCACCTCTGCGTGGCTCACAACAAGAACAATCTCATCTCCAAGGTG
CCCCGAGGAGCCCTGAGCCGCCAGACTCAACTCCGTGAGCTCTACCTCCAGCACACCAGCTGACAGACA
GTGGCCTGGATGCCACCACCTTCAGCAAGCTGCATAGCCTTGAATACCTGGATCTCTCCACAACCAGCT
GACCACAGTGCCCGCCGGCCTGCCCGGACCCTGGCTATCCTGCACCTGGGCCGCAACCGCATCCGGCAG
GTGGAGGCGGCTCGGCTGCACGGGGCGCGTGGTCTGCGCTATTTGTTGCTGCAGCACACCAGCTGGGGA
GCTCAGGGCTGCCCGCCGGGCTCTGCGGCCGCTGCGGGCCCTGCACACGCTGCACCTCTATGGCAATGG
GCTGGACCGCGTGCCTCCAGCCCTGCCCGCCGCTGCGTGCCTGGTGTGCCCAACAACACGTGGCC
GCGCTGGGTGCCCGTGACCTGGTGCACACCGGGCCTGACGGAGCTTAACCTGGCCTATAACCGCCTGG
CCAGCGCCCGTGTGCACACCAGGGCCTTCCGCCGTTGCGTGCCCTGCGCAGCCTCGACCTGGCAGGGA
TCAGCTAACCCGGCTGCCATGGGCTGCCACTGGCCTGCGCACCTGCAGCTGCAACGCAACCAGCTG
CGGATGCTCGAGCCGAGCCTCTGGCCGGCTGGACCAACTGCGGGAGCTCAGCCTGGCGACAACCGGC
TCCGGGTGCGGACATCGGGCCAGGCACCTGGCATGAGCTCCAAGCCCTCAGGTGAGGCACAGGCTGGT
TAGCCACACTGTCCCAGGGCCCTCCATCCCCCTGCCTGCCCTGCCACGTCCCAAACATTCTAGTTAGC
TGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG228336 representing NM_001146255
Red=Cloning site Green=Tags(s)

MGRPTQWPSLLLLLLPGPPVAGLEDAAFPHLGESLQPLPRACPLRCSCPRVDTVDCDGLDLRVFPDNI
TRAAQHLSLQNNQLQELPYNELSRLSGLRTLNLHNNLISSEGLPDEAFESLTQLQHLCAHNKNNLISKV
PRGALSRQTQLRELYLQHNQLTDSGLDATTF SKLHSLEYLDL SHNQLTTVPAGLPRTLAILHLGRNIRQ
VEAARLHGARGRLRYLLQHNQLGSSGLPAGALRPLRGLHTLHLYGNGLDRVPPALPRRLRALVLPNHVA
ALGARDLVATPGLTELNLAYNRLASARVHHRAFRRLRALRSLDLAGNQLTRLPMGLPTGLRTLQLQRNQL
RMLEPEPLAGLDQLRELSLAHNRLRVGDIGPGTWHELQALQVRHRLVSHTVPRAPSPCLPCHVPNILVS
W

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001146255

ORF Size: 1263 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001146255.1](#), [NP_001139727.1](#)

RefSeq Size: 3041 bp


RefSeq ORF: 1266 bp

Locus ID: 79883

UniProt ID: [Q6PEZ8](#)

Cytogenetics: 19p13.12

Product images:



Circular map for RG228336