

Product datasheet for **RG228344**

DISC1 (NM_001164550) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: DISC1 (NM_001164550) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: DISC1
Synonyms: C1orf136; SCZD9
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG228344 representing NM_001164550
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCAGGCGGGGTCTCAGGGCGCCCCAGCCGCCCGCGGGCGGGCGGTGAGCCACCGCAGGCA
 GCCGGATTGCTTACCACCTGCAGCGTGTCTTCGGAGGCGGGTGGCACGGAGCCGGGTACATGAG
 AAGCTCGACAGGGCTGGGATCGGGTTCCTTTCCCGAGCAGTGGGCACACTGTTCCGGTCCAGGAGGG
 GTGTCTGGCGAGGAGTCCACCCTCGGAGTCCAGGGCCAGACAGTGTGGCCTTGACTCGAGAGGCCTCT
 TGGTCCGGAGCCCTGTTCCAAGAGTGCAGCAGCCCTACTGTGACCTCTGTGAGAGGAACCTCGGCGCA
 CTTTGGGATTCAGCTCAGAGGTGGCACCAGATTGCCTGACAGGCTTAGCTGGCCGTGTGGCCTGGGAGT
 GCTGGGTGGCAGCAAGAGTTTGCAGCCATGGATAGTTCTGAGACCCTGGACGCCAGCTGGGAGGCAGCCT
 GCAGCGATGGAGCAAGGCGTGTCCGGGCAGCAGGCTCTCTGCCATCAGCAGAGTTGAGTAGCAACAGCTG
 CAGCCCTGGCTGTGGCCCTGAGGTCCCCCAACCCTCTGGCTCTCACAGTGCCTTTACCTCAAGCTTT
 AGCTTTATTCGGCTCTCGCTTGGCTGTCCGGGGAACGTGGAGAAGCAGAAGGCTGCCACCATCCAGAG
 AGGCTGAGTCCCATTGCCAGAGCCCCAGGAGATGGGAGCCAAAGCTGCCAGTTGGACGGGCTCACGA
 GGACCCGCGATGTCTCTCGGCCCTTCAGTCTCTTGGCTACACGGGTCTCTGCAGACTTGGCCCAGGCC
 GCAAGGAACAGCTCCAGGCCAGAGCGTGACATGCATTCTTTACCAGACATGGACCCTGGCTCCTCCAGTT
 CTCTGGATCCCTCACTGGCTGGCTGTGGTGGTGTGGGAGCAGCGGCTCAGGGGATGCCACTCTTGGGA
 CACCCTGCTCAGGAAATGGGAGCCAGTGTGCGGGACTGCCTGCTGAGAAACCGGAGGCAGATGGAGGTA
 ATATCCTTAAGATTAACCTTTCAGAACTTCAGGAAGTGCAGTTGAGAAATGATGATTATGATAAAGGAA
 GTCAGTGAAGGGTTACATTTTCATTTGGGGTGAACAGCAACTGTGGATACGGATCATGAAAATTGACAA
 TGGAAAATGGGCATGCCATTCTGGAACTTTTCCAGCTTTTTTCCAGAGCCAGCTGGAACCAATTGCTTT
 GGATCCACCATGGAAGCCTCGACATCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG228344 representing NM_001164550
 Red=Cloning site Green=Tags(s)


MPGGGPQGAPAAAGGGGVSHRAGSRDCLPPAACFRRRRLARRPGYMRSSSTGPGIGFLSPAVGTLFRFPGG
 VSGEESHSESRARQCGLDSRGLLVRSVPVSKSAAAPTVTSVRGTSAHFGIQLRGGTRL PDRLSWPCGPGS
 AGWQQEF AAMDSSETLDASWEAACSDGARRVRAAGSLPSAELSSNSCSPGCGPEVPPTPPGSHSAFTSSF
 SFIRLSLGSAGERGEAEGCPPSREAE SHCQSPQEMGAKAASLDGPHEDPRCLSRPF SLLATRVSADLAQA
 ARNSSRPERDMHSLPMDPGSSSLDPSLAGCGDGS SGGDAHSWDTLLRKWEPVLRDCLLRNRRQMEV
 ISLRKLKQLQEDAVENDDYDKGSHWKGYIFIWGEQQLWIRIMKIDNGKWACHSGTFPSFFPEPAGTNCF
 GSTMEASTS

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001164550

ORF Size:	1287 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001164550.2
RefSeq Size:	2832 bp
RefSeq ORF:	1290 bp
Locus ID:	27185
Cytogenetics:	1q42.2
Gene Summary:	This gene encodes a protein with multiple coiled coil motifs which is located in the nucleus, cytoplasm and mitochondria. The protein is involved in neurite outgrowth and cortical development through its interaction with other proteins. This gene is disrupted in a t(1;11) (q42.1;q14.3) translocation which segregates with schizophrenia and related psychiatric disorders in a large Scottish family. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]