

Product datasheet for **RG228351**

CCDC166 (NM_001162914) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	CCDC166
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG228351 representing NM_001162914
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGCCCAAGAAAAAGCGGGGCCGAGCGCAGGGAGCCAGCCAGGTGGTGCGGCGGCCGCGGGTGCCG
AGCAGCCGCTATCGGAGCGCGCAGTACCTGCAACGCGAACACGCGCTGCTCTCGGAGCAGCTGGACAC
CTGCGAGGAGAGCGTTGACCAGGTGCTGCGAGAGAACGCCTTCCTGGACCGCGAGGCGCTGCGCCTGCGC
GAGGAGAACCGGCTCTACGCCAGCTACGTGAGCGCGCGCGCCAGCGCTGCGCCAAAGCCATCGTCCGGC
TGGACGAGCAGAACCAGCGTGGACCTAGCGCAGATCCACTGGCAGCGGGCGGAAGTGGCTCGCTCTACCA
CGGGCGCGAGGACGGGTGCGCGCGCAGCTGTTGGAGATGGAGCGCGCGGCAGATGGCACAGCAG
GTGCAAGAGCTGCAGCCCTACAAGGTGCTGCAGCTGGAGCAGCTGGCCCCGATCCGGGCGCTGGAGCGCG
AGCTGCTGCATATGCGCGTGGAGCACAGCAGCTGCTCCACGCGTGAAGCGGCGCTTCCTGGAGGACAA
GGCGGCCTTCGAGCGCGAGGCGCGTCAAGCGCTGCACTGGCGCGCGCGCGGAGCGGGAGGCGGTG
CGCGCGCTTGTGGCGCACAGCAGGCCATCAAAGCGGACAAAGGACGCTGCGGCGAGGAGCTGCTGCTGC
TGCTCCGCCGACCCAGCTGCTGCACCACAGCGCGCGCAGCTGCTGGAACAGCGCGAGCAACTGCACCG
CGAGCACGAGGACAGCGGGACCTGGCGCGTGTGCACGGCTGGTGCAGGGGCCCCGAGGCCCCGCCG
CTCTGGGAGCGGCCGCGCTTCTCCAGCCACCTCGCGCCCCGGGTCAATAGCGCCCCCATAAGCCCGT
CGCGCGCGGCTCCAGACCCGCTGTGGTACCGTGCAGCGCGGCCCCCGGGCTCGTCCGTGGTCCC
GTCGCGGAGGCTTCCCGGGTCCCCTCGTTGGTGTATCGAGCATGGACTCCAGGGTCCCATCGCTGGCC
ACTTCAAAGGTGGGATCCCGGATGCCATCCCTGACCGCTTCGCGCGCAGGCTCCCGGGCCCTGTGCTGG
TGCAGTCACTCGAGGGCTCTGGGATCTCTCGGGTCTCACCAGGGTGTCTCACAGGACACTCTCCG
CTCCACGAAGTCTGGCCCCAAGCTTCTCTGGCTGTCCCGGGATCGGGATCCGGCTTTCTCCCCCA
CAGTCGGAGGACAGCGTGAACGCTGAAGCTGCGGCAGAAGCTCCCGGGTAGAGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG228351 representing NM_001162914

Red=Cloning site Green=Tags(s)


MAPKKKRGPSAGSQPGGAAAAGAEQPLSERAQYLQREHALLSEQLDTCESVDQVLRENAFLDREALRLR
EENRLYASYVSARAQRCAKAIIVRLDEQNRVDLAQIHWQRAELASLYHGREDGVRAQLLEMEARAAQMAQQ
VQELQPYKVLQLEQLARIRALERELLHMRVEHTQLLHRVKRRFLDKAAFEREARQVRVQSLARRAEREAV
RALVAHTQAIKADNGRLRQELLLLLRRRTQLLHHTRRQLLEQREQLHREHEDTRDLARVHGWLRRGPGGPP
LWERPAFSQPTSRPGSLAAPISPSRAASQTPSVVPSRAAPRASSVPSREASRVPSLVLSMDSRVPSLA
TSKVGSRMPSLTASRAGSRALSLVQSLEGGSSISSGSSPRVSSQDTRLRSTKSGPKLLSGLSRDRDPALLPP
QSEDSVNAEAAEASPGRA

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001162914

ORF Size:

1317 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001162914.1](#), [NP_001156386.1](#)

RefSeq Size: 1320 bp


RefSeq ORF: 1320 bp

Locus ID: 100130274

UniProt ID: [P0CW27](#)

Cytogenetics: 8q24.3

Product images:



Circular map for RG228351